ETCHING GEL

Etching Gel For Enamel And Dentin

en

Indications:
Etching of enamel surface or dentine before the preparation of a composite or bonded amalgam filling, before the preparation of a sandwich composite filling with a glass ionomer cement, before adhesive fixation of an aesthetic facet, either ceramic or composite, before the fixation of an inlay or onlay, the root pin, a crown or a bridge with a composite cement.

Application:
Connect the application tips to the syringe and check that it moves freely. The etching gel is applied to the enamel edges directly with the tip and left to act for 20-30 seconds. If the gel is used with the Retensin adhesive, only the enamel is etched. If dentine is to be etched ("total etching"), the etching gel is applied to the dentine 15-20 seconds later so that the overall action time does not exceed 10 seconds. The etching gel is removed with a sufficient stream of water (for 15 seconds) and the cavity surface is dried with a stream of air which is free of oil particles, so that the enamel edges become white and the dentine is not overried, i.e. does not lose its natural moisture. If the etched surface is contaminated e.g. with saliva, the etching must be repeated for 10 seconds and the surface washed in water and dried in the usual manner. Extending the etching time is recommended for teeth with a higher content of fluorine or for temporary teeth. The selected adhesive system is applied according to the instruction sheet.

Notice:
The use of protective glasses and contact-free technique is recommended when working with the etching gel. The surrounding soft tissues must be protected.

Warnings:
H314 Causes severe skin burns and eye damage.
P280 Wear protective gloves/protective clothing/eye protection/face protection. P301+310 IF SWALLOWED: Immediately call a POISON CENTER or doctor/physician. P305+351+338 IF IN EYES: Rinse cautiously with water for several minutes. Remove contact lenses, if present and easy to do. Continue rinsing.

Desinfection:
The used tips should be removed and disposed. After closing the syringe with a stopper, disinfect its surface with a solution of a common surface disinfectant.

Storage:
To be stored in a dry place, with temperatures below 23 °C, away from heat sources.

Expiry period:
2 years

Packaging:
3.5 g PE syringe + 10 application tips (supplied as a separate package)
30 g Etching Gel in a syringe + 20 application tips + 2 empty syringes

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