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EN Instructions for use IFU 2311 READ CAREFULLY!

Carbide bur

SAFETY PRECAUTIONS

The product is intended for professional use by dentists with the appropriate qualification.

TARGET PATIENT POPULATION

There is no age or other restriction for this instrument, provided that the indication criteria are met.

1. Description

- Carbide burs are dental rotary instruments. They are either manufactured from a single piece of tungsten carbide and steel or from a tungsten carbide working part brazed onto a steel shank. A shank of some variants of carbide burs is gold plated, which increases corrosion resistance. They are available with FG (friction grip), HP (hand piece) and RA (right angle) shanks which are designed to fit into the chuck of dental handpieces. The working part is equipped with 6 to 32 blades. A coloured groove on the shank makes identification of the adequate the number of blades on the working part easier.

2. Material

- Material that come into direct contact with patient:
 - tungsten carbide or stainless steel 1.4034

Material	Graphic sign	ISO Code
Tungsten carbide	TC	500
Ocel 1.4034	SSt	330

3. Intended Purpose

- Carbide Bur is intended for cutting and polishing of hard structures in the oral cavity such as teeth, bone, amalgam, composite, hard metals, polymers and glass ionomer cement.

4. Indications

- Cavity preparation.
- Hard fillings removal.
- Finishing and polishing of fillings.
- Cutting of bones, metals, polymers, filling.
- Finishing of polymers.
- Trepanation of crown restorations.

5. Contraindications

- Nickel allergy. Carbide bur contains a small amount of nickel, which may cause an allergic reaction in hypersensitive patients.

6. Instrument selection

- The dentist selects the device based on the required application.
- Colour marking of the devices indicated the number of blades of working part.

Colour	Number of blades
Blue	8
Red	8–12
Yellow	16–20
White	30

7. Compatibility

- The instruments are used in conjunction with a dental handpieces.
- The types of carbide burs according to the shank are as follows:

Shank type 1: Carbide burs RA (right angle), RA XL (right angle extra long). Shank length is standard or extra long.

Shank type 2: Carbide burs HP (hand piece) handpieces. Shank length is standard.

Shank type 3: Carbide burs SHORT FG (short friction grip), FG (friction grip), FG L (friction grip long), FG XL (friction grip extra long). The different shank lengths are distinguished by the markings SHORT, standard without marking, L and XL.

Shank type	Diameter [mm]	Length [mm]	ISO Code
FG	1.6	19	314
FG SHORT	1.6	16.5	313
FG L	1.6	21	315
FG XL	1.6	25	316
RA	2.35	25	204
RA XL	2.35	34	206
HP	2.35	44.5	104

- The shank type of a carbon bur is indicated on its label. Marking of shank type is a part of the device name or can be identified from the second three digits of the 15-digit ISO code.

8. Procedure for use

- Remove the instruments from the stand or cassette.
- The selection of instruments (shape, size, number of blades, material) depends primarily on the planned application. Improperly selected instrument may lead to an unsatisfactory result.
- Insert the instrument as deeply as possible into the chuck.
- Clamp the instrument in the dental handpiece and check the clamping strength. An instrument that is not secured properly in the handpiece may cause injury to the patient or the dentist.
- Select suitable speed according to recommended speed table.

Size	Recommended speed (RPM)	Maximum speed (RPM)
005–006	60000–120000	450000
007–008	45000–50000	450000
009–010	35000–70000	450000
012–014	25000–50000	450000
016–018	20000–40000	450000
021–023	15000–30000	450000
025–027	12000–25000	160000

- Set the handpiece at the required working speed before applying the instrument to the treated area.
- Apply a maximum operating force of 0.4 [N].
- Proper cooling is required while using the device (min. 50 ml/min).
- Place instruments into the cleaning box after using.

9. Undesirable side-effects

- Inadequate cooling may generate undesirable heat and cause patient discomfort, tooth or tissue necrosis, or patient burns.
- Excessive pressure on the bur could cause undesirable heat or may cause the bur to fail and cause patient or dentist injury.

10. Precautionary measures

- The instruments are strictly intended for professional use by dentists.
- Eye protection or protection shield must be worn to protect against eject particles. Surgical masks must be worn to avoid inhalation of any aerosol or dust generated and possible cross-infection.

11. Information for the patient

- Not applicable.

12. Warnings

- Do not use the device if the packaging is damaged.
- Never exceed the maximum permissible speed. Not observing the maximum speed will result in the increased risk of injury
- Attention should be paid to the speed of work (RPM). Some types of carbide burs have maximally speed of work lower than is recommended in the speed table, this information is stated on the device label. Therefore always refer to the product label for the Maximum RPM.
- Clean and sterilize non-sterile burs in accordance with the validated procedures provided below prior to first use and prior to each reuse.
- Do not apply excessive pressure on the bur as this could cause undesirable heat or may cause the bur to fail and cause patient or user injury.
- Cleaning agents with chlorine or chloride as the active ingredient are corrosive to stainless steel and must not be used. Cleaning agents with neutral pH are recommended.
- Patients identified as at-risk for Creuzfeldt-Jacob disease (CJD) and related infection should be treated with single-use instruments. Properly dispose of instruments and implants used, or suspected of use, on a patient with CJD after surgery and/or follow current national recommendation.

13. Usage restriction

- Carbide burs are reusable.
- The device lifetime is limited. Maximally number of cycle cannot be set in advance. The device lifetime is given by its wear or damage. The signs of wear may include but they not limited on: corrosion, colour change, deformation, blades damage and breakage. Devices should be inspected prior each use to ensure full functionality.
- Worn or damaged devices must never be reused and must be dispose in compliance to local policies.

INSTRUCTIONS

14. Point of use

- Remove gross soil from the devices after use. A drying of gross soil onto the device surface should be avoided.
- Than place a devices into cleaning box as soon as possible.

15. Storage and transport

- Contaminated devices and unused devices should not come into contact during transport prevent contamination.

16. Preparation for cleaning

- Devices are supplied in protective packaging, which is not intended for sterilising. Remove devices out of packaging. Dispose a packaging.
- The device must be placed into disinfection solution as soon as is reasonable to avoid dry of gross soil on device.
- If gross soil cannot be removed in standard manner a soft brushes can be carefully used. Clean the device so that was not damage.
- Care should be taken in the handling and cleaning of sharp devices. These are recommended to be cleaned separately to reduce risks of injury.

17. Cleaning and disinfection – manual

- Wash the devices manually, under tap water and then in fresh cleaning and disinfecting solution in an ultrasonic bath, then rinse again under tap water.

- 17.2. A use of a specific detergent and disinfectant is not prescribed.
- 17.3. The washing agent and disinfectant have to be used in accordance with the instructions of the manufacturer. Follow prescribed concentration, temperature, water quality and exposure time.
- 17.4. Remove devices after recommended exposure time and wash under the tap water. **Note:** fresh solution is one that is newly mixed and clean.

18. Cleaning and disinfection – automatized

- 18.1. The validated automated cleaning and disinfection procedure must always be preferred to manual cleaning because of the higher safety of the process. A good cleaning is also a requirement for successful sterilization.
- 18.2. We recommend washing devices using automatic washers with thermal disinfection.
- 18.3. Place devices into the washing agent carefully to avoid moving devices freely during washing process or their overlapping which could damage their surface.
- 18.4. A use of a specific detergent and disinfectant is not prescribed.
- 18.5. The washing agent and disinfectant have to be used in accordance with the instructions of the manufacturer. Follow prescribed concentration, temperature, water quality and action time when using.

19. Thermal disinfection

- 19.1. Disinfect at a minimum temperature of 90 °C for at least 5 minutes in case of cleaning in washer (according to 15883).

20. Drying

- 20.1. If drying cycle is not included in the washer:
 - dry each device thoroughly.

21. Service

- 21.1. Not applicable.

22. Inspection

- 22.1. Devices have to be visually inspected after treatment and prior sterilizing.
- 22.2. Devices that have visibly damaged cutting edges, corrosion, cracks, deposits or tissue residues that cannot be removed must be discarded.

23. Packaging

- 23.1. The cleaned and dry devices put in their place in a stand or cassette intended for carbide bur sterilisation.
- 23.2. Stand-alone instrument may be packaged into sterilisation pouch to ISO 11607.

24. Sterilization

- 24.1. Sterilize cleaned and disinfected devices only.
- 24.2. Sterilize with wet heat in a steam sterilizer at a temperature of 134 °C and the time of sterilization of 7 min. and pressure of 304 kPa.
- 24.3. Sterilize stand-alone stand with instruments and cassette. Sterilize open cassette.

25. Storage

- 25.1. The products should be stored in their protective packaging in a clean and dry room, off the direct sunlight.

26. Other information

- 26.1. The cleaning process with thermal disinfection intended by the manufacturer has been validated. The manufacturer cannot guarantee the result if washing method other than those recommended is, used. When using another cleaning method, it is recommended that this procedure be validated by the supplier of the relevant cleaning equipment.
- 26.2. Instructions for cleaning and sterilization are provided in accordance with the standards of ISO 15883, ISO 17664, ISO 17665-1 and they are validated by the manufacturer MEDIN, a. s., as suitable for repeated preparation of a single-use medical device. If another procedure is used, the manufacturer cannot guarantee that the procedure has provided the required level of sterilization certainty of the device, therefore the user is responsible for achieving the required state. If the user does not follow the specified preparation procedure, he must ensure that the processing, as done, has achieved the desired result. This requires verification and / or validation and ongoing monitoring of the process.
- 26.3. If a serious adverse event is suspected in connection with the use of a medical device, report it to the manufacturer and to the national competent authority responsible for vigilance. You can find contact information for the manufacturer at **www.medin.cz** or communicate the suspicion directly to your contact person on behalf of MEDIN, a. s.

27. Device disposal

- 27.1. All worn and damaged devices, and devices that are contaminated with human blood, tissues, secretions or body fluids, pose a biological hazard and are considered infectious waste. Dispose these products according to local regulations for infectious biological medical waste.
- 27.2. Dispose demonstrably non-contaminated material, such as non-contaminated packaging, according to local regulations for medical waste, the collection and disposal of which are not subject to special requirements with regard to infection prevention.
- 27.3. Handle contaminated sharp objects with particular care. Injuries caused by contaminated objects may result in the transmission of particularly blood-borne pathogens to the patient or user. All contaminated sharps that carry a potential risk of injury and infection require special precautions to prevent injury during handling and must be placed in a rigid container designed for sharps waste. When disposing of them, follow the locally applicable regulations for the disposal of sharps from the healthcare sector.

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