- 16.4. Removal of hard-to-remove tissue residues from the instrument can be performed with compressed air and subsequent rinsing with clean distilled or demineralized water.

17. Cleaning and disinfection – manual

- 17.1. Wash the instruments manually, under running water and then in fresh cleaning and disinfecting solution in an ultrasonic bath, then rinse again under running water.
- 17.2. The use of a specific detergent and disinfectant is not prescribed.
- 17.3. The washing agent and disinfectant have to be used in accordance with the instructions of their manufacturer. Follow prescribed concentration, temperature, water quality and exposure time. High concentrations or prolonged exposure can damage the instrument.
- 17.4. Remove the instrument after recommended exposure time and wash under the tap water. **Note:** fresh solution is one that is newly mixed and clean.
- 17.5. Instruments can be disinfected with 5 % sodium hypochlorite solution NaClO, but not longer than 5 minutes.
- 17.6. Do not use products containing phenol and hydrogen peroxide H2O2 for disinfection.

18. Cleaning and disinfection – automatized

- 18.1. In the case of a choice, the validated automated cleaning and disinfection process takes precedence over the manual one, due to higher process safety. Quality cleaning is also a prerequisite for successful sterilization.
- 18.2. We recommend washing the instruments before sterilization using automatic dishwashers with the possibility of thermal disinfection.
- 18.3. Place the instruments in the dishwasher so that they cannot move freely or overlap during the washing process and their surface cannot be damaged.
- 18.4. The use of a specific detergent and disinfectant is not prescribed.
- 18.5. When using detergent and disinfectant, follow the instructions of the appropriate sterilizer manufacturer, observing the prescribed concentration, solution temperature, water quality and exposure time.

19. Thermal disinfection

- 19.1. Disinfect at a minimum temperature of 90 °C for at least 5 minutes in case of cleaning in washer.

20. Drying

- 20.1. If drying cycle is not included in the manual washer, dry each instrument thoroughly.

21. Service

- 21.1. Does not apply.

22. Inspection

- 22.1. The instruments must be visually inspected for condition and function before sterilization. Instruments should be undamaged without signs of corrosion and without residues of contaminant. Instruments showing wear or damage must be discarded and disposed of, re-use is not permitted.
- 22.2. Instrument with damaged surface shall be discarded.

23. Packaging

- 23.1. Store cleaned and dry instruments in bowls or sieves intended for the sterilization method.
- 23.2. Separate instruments can be packed in the sterilization packaging that meets the requirements of ISO 11607.

24. Sterilization

- 24.1. Sterilization parameters apply only to instruments that are sufficiently cleaned according to the above instructions.
- 24.2. Sterilize the instruments with wet heat in a steam sterilizer at a temperature of 134 °C and the time of sterilization of 7 min. and pressure of 304 kPa.
- 24.3. Sterilize instrument sieves and cassettes separately. Sterilize the cassettes open.

25. Storage

- 25.1. Store packaged instruments in a dry, clean environment out of direct sunlight.

26. Other information

- 26.1. The process of washing with thermal disinfection has been verified by the MEDIN company. In case of use of some other methods of the washing procedure than recommended, the MEDIN company does not guarantee the result. When using some other method of washing than recommended, the MEDIN company recommends this process to be validated by the supplier of such washing facility.
- 26.2. The above instructions for cleaning and sterilization are in accordance with the norms and standards of ISO 15883, ISO 17664, ISO 17665-1 and have been verified for the preparation of reusable non-sterile medical devices before use. The responsibility of the processor is achieve the desired result by cleaning disinfection and sterilization. In the event that the processor does not comply with the above procedure, it should properly assess the impacts in terms of effectiveness and possible adverse consequences. This requires verification and / or validation and ongoing monitoring of the process.
- 26.3. The user is responsible for providing the necessary materials, equipment and training of personnel in accordance with these prescribed procedures in order to achieve the desired result.
- 26.4. If a serious adverse event is suspected in connection with the use of a medical device, report it to the manufacturer and to the national competent authority responsible for vigilance. You can find contact information for the manufacturer at **www.medin.cz** or communicate the suspicion directly to your contact person on behalf of MEDIN, a. s.

27. Disposal

- 27.1. All contaminated sharp devices that can damage the skin and have a potential risk of injury and infection require special precautions to prevent injury during manipulation. These devices must be placed into in a solid container intended for sharpened waste and when disposing of them, follow the locally valid regulations for the disposal of sharps from healthcare.
- 27.2. Dispose of demonstrably uncontaminated material, such as uncontaminated packaging, in accordance with local regulations for medical waste, the collection and disposal of which are not subject to special requirements with regard to infection prevention.

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