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Vrták

1. Popis

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EN Instructions for use IFU 1111 READ CAREFULLY!

Drill

1. Description

- A non-sterile metal shaft with a V-shaped blunt cutting edge designed for rotation for cutting into a bone to form a hole of the same size as the shaft diameter. This device, commonly known as the "helical drill", has a cylindrical design with a two-edged cutting head which is formed from two spiral grooves extending along the introduction part of the drill bit. It is usually inserted into a mechanical surgical instrument or hand surgical drill at its proximal end. It is typically made of high quality stainless steel and is available in different sizes. It is a reusable medical device.

2. Material(s)

- Stainless steels according to ČSN EN ISO 7153-1.

3. Intended purpose for use

- Drills are rotary surgical instruments that are used to drill a hole into the bone and hard tissues during surgical procedures. The clamping part is designed for clamping into various types of drills. Some drills are provided with marks for the surgeon's orientation about the depth of the drilling.

4. Indications

- They are used in surgery or traumatology to:
 - create a hole in the bone for the self-tapping screw introduction
 - create a hole for the cannulated screw introduction
 - create a guide hole in the bone at the implant or bone graft introduction

5. Contraindications

- Are not determined.

6. Choice of the drill

- The surgeon has to choose an appropriate type end size of the drill according to the surgical technique.
 - The surgeon has to choose in the other cases an appropriate drill, which best corresponds with the surgical procedure needs.

7. Compatibility

- Drill bits are connected to surgical drills which allow attachment with round, triangular, hexagonal shank and AO quick coupling.

- Check before its use, if its shank corresponds with the chuck of a drill.

8. Procedure of use

- Connect the drill bit into the drill and create the hole needed using an appropriate pressure.Use the protection sleeves. The drill is reusable.

10. Precautionary measures

- The surgeon must be familiar with the surgical technique.
 - Protect the drills from deformation and scratches. Avoid contact with other metallic materials and chemicals when handling them.
 - The drills are provided with a blade protection before their expedition and for further handling before their use.
 - Always check the quality of the blades visually, never by touching it.
 - When the instrument is dull, let it sharpen.
 - Before its use always check visually if the drill is not bent.

11. Information for the patients

- Instruct the patient that he/she must strictly follow the instructions of the attending doctor.

12. Warning

- Caution! Drills are sharp!
 - When using a drill which has a drill bit attached, avoid touching the drill bit outside of the intended purpose. Stainless steel surgical instruments cannot be inserted into the physiological saline (NaCl solution), because a longer contact causes corrosion such as a pitting corrosion and stress corrosion.
 - Cleaning solutions in which hydrogen peroxide is added and/or solutions for cleaning with higher alkalinity may lead to a change in colour. This could cause damage to the identification of the instruments.
 - Do not use any abrasive cleaning materials such as steel wool or abrasive detergents to clean the instruments. If possible, always clean the instruments and sieves in the washing and disinfecting facility for instruments according to the standard ISO 15883.
 - The long, narrow cannulas, blind holes and complex parts require a special attention when cleaning.
 - Sterilising parameters only apply to products which are sufficiently cleaned.
 - Sterilisation parameters apply only to correctly assembled, maintained, calibrated and adequate equipment for processing in accordance with ISO 15883 and 17665 standards.
 - The products are delivered NON-STERILE.
 - Remove the instruments from their packaging before first use and remove covers from the sharp ends. These covers are not intended to be sterilized. Place the covers into the waste sorted in accordance with the Catalogue of waste – 18 01 07 – chemicals that do not contain hazardous substances. It is necessary to return the instrument to the manufacturer for repackaging and control in case of impaired packaging. Dispose the package according to the recycling symbols.
 - Perform cleaning, disinfection and sterilization according to the regulation before first use (points 14.–24. of these instructions).
 - Clean and sterilize the instruments as soon as possible after their use.

13. Restriction on repeated use

- The lifetime of a drill is determined by a wear and a damage due to their use (drilling). The manufacturer guarantees safe use for drilling of 30 holes. The drill must be resharpened according to point 22.2. or must be disposed according to point 27 after 30 holes drilled.
 - Repeated cycles of cleaning, disinfecting and sterilization have minimal impact to instruments.

INSTRUCTIONS

14. Place of use

- Wipe the blood and any remaining tissues from the device during the surgery before they are allowed to dry to the surface.
 - Rinse the instruments with openings by the sterile or purified water before the impurities and/or remaining tissues are allowed to dry inside.
 - Soiled instruments have to be separated from the clean ones to avoid contamination of the unused instruments, staff and surroundings. Contaminated instruments must be separated from medical waste.
 - The instruments should be covered with the drape dampened in the sterile or purified water to avoid drying of blood or debris.

15. Storage and transport

- Soiled devices have to be transported separately from the clean ones to avoid their contamination.

16. Preparation for cleaning

- Open the instruments with toothed or threaded latches and hinges. Separate the pointed and sharp instruments. Wipe the oil from the locks using a cloth at the new instruments. A larger amount of oil in the locks of the new surgical could lead to the formation of spots on the instrument surface after sterilization.
 - Remove the sharp removable parts of the devices prior the manual cleaning or place them onto the separated tray.
 - Disassemble or loosen all connections before processing of the instruments. You will get further detailed information about the dismantling from your local sales representative.

17. Cleaning and disinfection – manual

- Place MEDIN instruments into the washing agent carefully to avoid moving of them freely during the washing process or their overlapping which could damage their surface.
 - Wash the instruments manually, remove coarse dirt mechanically with brush in fresh cleaning and disinfection solution in an ultrasonic bath, and rinse under tap water.
 - The MEDIN company does not prescribe the use of a specific detergent and disinfectant.
 - The washing agent and disinfectant have to be used in accordance with the instructions of the manufacturer. Keep prescribed concentration, temperature, water quality and exposure time.
 - Remove the drill after recommended exposure time and wash under the tap water. **Note:** fresh solution is one that is newly mixed and dean.

18. Cleaning and disinfection – automatized

- The validated machine cleaning and disinfection procedure must always be preferred to manual cleaning because of the higher safety of the process. A good cleaning is also a requirement for successful sterilization.
 - We recommend washing the instruments before sterilization using automatic washers with thermal disinfection.
 - Place MEDIN instruments into the washing agent carefully to avoid moving of them freely during washing process or their overlapping which could damage their surface.
 - The MEDIN company does not prescribe the use of a specific detergent and disinfectant.
 - The washing agent and disinfectant have to be used in accordance with the instructions of the manufacturer. Keep prescribed concentration, temperature, water quality and exposure time.

19. Thermal disinfection

- Disinfect at a minimum temperature of 90 °C for at least 5 minutes in case of cleaning in washer.

20. Drying

- If drying cycle is not included in the manual washer:
 - dry each device thoroughly.

21. Maintenance	
21.1.	It is necessary to lubricate the instruments with movable parts, such as hinges or connections with springs or threads. Use solely the oil for healthcare suitable for sterilization for the instrument maintenance and lubrication.
22. Check	
22.1.	MEDIN instruments have to be checked after treatment and before sterilizing for the following: <ul style="list-style-type: none">Cleanliness. Any damage including corrosion pitting, change of colour, excessive and significant scratching, chipping, wear and cracks. Proper functioning including the sharpness of the cutting tools, flexion of the elastic devices, moving of the joint connections / threaded latches and other moving parts, such as handles, toothed latches and connections. Missing or removed (scraped) numbers of the parts and wear and tear. Discard the instruments with the aforementioned defects.
22.2.	Make sure that the instruments have an intact surface and that they are properly adjusted and functional. Do not use instruments that are severely damaged, have illegible marks, show signs of corrosion or have dull blades. Discard these instruments from their use. For more detailed instructions on the check of functionality, please contact your MEDIN sales representative. The manufacturer only is entitled to perform service interventions.
22.3.	A dismantled device should be assembled prior to sterilization unless otherwise indicated.
23. Packaging	
23.1.	Put the disinfected, cleaned, rinsed and dry instruments in their place in a sieve. In addition, use a suitable sterilization packaging or a reusable hard container such as a system of sterile barriers according to ISO 11607.
22.2.	It is necessary to pay attention to the protection of the instruments and ensure that the spiky or sharp instruments do not come into contact with other objects whose surface could be damaged.
23.3.	Single instruments, for which MEDIN storage is not available, pack into the sterilization bags that meet the requirements of ISO 11607.
24. Sterilization	
24.1.	Sterilize cleaned and disinfected medical devices only.
24.2.	It is recommended to use sterilisation with moist heat in a steam sterilizer at a temperature of 134 °C and the time of sterilization of 7 min. and pressure of 310 kPa for the medical devices.
25. Storage	
25.1.	Store the packed instruments in a dry clean environment without extreme temperature and humidity, away from direct sunlight, animals, insects, dust, fungi and chemicals.
26. Other information	
26.1.	The process of washing with thermal disinfection has been verified by the MEDIN company. In case of use of some other methods of the washing procedure than recommended, MEDIN company does not guarantee the result. When using some other method of washing than recommended, the MEDIN company recommends this process to be validated by the supplier of such washing facility.
26.2.	The MEDIN company has validated the method of sterilization with the moist heat at the temperature of 134 °C and the time of the sterilization of 7 min and pressure of 310 kPa for the representatives of the WORST-CASE MEDIN products.
26.3.	The instructions for cleaning and sterilization are stated in accordance with the norms and standards ISO 15883, ISO 17664, ISO 17665-1. The instructions listed above have been confirmed by the manufacturer of the medical device as appropriate for the preparation of MEDIN non-sterile medical devices. It is the responsibility of the provider to achieve the desired outcome by preparing a product with the use of monitoring of the procedure. Likewise every deviation of the provider from these recommendations should be properly assessed in terms of the efficiency and any possible undesired consequences.
27. Medical device disposal	
27.1.	A discarded instrument is considered as a hazardous waste. The user is responsible for the implementation of the measures to the safe handling and disposal of the product.
27.2.	The used implants are discarded after drying according to the Catalogue of waste – Group 18 – medical waste.

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