DE

Gebrauchsanweisung
PGA resoquick™ (PGA)
steriles, resorbierbares, synthetisches
Nahtmaterial

EN

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PGA resoquick™

EN

Instructions for use, please read carefully!

The information given in this package leaflet is updated regularly. Please read instructions carefully prior to use.

PGA

sterile, absorbable, synthetic
suture material

DESCRIPTION

PGA *resoquick*™ is a synthetic, sterile, absorbable surgical suture that consists of a polymer of glycolic acid. The braided thread is coated with resolactone (a mixture of calcium stearate and polycaprolactone). The coating is inert, nonantigenic and nonpyrogenic. A characteristic of PGA *resoquick* is its rapid loss of tensile strength.

PGA *resoquick* is undyed.

PGA *resoquick* is available in various thicknesses and lengths. PGA *resoquick* is also available as a needle-thread combination. PGA *resoquick* satisfies the requirements of the European Pharmacopoeia for sterile synthetic absorbable suture material.

INDICATIONS

PGA *resoquick* is indicated for coaptation of soft tissues when rapid absorption is required. The tensile strength profile of PGA *resoquick* renders it suitable for skin sutures, especially in paediatric surgery, for episiotomies, circumcisions, and in oral and maxillofacial surgery for subcutaneous and mucous-membrane sutures. PGA *resoquick* is also used effectively in ophthalmosurgery, e.g. for suturing the conjunctiva.

USE

The suture material should be chosen in accordance with the condition of the patient, the experience of the surgeon, the surgical technique, and the size of the wound.

ACTION

PGA *resoquick* causes a slight inflammatory tissue reaction initially, and is then slowly and progressively encapsulated by connective tissue. The gradual loss of tensile strength and the absorption of PGA *resoquick* suture material are due to hydrolytic breakdown to glycolic acid, which is then absorbed and metabolized in the body. Absorption leads initially to a loss of tensile strength and then to a reduction in mass.

In-vitro studies have shown that approximately 50% of initial tensile strength is still retained after seven days, and no tensile strength remains after 14 – 21 days.

CONTRAINDICATIONS

Due to its rapid loss of tensile strength, PGA *resoquick* should not be used for wound coaptations which will be maintained under traction for an extended period of time, nor if wound support will be necessary for

more than seven days. PGA *resoquick* is not intended for ligatures nor for use in the central circulatory system, nerve tissue or the heart.

WARNINGS / PRECAUTIONS / INTERACTIONS

As with all suture material, care must be taken to ensure that the thread is not damaged during handling. In particular, it must not be kinked or crushed by surgical instruments such as needle holders.

When tightening the suture always pull on the thread between the needle and the puncture channel. Do not pull the thread too strongly or over sharp objects. When stretching the thread avoid friction with the rubber glove, as this can damage the thread.

In order not to damage the needle, always grasp it 1/3 – 1/2 of the distance from the reinforced end to the point. Do not bend the needle, as this leads to a loss of stability. When using surgical needles the user should take particular care not to incur stab wounds. Used needles must be disposed of correctly (in order to avoid possible risks of infection).

The risk of wound dehiscence varies with the location of the wound and the suture material used, therefore the user should be familiar with the surgical techniques in which PGA *resoquick* is to be used. As with all foreign bodies, prolonged contact between PGA *resoquick* and saline solutions can lead to the formation of calculi (urinary and biliary tracts).

PGA *resoquick* is used where indicated in accordance with standard surgical suturing and knot-tying techniques and the experience of the user. Contaminated wounds should receive suitable surgical attention. Since PGA *resoquick* is an absorbable suture material, the surgeon may wish to consider the use of non-absorbable sutures to supplement sutures under tension or those subject to stretch or requiring additional support. Skin sutures which are required to remain in situ for more than seven days can cause local irritation. The visible part of the suture can therefore be cut short if necessary.

Caution is required when absorbable suture materials are used in poorly perfused tissue, as rejection of the thread and delayed absorption may occur. The suture material should be placed as deep as possible intracutaneously in order to minimize the erythema and induration that generally accompany the absorption process.

PGA *resoquick* is inappropriate in elderly, malnourished or debilitated patients or in patients with delayed wound healing.

ADVERSE REACTIONS

Inflammatory tissue reactions to the foreign body and transient local irritation may occur initially. In addition, erythema and tissue induration can develop during absorption of the intracutaneous sutures. The foreign-body effect can exacerbate pre-existing tissue inflammation initially.

STERILITY

PGA *resoquick* is gamma sterilised. Do not use opened or damaged packs! Do not resterilise!

STORAGE CONDITIONS

Store at not more than 25 °C protected from moisture and direct heat. Do not use after the expiry date!

DESCRIPTION OF THE SYMBOLS USED ON THE PACKAGE

REF = Reference number

LOT = Batch number

 = Use by year - month

 = Consult instructions for use

 = Do not reuse

 = Do not resterilize

 = Do not use if package is damaged

STERILE R = Sterilised using irradiation

 = Upper limit of temperature

 = undyed, braided, coated, absorbable

PGA = Polyglycolic acid

 = HIBC code

 = ligature pack:
sterile individual suture on a small roll

CE 0197 = CE marking and identification number of the notified body. Product conforms to the essential requirements of the Council Directive 93/42/EEC concerning medical devices.

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RESORBA Medical GmbH