DE

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GLYCOLON®

EN

Instructions for use, please read carefully!

The information given in this package leaflet is updated regularly. Please read instructions carefully prior to use.

PGA/PCL

sterile, absorbable, synthetic
suture material

DESCRIPTION

GLYCOLON® is a monofilament, synthetic, sterile, surgical suture material that is absorbable and consists of a copolymer of glycolic acid and ϵ -caprolactone. The characteristic feature of GLYCOLON® is the rapid loss of tear resistance.

GLYCOLON® is dyed or undyed.

GLYCOLON® is available in various gauge sizes and lengths and also as a needle and thread combination. GLYCOLON® complies with the requirements of the European Pharmacopoeia for sterile, absorbable, synthetic, monofilament suture material.

INDICATIONS

GLYCOLON® is suitable for soft tissue approximation and/or for ligatures, where the use of an absorbable suture is indicated.

USE

The suture material used should be selected according to the patient's condition, experience of the surgeon, the surgical technique and wound size.

ACTION

GLYCOLON® elicits a minimal inflammatory tissue reaction, and an ingrowth of connective tissue cells occurs. The gradual loss of tensile strength and absorption of GLYCOLON® are the consequences of hydrolytic processes. On absorption, there is first a reduction in tensile strength followed by a loss of mass.

In-vitro investigations showed that the dyed suture threads of GLYCOLON® still displayed about 55% of the original tensile strength after seven days, the undyed threads about 45%.

CONTRAINDICATIONS

Due to the rapid reduction in tensile strength, GLYCOLON® must not be used for long-lasting tissue approximation under stress or where wound support is necessary for more than seven days. GLYCOLON® is not intended for ligature or for use in the central circulatory system, neurological tissues or the heart. GLYCOLON® must not be used for abdominal closures or for approximation of fascia.

GLYCOLON® should not be used in cases of hypersensitivity to one of its constituents.

WARNINGS / PRECAUTIONS / INTERACTIONS

As with any suture material, care should be taken to avoid damage to the suture thread from handling, in particular it must not be crimped or crushed by surgical instruments such as needle holders.

When tightening the suture, always pull on the thread between the needle and the puncture channel. Do not pull the thread too vigorously or over sharp objects. When stretching the suture thread, avoid surface friction with the rubber glove because this can damage the thread.

To avoid damaging the needle, it should always be grasped 1/3 – 1/2 of the distance from the attachment end to the needle point. Reshaping the needle should be avoided because this may cause a loss of stability. Users of surgical needles must always take care to avoid needle stick injuries. Used needles must be disposed of properly (in order to avoid possible risks of infection).

The risk of wound dehiscence varies according to the site of application and suture material used and therefore the user should be familiar with the surgical procedures in which GLYCOLON® is used and take account of the in vitro properties.

GLYCOLON® may be inappropriate in elderly, malnourished or debilitated patients or patients whose disease may delay wound healing. As with any foreign body, on contact with salt solutions (such as those in the urinary or biliary tracts) GLYCOLON® can result in calculus formation.

GLYCOLON® is used when indicated according to the accepted surgical suturing methods and knotting techniques and according to the experience of the user.

Contaminated wounds should be cleaned by appropriate surgical methods. Since GLYCOLON® is an absorbable suture material, the surgeon should use supplemental nonabsorbable suture material for sutures that are under tension or are stretched or require additional support.

Skin sutures that have to remain in place for longer than seven days may cause localised skin irritation and the external part of the suture should therefore be snipped off if necessary.

Caution should be exercised in the use of absorbable suture material in poorly perfused tissues, because here suture extrusion and delayed absorption may occur. Subcuticular sutures should be placed as deeply as possible to minimise erythema and induration normally associated with the absorption process.

ADVERSE REACTIONS

In the beginning, inflammatory tissue reactions to the applied foreign body and transient local irritation can occur. In addition, erythema and tissue induration may develop during the absorption of subcuticular sutures. The effect of a foreign body may initially enhance existing tissue inflammation.

STERILITY

GLYCOLON® is sterilised by ethylene oxide. Already opened or damaged packs must not be used! Do not resterilise!

STORAGE CONDITIONS

Store below 25 °C, protected from moisture and direct heat. Do not use after the expiry date!

DESCRIPTION OF THE SYMBOLS USED ON THE PACKAGE

REF = Reference number

LOT = Batch number

 = Use by year - month

 = Consult instructions for use

 = Do not reuse

 = Do not resterilize

 = Do not use if package is damaged

STERILE EO = Sterilised using ethylene oxide

 = Upper limit of temperature

 = undyed, monofilament, absorbable

 = dyed, monofilament, absorbable

PGA/PCL = Poly(glycolide-co-caprolactone)

 = HIBC code

CE 0197 = CE marking and identification number of the notified body. Product conforms to the essential requirements of the Council Directive 93/42/EEC concerning medical devices.

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RESORBA Medical GmbH