PARASORB® Cone

- DE** Gebrauchsanweisung:
Bitte aufmerksam lesen!
- EN** Instructions for use:
Please read carefully!
- FR** Instructions d'utilisation:
À lire attentivement!
- ES** Instrucciones de utilización:
¡Rogamos leer atentamente!
- IT** Istruzioni per l'uso:
Leggere attentamente!
- EL** Οδηγίες χρήσεως:
να διαβάσετε με προσοχή!
- LT** Naudojimo instrukcija:
prašome įdėmiai perskaityti!
- RU** Инструкция по применению:
пожалуйста, прочитайте внимательно!
- RO** Instrucțiuni de utilizare:
Vă rugăm citiți cu atenție!
- NO** Bruksanvisning:
Les nøye!
- DA** Brugsanvisningen:
Skal læses grundigt!
- FI** Käyttöohjeisiin:
Lue huolellisesti.
- SV** Bruksanvisningen:
Läs noggrant.
- NL** Gebruiksaanwijzing:
Gelieve aandachtig te lezen!
- PT** Instruções de utilização:
Por favor, leia atentamente!
- CS** Návod k použití:
Čtěte pozorně!
- PL** Instrukcjach używania:
Należy uważnie przeczytać!

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 HIBC-Code

EN Instructions for use – Please read carefully!

The information given in this package leaflet is updated regularly. Please read instructions carefully prior to use.

Absorbable Collagen Cone
- haemostyptic, sterile -

QUALITATIVE AND QUANTITATIVE COMPOSITION

1 PARASORB® Cone contains 22.4 mg of native collagen fibrils of equine origin.

PHARMACEUTICAL FORM

Absorbable, local haemostyptic agent.

THERAPEUTIC INDICATIONS

PARASORB® Cone is used as a haemostatic agent after tooth extractions and oral surgery, respectively. Combination with fibrin glue is recommended in patients with an increased risk of bleedings. The corresponding instructions must be strictly observed. PARASORB® Cone is also used to minimize alveolar ridge atrophy after tooth extractions.

POSOLGY AND METHOD OF ADMINISTRATION

Haemostasis occurs on a physical basis. The contact of collagen with blood leads to aggregation of platelets which precipitate in large numbers on the collagen matrix, disintegrate and release clotting factors which, together with plasma factors, facilitate the formation of fibrin. The contact of collagen with blood leads to aggregation of platelets which precipitate in large numbers on the collagen matrix, disintegrate and release clotting factors which, together with plasma factors, facilitate the formation of fibrin. The collagen matrix gives additional reinforcement to the blood clot thus preventing soft tissue to grow into the extraction alveola. Consequently alveolar ridge atrophy is depleted.

Due to its structure PARASORB® Cone is capable of absorbing large quantities of fluid. During this purely mechanical process of absorption of secretion, rejected materials such as bacteria and fibrin clots are also removed.

Medical treatment of severe bleeding disorders is mandatory prior to the beginning of the procedure.

Unless otherwise prescribed, PARASORB® Cone is applied directly into the extraction alveola. PARASORB® Cone moulds itself into the site of extraction or excision wound respectively, stops bleeding and will be absorbed. The wounds can be sutured.

CONTRAINDICATIONS

PARASORB® Cone must not be used in patients with a known hypersensitivity to products of equine origin. The application in infected areas is an absolute contraindication. No data on the safe use of PARASORB® Cone during pregnancy and lactation are available. Therefore PARASORB® Cone should only be used with strict medical indications during pregnancy and lactation.

UNDESIRABLE EFFECTS

Literature reports occasional events suggesting that potentially present infections may be aggravated. Despite its proven low antigenicity, hypersensitivity reactions may occur due to the xenogenous origin of the collagen.

INTERACTION

Chlorine-releasing antiseptic agents, tannins and caustics should not be used in combination with collagen. Due to the specific impact of PARASORB® Cone on thrombocyte aggregation and plasma coagulation, a reduction of its efficacy has to be considered through substances inhibiting aggregation and coagulation.

EFFECTS ON ABILITY TO DRIVE AND USE MACHINES

To date none known.

PLEASE NOTE

Because of the affinity of collagen to bleeding surfaces, dry instruments and gloves should be used when applying PARASORB® Cone.

SPECIAL PRECAUTIONS FOR STORAGE

Do not store above 25°C.
Protect from moisture.
Keep out of the reach of children.

SHELF LIFE

The contents of unused, yet opened or damaged, packets must not be resterilised and should therefore be discarded. Expiry date: See imprint on outer package. PARASORB® Cone must not be used after the expiry date.

PACK SIZES

PARASORB® Cone
10 pieces
Ø: 1.2 cm; H: 1.6 cm
REF DK1010

DESCRIPTION OF THE SYMBOLS USED ON THE PACKAGE

 Reference number

 Batch number

 Use by year - month

 Consult instructions for use

 Do not reuse

 Do not resterilize

 Do not use if package is damaged

 Sterilised using ethylene oxide

 25°C
Upper limit of temperature

 Protect against humidity

 0197

CE marking and identification number of the notified body. Product conforms to the essential requirements of the Council Directive 93/42/EEC concerning medical devices.

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COMPOSITION QUALITATIVE ET QUANTITATIVE

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