DE

Gebrauchsanweisung
SUPOLENE (PET)
steriles, nicht-resorbierbares
Nahtmaterial

EN

Instructions for use
SUPOLENE (PET)
sterile, nonabsorbable
suture material

FR

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SUPOLENE

EN

Instructions for use, please read carefully!

The information given in this package leaflet is updated regularly. Please read instructions carefully prior to use.

POLYESTER

Polytetrafluoroethylene-coated polyester,
sterile, nonabsorbable **suture material**

DESCRIPTION

SUPOLENE is a braided and coated nonabsorbable surgical suture material of polyethylene terephthalate with the molecular formula $(C_{10}H_8O_4)_n$. The material has a coating of polytetrafluoroethylene that acts as a lubricant. This reduces the capillarity of the thread to a minimum, improves tissue passage, and facilitates knotting.

SUPOLENE is available both undyed (white) and dyed green. SUPOLENE is available in various thicknesses and lengths and also as a needle-thread combination. SUPOLENE satisfies the requirements of the European Pharmacopoeia for sterile suture material of polyethylene terephthalate.

INDICATIONS

Because of its hydrophobic properties and long-term stability, SUPOLENE is suitable for the attachment of implants and prostheses. It is also intended for use in coaptation of soft tissues and/or for ligatures, in particular in neurosurgery and ophthalmic and cardiovascular surgery.

USE

The suture material should be chosen in accordance with the condition of the patient, the experience of the surgeon, the surgical technique, and the size of the wound.

ACTION

SUPOLENE causes only a slight inflammatory tissue reaction initially and is then slowly and progressively encapsulated by connective tissue.

CONTRAINDICATIONS

None known.

WARNINGS / PRECAUTIONS / INTERACTIONS

As with all suture material, care must be taken to ensure that the thread is not damaged during handling. In particular, it must not be kinked or crushed by surgical instruments such as needle holders.

When tightening the suture always pull on the thread between the needle and the puncture channel.

In order not to damage the needle, always grasp it 1/3 – 1/2 of the distance from the end to the point. Do not bend the needle, as this leads to a loss of stability. Because of the risk of infection, the user should take particular care not to incur stab wounds when using surgical needles. Used needles must be disposed of correctly (in order to avoid possible risks of infection).

The risk of wound dehiscence varies with the location of the wound and the suture material used, therefore the user should be familiar with the surgical techniques in which SUPOLENE is to be used. As with all foreign bodies, prolonged contact between SUPOLENE and saline solutions can lead to the formation of calculi (urinary and biliary tracts).

SUPOLENE is used where indicated in accordance with standard surgical suturing and knotting techniques and the experience of the user.

ADVERSE REACTIONS

Inflammatory tissue reactions and transient irritation in the area of the wound may occur initially. The foreign-body effect can exacerbate pre-existing tissue inflammation initially.

STERILITY

SUPOLENE is sterilised with ethylene oxide. Do not use opened or damaged packs! Do not resterilise!

STORAGE CONDITIONS

Store at not more than 25 °C protected from moisture and direct heat. Do not use after the expiry date!

DESCRIPTION OF THE SYMBOLS USED ON THE PACKAGE

REF = Reference number

LOT = Batch number

 = Use by year - month

 = Consult instructions for use

 = Do not reuse

 = Do not resterilize

 = Do not use if package is damaged

STERILE EO = Sterilised using ethylene oxide

 = Upper limit of temperature

 = undyed, braided, coated, nonabsorbable


 = dyed, braided, coated, nonabsorbable

PET = Polyester

 = HIBC code

 = control release needle  = loop suture

 = ligature pack:
sterile individual suture on a small roll

 = pledget

CE 0197 = CE marking and identification number of the notified body. Product conforms to the essential requirements of the Council Directive 93/42/EEC concerning medical devices.

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RESORBA Medical GmbH