**DE**

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**PGA RESORBA®**

**EN**

**Instructions for use, please read carefully!**

**The information given in this package leaflet is updated regularly. Please read instructions carefully prior to use.**

## **POLYGLYCOLIC ACID**

sterile, absorbable, synthetic  
suture material

### **DESCRIPTION**

PGA RESORBA® is a synthetic sterile absorbable surgical suture material that consists of a polymer of glycolic acid. The braided thread is coated with resolactone (a mixture of calcium stearate and polycaprolactone). The coating is inert, nonantigenic, and nonpyrogenic. For better visualization PGA RESORBA® is dyed by addition of D & C Violet no. 2, but is also available undyed. PGA RESORBA® is available in various thicknesses and lengths and also as a needle-thread combination. PGA RESORBA® satisfies the requirements of the European Pharmacopoeia for sterile synthetic absorbable suture material.

### **INDICATIONS**

PGA RESORBA® is intended for coaptation of soft tissues and/or for ligatures, including use in ophthalmology. Insufficient data are available as yet on use of PGA RESORBA® in the heart and central vascular system.

### **USE**

The suture material should be chosen in accordance with the condition of the patient, the experience of the surgeon, the surgical technique, and the size of the wound.

### **ACTION**

PGA RESORBA® causes a slight inflammatory tissue reaction initially and is then slowly and progressively encapsulated by connective tissue. The gradual loss of tensile strength and the absorption of PGA RESORBA® suture material are due to hydrolytic breakdown to glycolic acid, which is then absorbed and metabolized in the body. Absorption leads initially to a loss of tensile strength and then to a reduction in mass. The standard in-vitro test revealed the following **absorption parameters**:

Time of measurement	Remaining tensile strength (% of original tensile strength)
7 days	80 %
14 days	50 %
21 days	25 %

### **CONTRAINDICATIONS**

As PGA RESORBA® is absorbable, it should not be used for wound coaptations in which the wound edges are subject to considerable tension.

### **WARNINGS / PRECAUTIONS / INTERACTIONS**

As with all suture material, care must be taken to ensure that the thread is not damaged during handling. In particular, it must not be kinked or crushed by surgical instruments such as needle holders.

When tightening the suture always pull on the thread between the needle and the puncture channel. Do not pull the thread too strongly or over sharp objects. When stretching the thread avoid friction with the rubber glove, as this can damage the thread.

In order not to damage the needle, always grasp it 1/3 – 1/2 of the distance from the reinforced end to the point. Do not bend the needle, as this leads to a loss of stability. Because of the risk of infection, the user should take particular care not to incur stab wounds when using surgical needles. Used needles must be disposed of correctly (in order to avoid possible risks of infection).

The risk of wound dehiscence varies with the location of the wound and the suture material used, therefore the user should be familiar with the surgical techniques in which PGA *RESORBA*® is to be used. As with all foreign bodies, prolonged contact between PGA *RESORBA*® and saline solutions can lead to the formation of calculi (urinary and biliary tracts).

PGA *RESORBA*® is used where indicated in accordance with standard surgical suturing and knotting techniques and the experience of the user. Contaminated wounds should receive suitable surgical attention. As PGA *RESORBA*® is an absorbable suture material, the surgeon should if necessary also use nonabsorbable suture material for sutures that are under tension, are stretched, or need additional support.

Skin sutures that must remain in situ for more than seven days can cause local irritation. The external part of the suture should therefore be cut off if necessary.

In many cases, especially in orthopaedic operations, the surgeon has the option of resorting to external immobilization of joints. Caution is required when absorbable suture materials are used in poorly perfused tissue, as rejection of the thread and delayed absorption may occur. The suture material should be placed as deep as possible intracutaneously in order to minimize the erythema and induration that generally accompany the absorption process.

PGA *RESORBA*® may be inappropriate in elderly, malnourished or debilitated patients, or in patients with delayed wound healing.

### **ADVERSE REACTIONS**

Inflammatory tissue reactions to the foreign body and transient local irritation may occur initially. In addition, erythema and tissue induration can develop during absorption of the intracutaneous sutures. The foreign-body effect can exacerbate pre-existing tissue inflammation initially.

### **STERILITY**

PGA *RESORBA*® is sterilised by ethylene oxide. Already opened or damaged packs must not be used! Do not resterilise!

### **STORAGE CONDITIONS**

Store at not more than 25 °C protected from moisture and direct heat. Do not use after the expiry date!

## DESCRIPTION OF THE SYMBOLS USED ON THE PACKAGE

**REF** = Reference number

**LOT** = Batch number

 = Use by year - month

 = Consult instructions for use

 = Do not reuse

 = Do not re-sterilize

 = Do not use if package is damaged

**STERILE EO** = Sterilised using ethylene oxide

 = Upper limit of temperature

 = undyed, braided, coated, absorbable

 = dyed, braided, coated, absorbable

**PGA** = Polyglycolic acid

 = HIBC code

 = control release needle

 = loop

 = loop suture

 = ligature pack:  
sterile individual suture on a small roll

**CE** 0197 = CE marking and identification number of the notified body. Product conforms to the essential requirements of the Council Directive 93/42/EEC concerning medical devices.

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**RESORBA Medical GmbH**