CHIRMAX

VITREX

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CHIRLAC BRAIDED

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Instructions for Use

CHIRLAC BRAIDED

Sterile braided synthetic absorbable surgical suture composed of a polymer of glycolic acid. This suture is coated.

GB

Product description

Chirlac braided is synthetic absorbable braided surgical suture composed of a polymer of glycolic acid. The substances contained in the coating and suture are nonantigenic and noncollagenous. Sterile Chirlac braided is available in dyed violet or undyed natural. The product elicits a minimal acute inflammatory reaction in tissues. Progressive loss of tensile strength and absorption of Chirlac braided occurs by hydrolysis, where the polymer degrades to glycolic acid, which is subsequently absorbed and metabolized by the body. Absorption begins as a loss of tensile strength without appreciable loss of mass. Implantation studies in animals indicate that sterile Chirlac braided retains approximately 75 % of its original tensile strength at two weeks post implantation, with approximately 50 % remaining at three weeks. The absorption of sterile Chirlac braided is completed within 80-90 days.

Indication

Sterile Chirlac is designed for use in general soft tissue approximation and/or ligation, for the surgical adaptation of muscles and fasciae, for suturing of parenchymal organs, peripheral nerves, for ligatures of small vessels and for intradermal suturing.

Reaction in tissue

Slight initial inflammatory reaction in tissue.

Interaction

Sterile Chirlac should not be used in the tissues which may undergo expansion, stretching or distension or which may require a long-term mechanical support.

Warnings

This suture may be inappropriate in elderly, malnourished or debilitated patients, or in patients whose bad conditions may delay wound healing. The use of supplemental non-

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absorbable sutures should be considered by the surgeon. Skin sutures which must remain in place longer than 7 days may cause local irritation and so the part of the suture which is showing should be snipped off or removed as indicated. Prolonged contact of this suture with salt solution, such as those found in urinary or biliary tracts, may result in calculus formation. It is not possible to use the product repeatedly; it will not be sterile. **Handling instruction**
Product is intended for single usage!
Do not use if package is opened or damaged!
Check the expiry date printed on the packing before use!
Do not use the product after the expiry date!
Do not resterilize!
Only trained staff (doctors, nurses) may apply surgical sutures.

How to open the surgical suture packaging: take the bag out of the box and take it in both hands, open the bag by pulling the two loose ends of the package so that the sterility of the surgical suture is not affected. Do not bend the bag with the suture. Sutures in atraumatic needles: after having separated the upper part of the paper spool, grasp the needle with the needle holder in an area one-third (1/3) to one half (1/2) of the distance from the attachment end to the point and pull it out of the bag. Do not break the sterility of the suture during unpacking, when removing the needle from the case, always use a needle holder. The suture is sterile. Use the suture taken out of the bag immediately. Put the remaining material into a separate waste container. This will prevent injuries of the medical staff. Upon correct handling the producer guarantees absolute quality of the product. **Sterilization**
Chirlac is sterilized by ethylene oxide. **Storage**
The sutures have to be stored in their original packaging in a clean, dry room at the temperature of 10° C to 25° C, protected from direct light and all kinds of moisture.

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