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Instructions for Use

CHIRLAC RAPID BRAIDED

Sterile braided synthetic absorbable surgical suture composed of a polymer of glycolic acid. This suture is coated.

GB

Product description

Chirlac rapid is synthetic absorbable braided surgical suture composed of a polymer of glycolic acid. The substances contained in the coating and suture are nonantigenic and noncollagenous. Sterile Chirlac rapid is available dyed or undyed. Product elicits a minimal acute inflammatory reaction in tissues. Progressive loss of tensile strength and absorption of Chirlac rapid occurs by hydrolysis, where the polymer degrades to glycolic acid, which is subsequently absorbed and metabolized by the body. Absorption begins as a loss of tensile strength without appreciable loss of mass. Five days post implantation approximately 50 % of its original tensile strength remains. All of the original tensile strength is lost by approximately 10 to 14 days post implantation. Complete absorption at 42 days.

Indication

Sterile Chirlac rapid is designed for the use in general soft tissue approximation where only short term wound support is required and where the rapid absorption of the suture would be beneficial. Chirlac rapid is indicated for closure of skin, particularly in pediatric surgery, gynecology (e.g. episiotomies), ophthalmic surgery (e.g. conjunctiva suturing) and for circumcison of oral mucosa.

Reaction in tissue

Slight initial inflammatory reaction in tissue.

Interaction

Sterile Chirlac rapid should not be used in the tissues which may undergo expansion, stretching or distension or which may require a long-term mechanical support longer than seven days. Chirlac rapid suture is not for use in ligation, cardiovascular and neurological tissues.

Warnings

This suture may be inappropriate in elderly, malnourished or debilitated patients, or in patients whose bad conditions may delay wound healing. The use of supplemental non-absorba-

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tle sutures should be considered by the surgeon. Skin sutures, which must remain in place longer than 7 days may cause local irritation and so the appearing part of the suture should be snipped off or removed as indicated. Prolonged contact of this suture with salt solution, such as those found in urinary or biliary tracts, may result in calculus formation. It is not possible to use the product repeatedly; it will not be sterile.

Handling instruction

Product is intended for single usage!

Do not use if package is opened or damaged!

Check the expiry date printed on the packing before use!

Do not use the product after the expiry date!

Do not resterilize!

Only trained staff (physicians, nurses) may apply surgical sutures. **How to open the surgical suture packaging:** take the bag out of the box and take it in both hands, open the bag by pulling the two loose ends of the package so that the sterility of the surgical suture is not affected. Do not bend the bag with the suture.

Sutures in atraumatic needles: after having separated the upper part of the paper spool, grasp the needle with the needle holder in an area one-third (1/3) to one half (1/2) of the distance from the attachment end to the point and pull it out of the bag. Do not break the sterility of the suture during unpacking, when removing the needle from the case, always use a needle holder. The suture is sterile. Use the suture taken out of the bag immediately. Put the remaining material into a separate waste container. This will prevent injuries of the medical staff. Upon correct handling the producer guarantees absolute quality of the product.

Sutures without needles: after having separated the upper part of the paper spool, a paper fixture appears which allows you to take the sutures out of the bag easily. The suture is sterile. Use the filament taken out of the bag immediately. Put the remaining material into a separate waste container. This will prevent injuries to the medical staff. If used properly the producer guarantees absolute quality of the product.

Sterilization

Product is sterilized by ethylene oxide or by irradiation – Gamma rays..

Storage

The sutures have to be stored in their original packaging in a clean, dry room at the temperature of 10° C to 25° C, protected from direct light and all kinds of moisture.

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