IPS Empress® Direct is a light-curing, radiopaque, nano-hybrid composite based on state-of-the-art technology, which is suitable for the fabrication of restorations that fulfill the highest demands. IPS Empress Direct cures with light in the wavelength range of 400-500 nm (blue light).

IPS Empress Direct is offered in 32 different shades and 5 different degrees of translucency. The following shades are available:

A1, A2, A3, A3.5, A4, B1, B2, B3, B4, C1, C2, C3, D2, D3, Bleach-L, Bleach-XL

A1, A2, A3, A3.5, A4, IVA5, IVA6, B1, B2, C3, D2, Bleach-L, Bleach-XL

Trans 20, Trans 30, Opal*

*Not radiopaque due to its high transparency

The monomer matrix is composed of dimethacrylates (20-21.5 wt%, opalescent shade 17 wt%). The fillers contain barium glass, ytterbium trifluoride, mixed oxide, silicon dioxide and copolymer (77.5-79 wt%, opalescent shade 83 wt%). Additional contents: additives, catalysts, stabilizers and pigments (<1.0 wt%). The total content of inorganic fillers is 75-79 wt% or 52-59 vol% (opalescent shade 60.5 wt% or 45 vol%). The particle size of the inorganic fillers is between 40 nm and 3000 nm with a mean particle size of 550 nm.

- Anterior restorations (Class III, IV)
- Restorations in the posterior region (Classes I and II)
- Tooth neck restorations (Class V, e.g. in the case of cervical caries, root erosion, wedge-shaped defects)
- Corrections of tooth positions and tooth shape (e.g. diastema closure, closure of interdental "black triangles", lengthening of the incisal edge)
- Direct veneers

Placement of IPS Empress Direct restorations is contraindicated

- if a dry working field cannot be established, or if the stipulated technique cannot be applied
- if a patient is known to be allergic to any of the ingredients in IPS Empress Direct

In rare cases, components of IPS Empress Direct may lead to sensitization. IPS Empress Direct should not be used in such cases. To avoid possible irritation of the pulp, areas close to the pulp should be protected with a suitable pulp/dentin protector (selectively apply a calcium hydroxide-based preparation in areas close to the pulp and coat with a suitable liner).

Substances containing eugenol/oil of cloves may inhibit the polymerization of composite materials. Consequently, the application of such materials together with IPS Empress Direct must be avoided. Discoloration may occur in combination with cationic mouthwashes, plaque disclosing agents and chlorhexidine.

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— Light-curing resin-based dental restorative material

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Application

1. Shade determination

Clean the teeth prior to shade determination. The enamel shade is selected with the tooth still moist. The dentin shade can be taken either by placing the tab next to the prepared tooth or next to the cervical area of the canine.

2. Isolation

Appropriate isolation is required, best with a rubber dam (e.g. OptraDam® Plus).

3. Cavity preparation

Cavity preparation is carried out according to the requirements of the adhesive technique, i.e. protecting the tooth structure. Enamel edges in the anterior region are bevelled. In the posterior region, only sharp enamel edges should be rounded (finishing diamonds, 25-40 µm). Caries-free cervical defects are not prepared, only cleaned with pumice or other suitable cleaning pastes with the help of rubber cups or rotary brushes. In case of sclerotic dentin, superficial reconditioning with a rose bur can improve the adhesive bond. Subsequently, all residue in the cavity is removed with water spray. The cavity is dried with water- and oil-free air.

4. Pulp protection / Base

Do not apply a base material when using an enamel/dentin bonding agent. Only cover very deep areas close to the pulp with a calcium hydroxide material (e.g. ApexCal®) and subsequently use pressure-resistant cement (e.g. glass ionomer cement, such as Vivaglass® Liner). Do not cover other cavity walls, since they can be used to support the bond with an enamel/dentin adhesive.

5. Placement of a matrix / interdental wedge

Use a wrap-around matrix for cavities affecting the proximal area or a sectional matrix and wedge.

6. Conditioning / Application of the bonding agent

Condition the tooth structure and apply the bonding agent according to the Instructions for Use of the product in use. We recommend using Syntac® or ExciTE® F (with phosphoric acid etching) or AdheSE® or AdheSE® One F self-etching adhesive.

7. Application of IPS Empress Direct

- Since a large number of IPS Empress Direct Dentin and Enamel shades in different translucencies for the creation of anatomical detail are available, any layering technique can be creatively used to achieve highly esthetic restorations. For an optimum result, apply IPS Empress Direct in layers of max. 2 mm (or 1.5 mm Dentin shades) thickness and adapt with e.g. OptraSculpt®/ OptraSculpt® Pad.
- Sufficient exposure to the curing light prevents incomplete polymerization. Recommended curing times (exposure time) and light intensities are listed in Table 1. Hold the light emission window as closely as possible to the surface of the restorative material.
- When using a metal matrix, additionally polymerize the restoration from the buccal or lingual/palatal aspect after removing the matrix, if no Bluephase® curing light was used or if the light probe could not be optimally positioned (e.g. far away from the material or with divergent radiation).
- In many cases, a flowable composite is now used for the initial layer in order to create an even cavity floor and to simplify the adaptation of the subsequently used restorative material. The use of a flowable composite (e.g. Tetric EvoFlow®) as an initial thin layer is optional. This layer must be light-cured in a separate step (the instructions of the respective Instructions for Use must be observed).

8. Finishing / Checking the occlusion / Polishing

Remove excess material with suitable finishers (e.g. Astropol® F) or fine diamonds after polymerization. Remove proximal excess with diamond finishers, tungsten carbide finishers, finishing strips, or flexible finishing discs. Check the occlusion and articulation and apply appropriate grinding corrections to prevent premature contacts or undesired articulation paths on the surface of the restorations. Use silicone polishers (e.g. Astropol® P, Astropol® HP, Astrobrush or OptraPol® Next Generation) as well as polishing discs and polishing strips to polish the restorations to a durable high gloss.

Additional information

1. If necessary, additional IPS Empress Direct material can be directly applied to polymerized material. If the IPS Empress Direct restoration has already been polished, it must first be roughened and wetted with Heliobond® before fresh IPS Empress Direct is applied.
2. IPS Empress Direct should be at room temperature when it is applied. Cool temperatures render the material difficult to extrude.
3.  For single use only. If IPS Empress Direct is applied from the Cavifil directly in the mouth of the patient, the Cavifil must be discarded after use to prevent cross-contamination between patients.
4. Syringes or Cavifils should not be disinfected with oxidizing disinfection agents.

Warning

Unpolymerized IPS Empress Direct should not come in contact with the skin, mucous membrane, or eyes. Unpolymerized IPS Empress Direct may have a slight irritating effect and may lead to a sensitization against methacrylates. Commercial medical gloves do not provide protection against sensitivity to methacrylates.

Storage

- Storage temperature 2–28 °C / 36–82 °F.
- Close syringes/Cavifils immediately after use. Exposure to light causes premature polymerization.
- Do not use IPS Empress Direct after the indicated date of expiration.
- Shelf life: see information on Cavifils, syringes and packages.

Keep material out of children's reach.

For use in dentistry only.

The material has been developed solely for use in dentistry. Processing should be carried out strictly according to the Instructions for Use. Liability cannot be accepted for damages resulting from failure to observe the Instructions or the stipulated area of application. The user is responsible for testing the products for their suitability and use for any purpose not explicitly stated in the Instructions. Descriptions and data constitute no warranty of attributes and are not binding.

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Table 1

Unit	Bluephase® C8	Bluephase® Style	Bluephase®	Bluephase® 20i
Program				
Turbo	–	–	–	5 s
High Power	15 s	10 s	10 s	10 s
Soft Start	20 s	–	15 s	15 s

Light Intensity	Exposure Time
≥ 500 mW/cm²	20 s
≥ 1'000 mW/cm²	10 s

Ivoclar Vivadent AG

Bendererstrasse 2, 9494 Schaan, Liechtenstein,

Tel. +423 235 35 35, Fax +423 235 33 60

www.ivoclarvivadent.com

Ivoclar Vivadent Pty. Ltd.

1 – 5 Overseas Drive, P.O. Box 367, Noble Park, Vic. 3174, Australia

Tel. +61 3 9795 9599, Fax +61 3 9795 9645

www.ivoclarvivadent.com.au

Ivoclar Vivadent Ltda.

Alameda Caiapós, 723, Centro Empresarial Tamboré, CEP 06460-110 Barueri – SP, Brazil

Tel. +55 11 2424 7400, Fax +55 11 3466 0840

www.ivoclarvivadent.com.br

Ivoclar Vivadent Inc.

1-6600 Dixie Road, Mississauga, Ontario, L5T 2Y2, Canada

Tel. +1 905 670 8499, Fax +1 905 670 3102

www.ivoclarvivadent.us

Ivoclar Vivadent Shanghai

Trading Co., Ltd., 2/F Building 1, 881 Wuding Road, Jing An District

200040 Shanghai, China

Tel. +86 21 6032 1657, Fax +86 21 6176 0968

www.ivoclarvivadent.com

Ivoclar Vivadent Marketing Ltd.

Calle 134 No. 7-B-83, Of. 520, Bogotá, Colombia

Tel. +57 1 627 3399, Fax +57 1 633 1663

www.ivoclarvivadent.co

Ivoclar Vivadent SAS

B.P. 118, F-74410 Saint-Jorioz, France

Tel. +33 4 50 88 64 00, Fax +33 4 50 68 91 52

www.ivoclarvivadent.fr

Ivoclar Vivadent GmbH

Dr. Adolf-Schneider-Str. 2, D-73479 Ellwangen, Jagst, Germany

Tel. +49 7961 889 0, Fax +49 7961 6326

www.ivoclarvivadent.de

Wieland Dental + Technik GmbH & Co. KG

Schwenninger Strasse 13, D-75179 Pforzheim, Germany

Tel. +49 7231 3705 0, Fax +49 7231 3579 59

www.wieland-dental.com

Ivoclar Vivadent Marketing (India) Pvt. Ltd.

503/504 Raheja Plaza, 15 B Shah Industrial Estate, Veera Desai Road, Andheri (West)

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www.ivoclarvivadent.com

Ivoclar Vivadent Ltd.

12 Omega St, Rosedale, PO Box 303011 North Harbour, Auckland 0751, New Zealand

Tel. +64 9 914 9999, Fax +64 9 914 9990

www.ivoclarvivadent.co.nz

Ivoclar Vivadent Polska Sp. z o.o.

Al. Jana Pawla II 78, 00-175 Warszawa, Poland

Tel. +48 22 635 5496, Fax +48 22 635 5469

www.ivoclarvivadent.pl

Ivoclar Vivadent Marketing Ltd.

Prospekt Andropova 18 korp. 6/ , office 10-06, 115432 Moscow, Russia

Tel. +7 499 418 0300, Fax +7 499 418 0310

www.ivoclarvivadent.ru

Ivoclar Vivadent Marketing Ltd.

Qlaya Main St., Siricon Building No.14, 2nd Floor, Office No. 204

P.O. Box 300146, Riyadh 11372, Saudi Arabia

Tel. +966 11 293 8345, Fax +966 11 293 8344

www.ivoclarvivadent.com

Ivoclar Vivadent Pte. Ltd.

171 Chin Swee Road, #02-01 San Centre, Singapore 169877

Tel. +65 6535 6775, Fax +65 6535 4991

www.ivoclarvivadent.com

Ivoclar Vivadent S.L.U.

C/ Ribera del Loira nº 46, 5ª planta, 28042 Madrid, Spain

Tel. + 34 913 757 820, Fax + 34 913 757 838

www.ivoclarvivadent.es

Ivoclar Vivadent AB

Dalvägen 14, S-169 56 Solna, Sweden

Tel. +46 8 514 939 30, Fax +46 8 514 939 40

www.ivoclarvivadent.se

Ivoclar Vivadent Liaison Office

: Tesvikiye Mahallesi, Sakayik Sokak, Nisantas' Plaza No:38/2, Kat:5 Daire:24

34021 Sisli – Istanbul, Turkey

Tel. +90 212 343 0802, Fax +90 212 343 0842

www.ivoclarvivadent.com

Ivoclar Vivadent Limited

Ground Floor Compass Building, Feldspar Close, Warrens Business Park, Enderby

Leicester LE19 4SE, United Kingdom

Tel. +44 116 284 7880, Fax +44 116 284 7881

www.ivoclarvivadent.co.uk

Ivoclar Vivadent, Inc.

175 Pineview Drive, Amherst, N.Y. 14228, USA

Tel. +1 800 533 6825, Fax +1 716 691 2285

www.ivoclarvivadent.us