BRILLIANT FLOW UNIVERSAL COMPOSITE

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Note <ul style="list-style-type: none">Light can lead to premature curing. Recap syringe after each use. To prevent run off, pull back the syringe plunger. For repairs, the treated surface must be roughened first.	
General information	
Preliminary preparation <p>Cleansing the tooth</p> <p>Cleanse tooth to be treated and the adjacent teeth using a brush and prophy paste which does not contain fluoride.</p>	
Selection of shade <p>Selection of shade should be made before isolation. The shade is ideally determined in broad daylight using the VitaTM Shade Guide. Shades reach their final tone 24 hours after curing.</p>	
Isolation of the cavity (rubber dam) <p>Adequate isolation is absolutely necessary for achieving optimum results. The use of HYGENIC®/Roeko Dental Dams are recommended.</p>	
Matrix and inter-dental wedges <p>When working approximally, apply a thin matrix band. Fix the matrix proximally using interdental wedges.</p>	
Adhesive system <p>Apply adhesive system (e.g. Coltène® One Coat Bond SL oder A.R.T. BOND) according to the corresponding instructions for use.</p>	
Working time <p>Brilliant Flow is light sensitive and should not be left too long under intensive light before curing, especially operatory light or sun light.</p>	
Polymerization <p>Brilliant Flow can be cured using all current curing lights (e.g. Coltolux® LED).</p>	
Recommended exposure times:	
Layer Thickness: 2 mm	> 500 mW/ cm²
Brilliant Flow A1 /B1	20 s
Brilliant Flow A2 /B2	20 s
Brilliant Flow A3 /D3	20 s
Brilliant Flow A3.5 /B3	20 s
Brilliant Flow A4 /M5	20 s
Brilliant Flow Super White	20 s
Brilliant Flow Transparent	20 s

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Finishing <p>For the tooth preparation, appropriate rotary instruments and polishers from the Diatech/Alpen range can be used. Fluoridate all surfaces after completion.</p>	
Application of Brilliant Flow for direct restorations	
Preparation of the cavity <p>When preparing the cavity, every effort possible should be made to use a technique that preserves tooth structure (principle of the Adhesive Restoration Technique). Prepare enamel and dentin using 80 µm preparation diamonds and finish using 25 µm finishing diamonds. Bevelled enamel margins are recommended to increase the area of adhesion between tooth and filling material and thereby optimize the marginal seal.</p>	
Pulp Protection <p>Covering the areas close to the pulp with hard-setting calcium hydroxide cements using the spot technique protects the pulp against bacterial infiltration.</p>	
Application of Brilliant Flow <p>Apply the material using the mounted applicator tip directly into the cavity.</p>	
Application of Brilliant Flow for indirect restorations	
Pretreat contact surfaces of cavity according to the corresponding manufacturer's instructions for use.	
Application of Brilliant Flow <p>Dispense the material using the mounted applicator tip directly into the restoration or cavity preparation. Position restoration in place using gentle pressure.</p>	
Removal of excess material <p>Rough excess material (e.g. brush or spatula). Afterwards, hold restoration with pressure in position and remove additional excess material with a spatula or a suitable instrument.</p>	
Polymerization <p>Recommended curing time: cure each side for 40 s.</p>	
Emergency measuress <p>In the case of direct contact with the oral mucosa, rinsing with tap water is sufficient. In case of contact with the eyes, rinse thoroughly with water (10 min), and consult an ophthalmologist showing these instructions for use.</p>	

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Flow is cured in a tooth cavity or on a model with all current light-curing lamps (e.g. Coltolux® LED).

Composition <p>Metacrylates Barium glass, silanized Amorphous silica, hydrophobic</p>	
Technical data <p>Complies with ISO 4049 Average filler particle size: 0.6 µm Range of particle size: 0.04–2.5 µm</p>	
Brilliant Flow <p>Percentage by volume of total inorganic filler: 42% Percentage by weight of total inorganic filler: 63%</p>	
Indications <ul style="list-style-type: none">Direct fillings of class V (cervical caries, root erosion, wedge-shaped defects) Repairs to anterior teeth (class III, IV) Small fillings in the posterior region Fillings of minimal cavities Repairs of composite and ceramic veneers Block out of undercuts Adhesive cementing of composite and ceramic restorations, if light penetration is possible. Preventive resin restorations	
Contra-indications <p>If any allergies exist to any components within Brilliant Flow. If the application field cannot be isolated during application and curing of Brilliant Flow. For patients with poor oral hygiene.</p>	
Side effects <p>No systemic side effects are known. Contact allergies with products of similar composition have been reported in isolated cases.</p>	
Interactions with other agents <p>Agents containing eugenol and/or oil of cloves may affect the polymerization of Brilliant Flow. The use of zinc-oxide-eugenol cements in combination with Brilliant Flow should therefore be avoided. Discolorations may occur when using cationic mouth rinses as well as plaque indicators or chlorhexidin.</p>	

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