Excite F

[en] Instructions for Use

Dental adhesives

Rx ONLY

C € 0123

Date information prepared:
2023-07-06 / Rev. 0

Manufacturer:
Ivoclar Vivadent AG
Bendererstrasse 2
9494 Schaan / Liechtenstein
www.ivoclar.com

ivoclar

English

1 Intended use

Intended purpose

Bonding of dental restorations to dental hard tissue

Patient target group

- Patients with permanent teeth
- Patients with deciduous teeth

Intended users / Special training

- Dentists
- No special training required.

Use

For dental use only.

Description

The product is a light-curing, single-component adhesive for dentin and enamel bonding in conjunction with the total-etch technique.

Indications

Missing tooth structure in anterior and posterior teeth

Areas of application

- Light-curing adhesive for direct composite restorations
- Light-curing adhesive for indirect all-ceramic and composite restorations

Contraindications

- The placement of the product is contraindicated if a patient is known to be allergic to any of the material's ingredients.
- Do not use the product for direct pulp capping

Limitations of use

- If the stipulated working technique cannot be applied.
- The product must not be used in combination with purely self-curing composite materials. Dual-curing materials always have to be light-cured.
- Do not use the product as a primer for ceramic restorative materials.
A suitable ceramic primer must be used.
- The applicators are for single use only and must not be reused due to hygienic reasons.
- The Soft Touch Single-Dose units are for single use only and must not be reused.

Side effects

In rare cases, components of the product may lead to sensitization. In such cases, the material should no longer be used.

Interactions

- Peroxide-based bleaching agents and disinfectants may inhibit the polymerization of adhesives. Such materials must not be used ≤ 14 days prior to the application of the product.
- Phenolic substances (e.g. eugenol) may inhibit the polymerization of the product.
- Before any treatment with the product, completely rinse off all iron- and aluminium-based styptic agents with water.
- Fluoridation agents must not be used prior to the application of the product.

Clinical benefit

- Reconstruction of chewing function
- Restoration of esthetics

Composition

Bis-GMA, ethanol, HEMA, phosphoric acid acrylate, GDMA, UDMA, phenylbis (2,4,6-trimethylbenzoyl)-phosphine oxide.

2 Application

2.1 Isolation

Ensure adequate relative or absolute isolation.

2.2. Pulp protection / Cavity lining

In very deep cavities, selectively coat areas close to the pulp with a calcium hydroxide liner and subsequently cover with a pressure-resistant cement (e.g. a glass ionomer cement).

2.3 Conditioning with phosphoric acid gel

Please observe the instructions for use of the phosphoric acid gel.

- Apply phosphoric acid gel onto the prepared enamel first, and then to the dentin. The etchant should be left to react on the enamel for 15–30 s and on the dentin for 10–15 s.
- Rinse off the phosphoric acid gel with a vigorous water spray for at least 5 seconds.
- Note: Do not overdry the dentin.
Remove any excess moisture, leaving the dentin surface with a glossy wet appearance (wet bonding). This can be done with either an air gun, a dry brush, a foam pellet or other lint-free absorbent wipes.

2.4 Dosage

Bottle

- Note: Protect the product from light. Exposure to light may lead to premature polymerization.
- Dispense the required amount of material into a mixing well.
- Pick up and apply the product using a customary disposable applicator.
- Seal the bottle tightly after each use.

Soft Touch single-dose units with pre-mounted applicator

Activate the Soft Touch single-dose unit to prepare it for use as follows:

- Place the single-dose unit on a firm surface.
- Press the pre-mounted applicator brush into the single-dose unit until a noticeable resistance is overcome.

2.5 Application of the adhesive

2.5.1 Direct restorations

- Note: Make sure that all the cavity walls are completely covered.
- Saturate enamel and dentin with a generous amount of Excite F and agitate the adhesive on the prepared surfaces for at least 10 seconds
- Disperse any excess of the product to a thin layer with a weak stream of water- and oil-free air. Pooling should be avoided. A uniform, glossy appearance shows that the tooth surface is completely sealed.

2.5.2 Indirect restorations

- Prepare the surfaces of the composite or ceramic restoration according to the instructions of the manufacturer.
- Note: Make sure that all the cavity walls are completely covered.
- Saturate enamel and dentin with a generous amount of Excite F and agitate the adhesive on the prepared surfaces for at least 10 seconds.
- Disperse any excess of the product to a thin layer with a weak stream of water- and oil-free air. Pooling should be avoided. A uniform, glossy appearance shows that the tooth surface is completely sealed.

2.6 Curing

A curing time of 10 s at a light intensity of between 500–1300 mW/cm² is sufficient.

If the diameter of the light guide tip is smaller than the restoration diameter, perform overlapping polymerization in order to ensure that all areas of the restoration are covered.

 **Please observe the instructions for use of the curing light employed.**

Safety note

Avoid direct exposure of the gingiva, mucous membrane or skin to the light emitted by the curing light.

2.7 Next steps

- Proceed according to the instructions for use of the composite filling material / composite luting material used.
- Seat the restoration using a light-/ dual-curing luting material and light-cure according to the instruction of the manufacturer.

Application notes

- Do not expose the product to intensive light, as this will considerably reduce the working time.
- Do not use the bottle intraorally. For hygienic reasons, suitable application aids must be used.
- Avoid pooling, as this may compromise the accuracy of fit of the restoration.
- Prevent any contamination of the applicator with blood, saliva, sulcus fluid or water during application. In the event of contamination, the cavity must be rinsed again and the procedure must be repeated with a new applicator.

3 Safety information

- The product is an irritant. Contact with the skin, mucous membrane or eyes should be avoided.
- In case of skin contact, immediately wash with copious amounts of water.
- If the material should come in contact with the eyes, immediately rinse with copious amounts of water and consult a physician or ophthalmologist.
- In rare cases, contact with the skin may cause sensitization to the ingredients.
- Commercial medical gloves do not provide protection against the sensitizing effects of methacrylates.
- In the case of serious incidents related to the product, please contact Ivoclar Vivadent AG, Bendererstrasse 2, 9494 Schaan/ Liechtenstein, website: www.ivoclar.com and your responsible competent authority.
- The current Instructions for Use are available in the download section of the Ivoclar Vivadent AG website (www.ivoclar.com).
- The Summary of Safety and Clinical Performance (SSCP) can be retrieved from the European Database on Medical Devices (EUDAMED) at <https://ec.europa.eu/tools/eudamed>.
- Basic UDI-DI: 76152082ABOND011HM

Warnings

- Observe the Safety Data Sheet (SDS) (available in the download section of the Ivoclar Vivadent AG website www.ivoclar.com).
- Warnings for phosphoric acid gel: Phosphoric acid gel is corrosive. Avoid contact with the eyes, mucous membrane and skin (goggles are recommended for both patients and operators). If the material should accidentally come in contact with the eyes, immediately rinse with copious amounts of water and consult a physician / ophthalmologist. Observe the Safety Data Sheet (SDS) (available at www.ivoclar.com).

Disposal information

Remaining stocks must be disposed of according to the corresponding national legal requirements.

Residual risks

Users should be aware that any dental intervention in the oral cavity involves certain risks. The following clinical residual risks are known:

- Failure of the adhesive bond
- Postoperative sensitivity

4 Shelf life and storage

- Close bottles immediately after use; exposure to light causes premature polymerization.
- Storage temperature 2–28 °C (36–82 °F)
- If the indicated storage temperature is exceeded, proper functioning of the product cannot be ensured.
- Do not use the product after the indicated date of expiration.
- Expiry date: See note on packaging.
- Before use, visually inspect the packaging and the product for damage. In case of any doubts, please contact Ivoclar Vivadent AG or your local dealer.

5 Additional information

Keep material out of the reach of children!

Not all products are available in all countries.

The material has been developed solely for use in dentistry. Processing should be carried out strictly according to the Instructions for Use. Liability cannot be accepted for damages resulting from failure to observe the Instructions or the stipulated area of application. The user is responsible for testing the materials for their suitability and use for any purpose not explicitly stated in the Instructions.

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