Prior to use, carefully read the instructions for use. **EN**

GC Fuji™ II LC (Improved)

RADIOPAQUE LIGHT CURED REINFORCED GLASS IONOMER RESTORATIVE

For use as a dental professional in the recommended indications.

RECOMMENDED INDICATIONS

- Class III and V restorations, particularly for cervical erosions and root surface caries.
- Restoration of primary teeth.
- Core build up.
- Cases where radiopacity is required.
- Esthetic applications.
- As a base or liner.

CONTRAINDICATIONS

- Pulp capping.
- Use of this product in patients with known allergies to glass ionomer cement, methacrylate monomer or methacrylate polymer.

DIRECTIONS FOR USE

Powder / Liquid Ratio (g/g)	3,2 / 1,0
Mixing Time (sec.)	20" 25"
Working Time (min., sec.)	3'45"
Light Curing Time (sec.)	20"
Depth of Cure (A2) (mm)	1,8

Test conditions: Temperatura (23+/- 1°C)
Relative humidity (50+/- 10%)

- POWDER AND LIQUID DISPENSING**
 - Select shade based on Vita® shade guide.
 - Vita is a registered trademark of Vita Zahnfabrik, Bad Säckingen, Germany.
 - The standard powder to liquid ratio is 3,2g / 1,0g. 1 level scoop of powder to 2 drops of liquid.
 - For accurate dispensing of powder, tap the bottle gently.
 - Do not shake the powder container.
 - Fill the liquid bottle vertically and squeeze gently.
 - Close bottles immediately after use.
- MIXING**
 - Fluff powder of bottle before dispensing (Fig. 1).
 - Divide one scoop of powder and two drops of liquid on pad.
 - Mix the powder and liquid ratio is 3,2g / 1,0g.
 - Spread liquid out into a thin layer (about the size of half dollar or 3 cm) with plastic spatula (Fig. 3).
 - Divide half of the powder and mix with lapping stoches (like impression material) for 10 to 15 seconds (Fig. 4).
 - Pull in remaining powder and mix thoroughly to a glossy consistency. Do not exceed 20-25 second total mixing time (Fig. 5).
- RESTORATIVE TECHNIQUE**
 - Prepare tooth and standard techniques. Extensive mechanical treatment is necessary. For pulp capping use calcium hydroxide.
 - Apply CAVITY CONDITIONER (10 seconds) or DENTIN BONDING COAT (20 seconds) to the bonding surfaces using a cotton pellet or sponge.
 - Rinse thoroughly with water. Dry by blotting with a cotton pellet or gently blowing with an air syringe. DO NOT DESICCATE. Best results are obtained when prepared surfaces appear moist (glistening).
 - Mix the required amount of cement. Working time is 45 seconds from the start of mixing at 23°C (73,4°F). Higher temperatures will shorten working time.
 - Transfer the cement to a syringe or suitable placement instrument. Avoid air bubbles.
 - Form the contour and place a transparent mat if required.
 - Light cure for 20 seconds with a light curing device (470 nm wavelength). Place light source as closely as possible to the cement surface.
 - For cavities deeper than 1,8 mm, use a layering technique.
- FINISHING**
 - Remove matrix and finish under water spray using standard techniques.
 - Note.
 - Apply a final coat of GC Fuji VARNISH (blow dry) or GC Fuji COAT LC (light cure for 10 sec) or G-COAT PLUS (light cure for 20 sec.)

STORAGE
Recommended for optimal performance, store in a cool and dark place (4-25°C / 39-77,0°F).

SHADE
(11 shades) A1, A2, A3, A3.5, A4, B2, B3, B4, C2, C4, D2.

PACKAGES
1. 3-2 package: 15 g powder (3) (A2, A3, B3), 8 g (6,8 ml) liquid (2), powder scoop (1), mixing pad (No.22) (1), plastic spatula (1) 1.
2. Bottle of 15 g powder with skin, 3.
3. Bottle of 8 g (6,8 ml) liquid

CAUTION
1. In case of contact with oral tissue or skin, remove immediately with a sponge or cotton soaked in alcohol. Flush with water.
2. In case of contact with eyes, flush immediately with water and seek medical attention.
3. Do not mix powder or liquid with other glass ionomer components.
4. Personal protective equipment (PPE) such as gloves, face masks and safety eyewear should always be worn.
5. In rare cases the product may cause sensitivity in some people. If any such reactions are experienced, discontinue the use of the product and refer to a physician.

Some products referenced in the present IFU may be classified as hazardous according to GHS. Always familiarize yourself with the safety data sheets available at:
<http://www.gceurope.com>
They can also be obtained from your supplier.

CLEANING AND DISINFECTING.
MULTI-USE DISINFECTING SYSTEMS: to avoid cross-contamination between patients this device requires mid-level disinfection. Immediately after use inspect device and label for deterioration. Discard device if damaged.
AUTOMATICALLY CLEANING: thoroughly clean device to prevent drying and accumulation of contaminants.
Disinfect with a mid-level registered healthcare-grade infection control product according to regional/national guidelines.

Undesired effects- Reporting:
If you become aware of any kind of undesired effect, reaction or similar events experienced by use of this product, including those not listed in this instruction for use, please report them directly through the relevant vigilance system, by selecting the proper authority of your country accessible through the following link:
https://ec.europa.eu/growth/sectors/medical-devices/contacts_en as well as to our internal vigilance system: vigilance@gcc.dental. In this way you will contribute to improve the safety of this product.
Last revised: 06/2019

GC 2797
MANUFACTURED BY
GC CORPORATION
76-1 Hasunuma-cho, Itabashi-ku, Tokyo 174-8585, Japan
DISTRIBUTED BY
GC CORPORATION
76-1 Hasunuma-cho, Itabashi-ku, Tokyo 174-8585, Japan
EU: GC EUROPE N.V.
Nieuwland 1, Nieuwland 2400, Interleuvenlaan 33,
8-3000 Leuven, Belgium TEL: +32 16 74 10 00

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