Speedcem® Plus

[en] Instructions for Use

Self-adhesive polymer-based dental luting material
(intraoral dual-curing)

Rx ONLY

C € 0123

Date information prepared:
2023-12-13 / Rev. 0

Manufacturer:
Ivoclar Vivadent AG
Bendererstrasse 2
9494 Schaan / Liechtenstein
www.ivoclar.com

ivoclar

English

1 Intended use

Intended purpose

Self-adhesive cementation of indirect anterior and posterior restorations

Patient target group

Patients with permanent teeth

Intended users / Special training

- Dentists
- No special training required.

Use

For dental use only.

Description

Speedcem® Plus is a self-adhesive, polymer-based dental luting composite (intraoral dual-curing) for the cementation of indirect restorations made of high strength all-ceramics, metal-ceramics and metal. Additional dentin bonding agents are not required. The special filler composition gives the luting material a very high radiopacity.

Shades

The product is available in three shades with different degrees of translucency:

- transparent
- yellow
- opaque

Working time

The working and setting times depend on the ambient temperature. Once the product has been dispensed from the automix syringe, the following times apply:

	At room temperature: approx. 23 °C / 73 °F	Intraorally
Working time	approx. 2 min	approx. 1 min
Setting time (incl. working time)	approx. 6 min	approx. 3 min

Note

After Speedcem Plus has been dispensed from the automix syringe, the curing process can be accelerated by intensive illumination (operating light, ambient light).

Indications

- Missing tooth structure in anterior and posterior teeth
- Partial edentulism in the anterior and posterior region

Areas of application:

- Permanent cementation on natural teeth with adequate retentive features in conjunction with indirect restorations made of:
 - oxide ceramics, e.g. zirconium oxide (crowns, bridges)
 - metal and metal-ceramics (inlays, onlays, crowns, bridges, endodontic posts)
 - lithium disilicate glass-ceramics (crowns, bridges)
 - fibre-reinforced composites (endodontic posts)
- Permanent cementation of crowns and bridges on the following implant abutment materials:
 - oxide ceramic (e.g. zirconium oxide)
 - metal (e.g. titanium)
 - lithium disilicate glass-ceramic

Contraindications

The use of the product is contraindicated if the patient is known to be allergic to any of its ingredients.

Limitations of use

- If the preparation/shape of the implant abutment does not provide adequate retention (e.g. veneers, short or severely tapered tooth preparations/implant abutments).
- In general, Speedcem Plus should not be applied to the exposed pulp or to dentin which is close to the pulp.
- Do not use if a dry working field cannot be established or the prescribed application technique cannot be applied.
- Use a suitable hygienic protective sleeve for the syringe as indicated by the manufacturer.
- The mixing tips and root canal tips are not intended to be reused.

Side effects

There are no known side effects to date.

In individual cases, allergic reactions to individual components have been reported.

Interactions

- Phenolic substances (e.g. eugenol, wintergreen oil) inhibit polymerization. Consequently, the use of products, e.g. mouth rinses and temporary cements, which contain these components should be avoided.
- Disinfectants with an oxidative effect (e.g. hydrogen peroxide) may interact with the initiator system, which in turn may impair the curing process.
Therefore, do not disinfect the preparation using oxidative agents.
- Alkaline jet media may compromise the effect of Speedcem Plus.

Clinical benefit

- Reconstruction of chewing function
- Restoration of esthetics

Composition

The monomer matrix is composed of dimethacrylates and acidic monomers. The inorganic fillers are barium glass, ytterbium trifluoride, co-polymer and highly dispersed silicon dioxide. Additional contents are initiators, stabilizers and colour pigments (< 1 %).

The primary particle size of the inorganic fillers is between 0.1 and 7 µm. The mean particle size is 5 µm.

The total content of inorganic fillers is approx. 40%vol.

2 Application

For more detailed information also refer to the Instructions for Use of the products used in conjunction with Speedcem Plus.

2.1 Pre-treatment of the prepared tooth or implant abutment

2.1.1 Pre-treatment of the prepared tooth

Removal of the temporary restoration and thorough cleaning of the prepared tooth

Remove possible residue of the temporary luting cement from the prepared tooth with a polishing brush and an oil- and fluoride-free cleaning paste (e.g. Proxylt® fluoride-free prophylaxis paste). Rinse with water spray. Subsequently, dry with oil- and moisture-free air. Avoid overdrying.

Note: Cleaning with alcohol can lead to dehydration of dentin.

Try-in of the restoration and isolation

Next, the shade, accuracy of fit and occlusion of the restoration can be checked. Care should be taken when checking the occlusion of fragile and brittle ceramic objects, as there is a risk of fracture before they are permanently cemented.

It is essential to ensure reliable isolation of the operating field

with e.g. OptraGate®, cotton rolls, saliva ejector and absorbent pads. Dental hard tissue that has become contaminated with saliva needs to be cleaned again.

2.1.2 Pre-treatment of the implant abutment

Removal of the temporary restoration

Remove the temporary restoration and, if present, the temporary abutment. Rinse the implant lumen and peri-implant gingiva.

Insertion and inspection of the final abutment

According to the manufacturer's directions.

A retraction cord may be placed to be better able to check the accuracy of fit of the restoration and remove excess luting composite.

Try-in of the restoration and isolation

Next, the shade, accuracy of fit and occlusion of the restoration can be checked. Care should be taken when checking the occlusion of fragile and brittle ceramic objects, as there is a risk of fracture before they are permanently cemented. It is essential to ensure reliable isolation of the operating field – preferably with OptraGate, cotton rolls, saliva ejector and absorbent pads – when using adhesive cementation with composites.

Pre-treatment of the abutment surface

According to the manufacturer's directions.

Note: Contamination with saliva or blood must be avoided during and after the pre-treatment of the abutment (try-in). If necessary, the abutment should be cleaned again in situ using water and fluoride-free cleaning paste (e.g. fluoride-free Proxylt).

A) Abutments made of zirconium oxide and titanium:

Clean or sandblast the abutment surface extraorally as described in sections 2.2.1/2.2.2.

B) Abutments made lithium disilicate glass-ceramic:

Pre-treat extraorally as described in section 2.2.3 using e.g. Monobond Etch & Prime®

Sealing of the abutment

- Seal the screw channel of the abutment using e.g. a temporary restorative material (e.g. Telio® Inlay / Onlay). It is recommended to isolate the screw head from the temporary restorative material with cotton pellets.
- Rinse the abutment with water spray.
- Dry the abutment.

2.2 Cleaning and pre-treatment of the restoration

According to the manufacturer's directions.

Note: Contamination with saliva or blood must be avoided during and after pre-treatment of the restoration (try-in).

2.2.1 Restorations made of oxide ceramics - (e.g. zirconium oxide, IPS e.max® ZirCAD)

IMPORTANT! In order to achieve an optimum bond, do not clean surfaces with phosphoric acid.

If the restoration has been sandblasted in the laboratory:

- Apply Ivoclean, as described in section 2.2.4, to clean the restoration.

If the restoration has not been sandblasted in the laboratory:

- Sandblast the inner surface of the restoration (use the blasting parameters stated in the manufacturer's instructions for use of the restorative material, e.g. IPS e.max® ZirCAD, max. 1 bar/ 14.5 psi, Al₂O₃).
- Clean the restoration e.g. in an ultrasonic unit for approx. 1 minute.
- Thoroughly rinse the restoration with water spray and dry with water-/oil-free air.

2.2.2 Metal or metal-supported restorations

IMPORTANT! In order to achieve an optimum bond, do not clean surfaces with phosphoric acid.

If the restoration has been sandblasted in the laboratory:

- Apply Ivoclean, as described in section 2.2.4, to clean the restoration.

If the restoration has not been sandblasted in the laboratory:

- Sandblast the inner restoration surface (use the blasting parameters stated in the manufacturer's instructions for use of the restorative material).
- Clean the restoration e.g. in an ultrasonic unit for approx. 1 minute.
- Thoroughly rinse the restoration with water spray and dry with water-/oil-free air.

If the restoration comprises precious metal, apply Monobond Plus to the pre-treated surfaces using a brush or microbrush, allow to react for 60 seconds and then disperse with a strong stream of air.

2.2.3 Restorations made of lithium disilicate glass-ceramic (e.g. IPS e.max® Press/CAD)

A) Procedure using Monobond Etch & Prime regardless of whether the restoration has been pre-treated in the laboratory:

- Following try-in, thoroughly rinse the restoration with water spray and dry with oil- and moisture-free air.
- Apply Monobond Etch & Prime onto the bonding surface using a microbrush and agitate it into the surface for 20 seconds using light pressure.
- Allow to react for another 40 s.
- Thoroughly rinse off Monobond Etch & Prime with water until the green colour has disappeared. If any residue remains in the microporosities, the restoration may be cleaned with water in an ultrasonic unit for up to 5 min.
- Dry the restoration with a strong stream of oil- and moisture-free air for approx. 10 s.

B) Procedure using Monobond Plus:

If the restoration has been pre-treated in the laboratory

- If the restoration has already been pre-treated in the laboratory, it should be cleaned after the try-in using Ivoclean as described above.
- Subsequently, re-apply Monobond Plus to the cleaned surface using a brush or a microbrush, allow to react for 60 seconds and then disperse with a strong stream of air.

If the restoration has not been pre-treated in the laboratory:

- Etch the restoration with 5% hydrofluoric acid (e.g. IPS® Ceramic etching gel) for 20 s or according to the instructions of the manufacturer of the restorative material being used.
- Thoroughly rinse the restoration with water spray and dry with water-/oil-free air.
- Apply Monobond Plus to the pre-treated surfaces with a brush or microbrush, leave to react for 60 seconds and then disperse with a strong stream of air.

2.2.4 Cleaning of restorations contaminated with blood or saliva

Regardless of any previous conditioning, clean contaminated restoration surfaces extraorally using Ivoclean as follows:

- Following try-in, thoroughly rinse the restoration with water spray and dry with oil-free air.
- Cover the entire bonding surface of the restoration with a layer of Ivoclean using a microbrush or brush.
- Allow 20 seconds for Ivoclean to take effect, then thoroughly rinse with water spray and dry with oil-free air.
- Restorations made of lithium disilicate glass-ceramic (e.g. IPS e.max Press/CAD) must be primed again with Monobond Etch & Prime or Monobond® Plus.

2.3 Application of Speedcem Plus into the restoration

Note: Make sure the product has reached ambient temperature before use. For each application, place a new automix tip on the double syringe. Extrude Speedcem Plus from the automix syringe, apply the desired quantity directly into the restoration and cover the entire bonding surface. As the luting material will cure in the used mixing tip, it may serve as a seal for the contents of the syringe until the next application.

2.4 Placement of the restoration and removal of excess cement

Note: As with all composites, Speedcem Plus is subject to oxygen inhibition. This means that the surface layer does not polymerize during curing, as it comes in contact with atmospheric oxygen. To prevent this, cover the restoration margins with glycerine gel/air block (e.g. Liquid Strip) immediately after excess removal. After complete polymerization, the glycerine gel/air block is rinsed off with water.

2.4.1 Solely self-curing

- Seat the restoration and retain it in place exerting uniform pressure until the excess material has been removed and the cement has completely cured (approx. 3 min.)
- Remove excess material e.g. with a microbrush/brush/foam pellet/dental floss or scaler. Make sure to remove excess material in time, especially in areas that are difficult to reach (proximal areas, gingival margins, pontics).
- To remove excess material from implant abutments, special implant scalers/curettes should be used to avoid scratching of the implant surface.

2.4.2 Self-curing with additional light-activation to tack-cure excess material or accelerate the curing process (quarter technique, indicated for cases with up to 2 bridge abutments = 3- to 4-unit bridge)

- Seat the restoration and retain it in place exerting uniform pressure until the excess material has been removed and the cement has completely cured.
- Excess cement is light-cured with a polymerization light (approx. 650 mW/cm²) from a distance of approx. 0-10 mm for 1 second per quarter surface (mesio-oral, disto-oral, mesio-buccal, disto-buccal).
- Following this procedure, the excess material can be easily removed with a scaler.
- Make sure to remove excess material in time, especially in areas that are difficult to reach (proximal areas, gingival margins, pontics).
- To remove excess material from implant abutments, special implant scalers/curttes should be used to avoid scratching of the implant surface.
- Subsequently, light-cure all cement lines again for 20 seconds (approx. 1,100 mW/cm²).
If the light output is lower, curing takes longer.

2.5 Finishing of the completed restoration

- Check occlusion and functional movements and make adjustments if necessary.
- Finish the cement lines with finishing diamonds if necessary.
 - Smooth out the cement lines using finishing and polishing strips and polish them with suitable polishing instruments (e.g. OptraGloss®).
 - If necessary, finish the restoration margins using suitable polishers (ceramics: e.g. OptraGloss).

Special instructions for the cementation of endodontic posts

- For the cementation of endodontic posts, carefully clean the root canal to remove any residual root canal filling material. (Residue of eugenol-based sealers may inhibit the polymerization of the luting composite.)
- Wet the endodontic post, which has been prepared according to the instructions of the post manufacturer, with the mixed Speedcem Plus and apply Speedcem Plus directly into the root canal using endodontic tips. The material should be applied in a sufficient quantity to provide excess cement.
- Insert the endodontic post.
- Remove excess cement from the occlusal preparation surface. Then, light-cure Speedcem Plus for 20 seconds.
- Pre-treat the occlusal preparation surface with an adhesive (e.g. Adhese Universal) according to the respective Instructions for Use.
- Apply the core build-up material directly onto the occlusal preparation surface and cure according to the manufacturer's instructions (light-curing).

3 Safety information

- In the case of serious incidents related to the product, please contact Ivoclar Vivadent AG, Bendererstrasse 2, 9494 Schaan/Liechtenstein, website: www.ivoclar.com, and your responsible competent authority.
- The current Instructions for Use are available on the website (www.ivoclar.com).
- Explanation of the symbols: www.ivoclar.com/elFU
- The Summary of Safety and Clinical Performance (SSCP) can be retrieved from the European Database on Medical Devices (EUDAMED) at <https://ec.europa.eu/tools/eudamed>.
- Basic UDI-DI: 76152082ACEME002DX

Warnings

- Observe the Safety Data Sheet (SDS) (available at www.ivoclar.com).
- Avoid contact of uncured Speedcem Plus with the skin/mucous membrane or eyes.
- Unpolymerized Speedcem Plus may cause slight irritation and may lead to a sensitization against methacrylates.
- Customary medical gloves do not provide protection against the sensitizing effects of methacrylates.

Disposal information

Remaining stocks must be disposed of according to the corresponding national legal requirements.

Residual risks

Users should be aware that any dental intervention in the oral cavity involves certain risks.

The following clinical residual risks are known:

- Failure of the adhesive bond

4 Shelf life and storage

- Storage temperature: 2-8 °C / 36-46 °F
- Store automix syringe with mixing tip attached after use.
- Do not use the product after the indicated expiry date.
- Expiry date: see information on syringes and packages.
- Before use, visually inspect the packaging and the product for damage. In case of any doubts, please contact Ivoclar Vivadent AG or your local dealer.

5 Additional information

Keep material out of the reach of children!

The material has been developed solely for use in dentistry. Processing should be carried out strictly according to the Instructions for Use. Liability cannot be accepted for damages resulting from failure to observe the Instructions or the stipulated area of application. The user is responsible for testing the products for their suitability and use for any purpose not explicitly stated in the Instructions. Descriptions and data constitute no warranty of attributes and are not binding.

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