VOCO
Ufi Gel® hard

EN Instructions for use
MD EU Medical device

Product description:
Ufi Gel hard is a cold-curing permanently hard relining material for dentures and is polymethacrylate-based. It is simple and quick to use for direct as well as indirect relinings.
The Ufi Gel hard liquid contains no methylmethacrylates, thus minimising the risk of allergies and irritations of the mucosa. The set additionally contains a conditioner to achieve a permanent and stable bond.

Indications:
- hard, permanent, total or partial relining to restore the function of partial and total dentures
- lengthening of denture margins

Contraindications:
Ufi Gel hard contains methacrylate and benzoyl peroxide. Ufi Gel hard should therefore not be used in patients with a known hypersensitivity (allergy) to these ingredients.

Patient target group:
Ufi Gel hard is suitable for use in all patients without any age or gender restrictions.

Performance features:
The product's performance features satisfy the requirements of the intended use and the relevant product standards.

User:
Ufi Gel hard should only be used by a professionally trained dental practitioner.

Use:

1. Preparation of the denture:
Check occlusion and carry out corrections, if necessary. Disinfect the prosthesis and dry thoroughly afterwards. Clean denture thoroughly with a brush and dry. Roughen all areas to be relined incl. buccal and labial margins with a suitable bur. Then clean and dry. Insulate all areas not to be relined (e.g. artificial teeth, labial and buccal areas beyond the roughened surfaces), e.g. with vaseline, to prevent unwanted contact with Ufi Gel hard.

Clasps, anchors and attachments of partial dentures have to be blocked out with thin-flowing silicone or wax with a low melting point.

2. Application of the conditioner
Apply conditioner with enclosed brush on all surfaces to be relined and let it dry in the air (approx. 30 s). First coat the labial and buccal surfaces, then the base of the denture. The conditioner remains effective for about 10 min after application.

Attention: Close the bottle tightly immediately after use because of high volatility.
The brush can be cleaned with alcohol, for example.

3. Dosing and mixing
Ufi Gel hard is mixed in a ratio of 1 ml liquid to 3 ml (≈ 1.8 g) powder. This corresponds to 2 graduation marks of the dropper to 1 graduation mark of the glass cylinder.

2 graduation marks liquid : 1 graduation mark powder

A thicker consistency, e.g. for lengthening the margins of the denture, is achieved by increasing the amount of powder.
Extraxt liquid with the dropper and drop it into a PP mixing cup. Shake the powder briefly and dispense it into the glass cylinder. For exact dosing smooth the surface of the powder by lightly tapping the side of the glass cylinder. Place powder into the liquid and mix carefully with a plastic spatula until homogeneous (approx. 10 s). Prevent formation of bubbles by stirring slowly and along the side of the mixing cup. Make any potential bubbles rise to the surface by tapping the mixing cup. Let the material soak until a workable consistency is achieved (approx. 1.5 min after commencement of mixing).

4. Application of Ufi Gel hard
Apply the material evenly onto the prepared margin and/or base of the denture with a plastic spatula. Prevent build-up of excess material, remove with a suitable instrument if necessary.
Re-insert denture and have the patient exert slight pressure in terminal occlusion for 1 min. Then carry out functional, chewing and swallowing movements for 2 min. Ensure that no material flows into the throat when relining an upper denture, especially at the transition from hard to soft palate (A-line).

4.1 Partial and total denture with undercuts
Remove denture 5 min after commencement of mixing and remove excess material immediately with scissors or a scalpel. Insert again into the mouth for a final occlusion check and allow to cure for about 2 - 3 min.
Instead of re-inserting the denture, curing can also be completed in warm water, e.g. in a pressure pot at approx. 40 °C. Do not let the material cure in contact with air since atmospheric oxygen will cause an uncured inhibition layer on the surface.

4.2 Total dentures without undercuts
Excess material can be removed intra- or extraorally before final polymerisation.
Intraorally: remove excess material after 5 min (from beginning of mixing) with a suitable instrument. Let the denture cure for further 2 - 3 min in the mouth until Ufi Gel hard is completely cured.
Extraorally: see section 4.1

5. Finishing the relining
The relined denture can be finished and polished with the usual instruments (tungsten carbide bur, rubber polisher, polishing disc).

Table with 3 columns: Action, Time, and Note. Rows include: beginning of mixing (2 min, mix and apply), insert denture into the mouth (1 min exert pressure, 2 min functional movements), remove denture if necessary (1.5 min, remove excess material), insert denture into the mouth (2 min, curing), remove denture (1 min, finishing, polishing), and a final row for 9.5 min.

extraorally intraorally intra- or extraorally

Warnings, precautionary measures:
- Store bottles of liquid and conditioner carefully closed and in upright position.
- Avoid contact of liquid or conditioner with skin. Rinse sites of contact thoroughly with water and soap.
- In case of eye contact, rinse well with plenty of water and consult an ophthalmologist.
- The indicated time periods refer to a room temperature of 25 °C and a temperature of 35 °C in the mouth. Polymerisation will be slower at lower temperatures, higher temperatures will accelerate polymerisation.
- Excessively long or violent mixing may lead to air bubbles and an un-homogeneous consistency, which produces a rough surface.
- Use the mixing cup several times.
- The enclosed dropper are intended for a single use and must be disposed of after only one use.
- Our information and/or advice do not relieve you of the obligation of checking that the products supplied by us are suitable for their intended purpose.

Storage:
Storage at 4°C-23°C. Do not use after the expiry date.

Cleaning indications:
Dentures relined with Ufi Gel hard can be cleaned with the usual cleansers and procedures. This applies to domestic as well as professional cleaning.

Disposal:
Dispose of the product in accordance with local regulations.

Reporting obligation:
Serious events such as death, temporary or permanent serious deterioration of a patient's, user's or other person's health condition and a serious risk to public health that arise or could have arisen in association with the use of Ufi Gel hard must be reported to VOCO GmbH and the responsible authority.

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