Product description:

Final Varnish LC seals and smoothenes the surface of glass ionomer cements. The light-curing covering varnish **Final Varnish LC** protects the cement surface from damage due to drying and from the ingress of moisture or saliva until the cement is cured completely, and also produces a smooth surface.

Indications:

- Protection of glass ionomer cement surfaces from the effects of humidity and dryness immediately after placing the filling
- Sealing and smoothing of glass ionomer cement surfaces after polishing
- Sealing of the cervical margins of cementations immediately after removal of excess material
- Sealing of linings/build-up fillings before taking impressions

Contraindications:

Final Varnish LC contains methacrylate and BHT. **Final Varnish LC** should therefore not be used in patients with a known hypersensitivity (allergy) to these ingredients.

Patient target group:

Final Varnish LC is suitable for use in all patients without any age or gender restrictions.

Performance features:

The product's performance features satisfy the requirements of the intended use and the relevant product standards.

User:

Final Varnish LC should only be used by a professionally trained dental practitioner.

Use:

Dispense 1-2 drops onto a pad or the application instrument and close the bottle immediately. Apply **Final Varnish LC** with a brush or other appropriate instrument immediately after placing the filling or finishing the glass ionomer cement. Remove any air bubbles (with probe). Cover all glass ionomer cement surfaces with the varnish (including approximal surfaces and surfaces after removal of matrices or similar). If the glass ionomer cement surfaces require additional treatment within 24 h after the initial application, they must be resealed once again.

Light-curing:

Cure **Final Varnish LC** with halogen/LED light (10 s per surface). Conventional polymerisation devices are suitable for light-curing the material. The light output should be a minimum of 500 mW/cm² for halogen polymerisation devices and 300 mW/cm² for LED devices.

Warnings, precautionary measures:

- Expose **Final Varnish LC** to the operating light for a short time only, as otherwise polymerisation will start (viscosity increases).
- Check sulcus after application, remove any excess material.
- Avoid contact with the mucosa as much as possible.
- Apply **Final Varnish LC** again if the seal is damaged prior to final curing.
- Unpolymerised or incompletely dried varnish may impede the setting of addition-curing silicones.
- Our information and/or advice do not relieve you of the obligation of checking that the products supplied by us are suitable for their intended purpose.

Storage:

Carefully close the **Final Varnish LC** bottle immediately after use. Store **Final Varnish LC** at 4°C-23°C and protected from light. Do not use after the expiry date.

Disposal:

Dispose of the product in accordance with local regulations.

Reporting obligation:

Serious events such as death, temporary or permanent serious deterioration of a patient's, user's or other person's health condition and a serious risk to public health that arise or could have arisen in association with the use of **Final Varnish LC** must be reported to VOCO GmbH and the responsible authority.

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VC 60 BB1241 E1 0622 99 © by VOCO