Start-X®

EN	Instructions For Use	2
FR	Instructions d'utilisation	7
IT	Istruzioni per l'uso	12
DE	Gebrauchsanweisung	17
ES	Instrucciones de uso	23
PT	Instruções de utilização	28
NL	Gebruiksaanwijzing	33
SV	Användarinstruktioner	39
DA	Brugsanvisning.....	44
NO	Brukerinstruksjoner.....	49
LT	Naudojimo instrukcija	54
LV	Lietošanas norādījumi	59
ET	Kasutusjuhised	64
CS	Návod k použití.....	69
SK	Návod na použitie	74
HU	Használati utasítás	79
RO	Instrucțiuni de utilizare	84
SR	Uputstva za upotrebu	89
BS	Uputstva za upotrebu	94
CNR	Uputstva za upotrebu	99
SQ	Udhëzimet për përdorim	104
EL	Οδηγίες χρήσης.....	109
TR	Kullanım Kılavuzu	114
PT-BR	Instruções de uso.....	119
KO	사용 지침.....	124

Start-X®

FOR DENTAL USE ONLY

FOR PROFESSIONAL USE ONLY

NON-STERILE – MANDATORY PROCESSING BEFORE FIRST USE

REUSABLE – LIABLE FOR REPROCESSING

INSTRUCTIONS FOR USE

REUSABLE NON-STERILE ULTRASONIC DEVICES

REF. A0660 (Start-X® Compatible SATELEC), A0661 (Start-X® Compatible EMS), A9660 (Start-X® Trial Kit Compatible SATELEC), A9661 (Start-X® Trial Kit Compatible EMS)

1 PRODUCT DESCRIPTION

- Product description:

Start-X® ultrasonic inserts are stainless-steel devices composed of:

- A part to be screwed in an ultrasonic generator handpiece thanks to a stainless-steel wrench. This part is marked with the insert n°.
- An active part with an end tip whose shapes depend on the insert n° and clinical application. Tips are driven by the handpiece of an ultrasonic generator which transmit vibration.
- A water port, to cool down the insert and avoid overheating of the treated tooth.

Five types of Start-X® inserts are available, depending on their clinical application as indicated in the below table:

Tip n°	Tip diameter (mm)	Active part diameter (mm)	Active part length (mm)	Total length (mm)	Clinical application
Tip n°1	0.84	1.46	12.0	26.65	Refinement of access cavity walls
Tip n°2	1.04	1.40	8.00	26.70	MB2 canal scouting (2nd mesio buccal root)
Tip n°3	0.65	0.83	8.00	26.40	Removing obstructions to the canal orifice
Tip n°4	1.40	0.80	/	26.20	Metal post removal
Tip n°5	1.00	1.00	10.0	26.50	Re-locating the original pulp chamber floor anatomy (calcification, removal of old obturation material)

- Clinical benefits:

The clinical benefit for the patient being treated with Start-X® ultrasonic inserts is the restoration of mastication function and the conservation of the natural tooth.

2 INDICATIONS FOR USE

- Indication for use:
The products are indicated for treatment of endodontic disease.
- Intended purpose:
Ultrasonic driven devices intended to refine the access cavity, to remove obstructions and dental material (e.g., post, obturation), locate the canal orifice.
- Intended patient population:
All patients who need root canal treatment can be treated with the device (no restriction in terms of age, gender and origin). It is nevertheless recommended not to use these ultrasonic inserts on patients who have a Pacemaker.
- Intended user and environment of use:
Endodontic instruments are to be used only in a clinical or hospital environment, following good dental practice, by qualified dental professionals such as general practitioners as well as endo specialists (endodontist) and dental assistants.

3 CONTRAINDICATIONS

None known when used as recommended. It is nevertheless recommended not to use these ultrasonic inserts on patients who have a pacemaker.

4 DELIVERY FORM

Start-X® ultrasonic inserts are delivered in non-sterile plastic boxes and the trial kit is composed by an assortment sold with the DVD in a paperboard box.

Start-X® ultrasonic range comprises the following endodontic instruments:

SKU x: Thread y: Tip n°	Product name	Available Threads (x)	Available Tip n° (y)	Packaging and content
A066x00000y00	Start-X®	0: type M3 x 0.6 mm 1: type M3 x 0.5 mm	1: Tip n°1 2: Tip n°2 3: Tip n°3 4: Tip n°4 5: Tip n°5	Plastic box containing 1 insert
A066000090000	Start-X® Compatible SATELEC Endodontic Instrument	0: type M3 x 0.6 mm		Assortment with the following components: Insert 1, insert 2, insert 3, insert 4, insert 5, the wrench
A066100090000	Start-X® Compatible EMS Endodontic Instrument	1: type M3 x 0.5 mm		Assortment with the following components: Insert 1, insert 2, insert 3, insert 4, insert 5, the wrench
A966000090100	Start-X® Trial Kit Compatible SATELEC	0: type M3 x 0.6 mm		Trial kit with the following components: Insert 1, insert 2, insert 3, insert 4, insert 5, the wrench, a DVD
A966100090100	Start-X® Trial Kit Compatible EMS	1: type M3 x 0.5 mm		Trial kit with the following components: Insert 1, insert 2, insert 3, insert 4, insert 5, the wrench, a DVD

5 COMPOSITION

Start-X® ultrasonic inserts are stainless-steel devices. All parts described in section 1 are made from the same stainless steel.

The Wrench for Start-X, used to screw the insert in the handpiece, is made from the same stainless steel.

6 COMPATIBLE DEVICES

Start-X® are designed to be exclusively used with Piezo ultrasonic generators as follows:

- **SATELEC ultrasonic generator**, for type M3 x 0.6 mm threads (**REF. A0660, A9660**).
- **EMS ultrasonic generator**, for type M3 x 0.5 mm threads (**REF. A0661, A9661**).

Power setting shall be adjusted, depending on the specific clinical application and insert n°:

	REF. A0660, A9660		REF. A0661, A9661	
	SATELEC / P5 Newtron		EMS Minipiezon	
	Min	Max	Min	Max
Insert n°1	16	18	8	9
Insert n°2	16	18	8	9
Insert n°3	10	18	4	7
Insert n°4	18	20	9	10
Insert n°5	7	18	3	8

Please refer to the generators IFUs for their Instructions For Use.

7 SAFETY NOTES

7.1 Warnings

- Strictly follow the Instructions For Use and the General Processing Instructions For Endodontic Products (see section [9]) to minimize the following risks to the devices, the patient and/or the user:
 - Serious health consequences (infection, toxic or allergic reactions)
 - Soft tissue damage
 - Swallowing / inhalation of a part of the device
 - Tooth extraction
 - Tissue necrosis
 - Treatment delay
- The use of inappropriate power settings, pressure and contact duration would speed up the wear of Start-X® inserts and increase the risk of breakage and heat generation.

7.2 Precautions

- Safety and effectiveness of use have not been established in pregnant or breastfeeding women or in children.
- Use a rubber dam system during the endodontic procedure.
- For your own safety, wear personal protective equipment (gloves, glasses, mask).
- All inserts should be used with minimal pressure and a light brushing motion.
- Inspect the packaging before use and do not use the instruments if the packaging is damaged.
- Inspect the instrument before each use for signs of defects such as deformations, breakage, corrosion or loss of marking. With these indications the devices are not able to fulfil the intended use with the required safety level, instruments should be discarded.
- Do not start the generator with damaged ultrasonic inserts.

- Processing step is mandatory before use according to section 8.
- The insert must be tightened onto the ultrasonic generator handpiece with the enclosed wrench prior to use.
 - Be careful to screw the selected insert in the right axis with the right angle thanks to the wrench to prevent the insert from under tightening or thread deterioration. Make sure it is well connected before starting the generator.
 - Do not apply an excessive torque when screwing the insert thanks to the wrench. Pay attention to over tightening.
- Do not apply excessive pressure when using the device.
- Irrigate during the use of the device to prevent overheating.
- The use of visual aid is recommended for improved working comfort and safety.

7.3 Adverse reactions / Undesirable side effects

Temporary post-operative discomfort could arise generally as a result of receiving a root canal treatment.

7.4 Storage conditions / Handling

Keep products in a dry and clean environment, away from sources of moisture and direct sunlight.

8 STEP BY STEP INSTRUCTIONS

- 1) Before use, check if the instrument is sterile or if processing is needed (refer to Processing Instructions).
- 2) Adjust the power of the ultrasonic generator according to the selected ultrasonic insert thread and n° as described in section [2].
- 3) Use the wrench to screw the selected insert in the handpiece.
- 4) Start with lower power settings and only increase the power as required to accomplish the clinical task. All tips should be used with minimal pressure and a light brushing motion.
- 5) Irrigation is mandatory to avoid overheating both the insert and the tooth.

9 PROCESSING: HYGIENE, CLEANING, DISINFECTION AND STERILIZATION INSTRUCTIONS

For all products listed in this IFU, please follow the General Processing Instructions For Endodontic Products (dentsplysirona.com/ifu) for the disinfection, cleaning and sterilization steps.

10 SAFE DISPOSAL OF THE DEVICE

Products shall be disposed of according to local regulations for the safe disposal of contaminated devices.

11 EXPIRATION DATE

None.

12 SERIOUS INCIDENT REPORTING

Any serious incident in relation to the product should be reported to the manufacturer and the competent authority according to local regulations.

13 WARRANTY AND LIMITED REMEDY / LIMITATION OF LIABILITY

Disclaimer: Users are solely responsible for any deviation from these instructions, and/or the use of alternative methods for processing. Dentsply Sirona accepts no liability for damage, injury, or any legal responsibility incurred directly or indirectly by the user due to a deviation from the Instructions For Use set forth above. The user shall observe safe and lawful practices including, but not limited to, those set forth in this document.

14 ADDITIONAL INFORMATION

See webpage dentsplysirona.com/ifu to:

- Get a free print copy of the IFU.
- Get an explanation of non-harmonized symbols for IFUs and labels.

Authorised representatives

DENTSPLY DeTrey GmbH
De-Trey-Straße 1
78467 Konstanz
Germany

DENTSPLY IH Limited
Brunel Way
Stonehouse
Gloucestershire
GL10 3GB

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Dentsply Sirona Europe GmbH