Gnathometer M with UTS Adapter

EN Operating Instructions

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
ivoclar

**List of Parts**

- 1 Lower basic arch
- 2 Maxillary registration plate with UTS adapter
- 3 Fixing plate
- 4 Bite rim supports
- 5 Mounting jig GN
- 6 Registration pin
- 7 Mandibular registration plate
- 8 Maxillary registration plate
- 9 Fixation screw
- 10 Washer
- 11 Guide component
- 12 Funnel hole
- 13 Stylus
- 14 UTS adapter
- 15 Recess for the guide component
- 16 Retentions
- 17 Snap-action mechanism
- 18 Registration pin holder

Preface

Dear Customer,

Thank you for having purchased the Gnathometer M. Inappropriate use may damage the equipment and result in personal harm. Please observe the relevant notes and read the Operating Instructions. We hope you enjoy working with Gnathometer M.

1. Intended use**Intended purpose**

Intraoral registration

Use

For dental use only!

Description

Gnathometer M is an intraoral support pin registration device for the individual determination of the centric TMJ position in edentulous patients.

Indications

Determining the TMJ relations in edentulous patients.

Contraindications

If the device is used according to the instructions, no contraindications are known to date.

Side effects

There are no known side effects to date.

Composition

The following components of the product are used intraorally:

- Position 1, lower basic arch, stainless steel 1.4301, (X5CrNi18-10)*
- Position 2, upper registration plate with UTS adapter, stainless steel 1.4301, (X5CrNi18-10)*
- Position 3, fixing plate, PMMA
- Position 4, bite rim supports, PA6
- Position 6, registration pin, stainless steel 1.4301, (X5CrNi18-10)*
- Position 7, mandibular registration plate, stainless steel 1.4301, (X5CrNi18-10)*
- Position 8, maxillary registration plate, stainless steel 1.4301, (X5CrNi18-10)*
- Position 9, fixation screw, stainless steel 1.4301, (X5CrNi18-10)*
- Position 10, washer, stainless steel 1.4301, (X5CrNi18-10)*

*contains nickel

Signs and symbols

The signs and symbols in these Operating Instructions facilitate the finding of important points and have the following meanings:

Risks and dangers

Contraindication

Important information

Target group

Dentists, dental staff

2. Application

The basic principle of the intraoral support pin method is based on the idea of McGrane, Gysi, Gerber etc. of the threepoint rest of the two jaw joints and the registration pin, which is placed in the balance point of the mandible. The support pin draws the typical image of the Gothic Arch on a coloured registration plate if the patient performs movements towards the right lateral and the left lateral out of the dorsal position.

The tip of the Gothic Arch marks the retral contact position, while the centric occlusal position is marked by the crossing point. The desired centric position can then be fixed using the fixing plate. Subsequently, the entire registration is three-dimensionally locked with a suitable material (e.g. registration silicone or stone). The registration can then be passed on to the laboratory for further processing.

The special feature of this registration device is that the registration plates are equipped with a snap-action mechanism. With this snap-action mechanism, the bite rim supports and the registration plates can be quickly and easily attached and removed in order to conduct the different working steps. In a first step, the bite rim supports are used for the functional impressions. After that, the bite rim supports are removed and replaced with the registration plates. Then, needle point tracing is performed to determine the centric function. Also, Gnathometer M can be connected with the UTS transferbow by means of the UTS adapter and then used for the skull-related transfer of the models into the Stratos articulator.

2.1. Installation and initial start-up

2.1.1. Unpacking and checking the delivery contents

Please make sure that the contents are complete.

- 1 Upper registration plate with UTS adapter
- 1 Mandibular registration plate
- 2 Bite rim supports
- 1 Registration pin
- 1 Fixing plate
- 1 Fixation screw
- 1 Washer
- 2 Guide components
- 1 Mounting jig

2.1.2. Assembling the components

Important information

Clean the components before they are used for the first time (see Cleaning in Chapter 3). Gnathometer M is already partly assembled.

2.2. Handling and mounting

2.2.1. Mounting Gnathometer M in the laboratory

Align Gnathometer M with the help of the mounting jig as follows:

- dorsally on the retromolar pads
- ventrally with the front edge over the mucolabial fold (the result should be a vertical anterior arch)
- horizontally halfway between the alveolar ridges
- showing bilateral symmetry from the centre of the maxilla

The Gnathometer M can be moved in a sagittal direction on the mounting jig. In the mandible, the stylus should be located in the region of the premolars. The mounting jig is tapered in the dorsal region and can be adjusted to adapt to the retromolar pads. Please block out any undercuts on the model with wax.


Isolate models with Separating Fluid. Then fabricate the trays in the usual manner using SR Ivolen.

Apply SR Ivolen on the mandibular tray. Set Gnathometer M in the corresponding position into the still soft material and secure.

Next, the maxilla can be mounted. To do so, also apply SR Ivolen to the upper model and close the occludator or articulator. The Gnathometer M is now mounted in the correct position.

To facilitate finishing, remove the mounting jig and registration plate.

Important information
Ivoclar Vivadent recommends securing Gnathometer M directly on the tray with resin. This saves space and ensures stability during the clinical procedure. Before mounting, it is essential to check the fit of the registration plates in the arches.

2.3. Practical use on the patient

2.3.1. Taking the functional impression


Check the marginal design of the individual trays with bite rim supports (4) and adjust, if required. The bite rim supports (4) close flush with an even distribution of pressure on the denture bases.

Now, the functional impression can be taken. Whilst the impression is being taken, the antagonist tray remains in the mouth to provide the correct jaw relationship for the subsequent registration.

2.3.2. Mounting the registration components

Remove bite rim supports (4).

Replace the mandibular registration plate (7) and the maxillary registration plate (8). Use suitable pliers for this purpose, as can also be seen in the image, to secure the registration plate (7) in the snap-action mechanism (17) of the basic arch (1). This also ensures that the impression margins are not damaged when the plates are inserted using needle-nosed pliers.

Remove the fixing plate (3) by loosening the fixation screw (9) with a suitable screwdriver.

Warning
To prevent injury by the registration pin, the patient must be instructed to slide the tongue below the registration plate or towards the dorsal region.

Now reposition the registration components in the mouth. When the patient closes the mouth, only one point of contact will be on the registration plate.

2.3.3. Adjusting the occlusal height

Check vertical relation. If necessary, adjust the height by turning the registration pin (6) with a screwdriver. Fine adjustments may also be carried out with the help of the cross hole in the registration pin and a suitable instrument.


Important information
One full rotation results in a vertical change of 0.8 mm, which increases to 1mm in the anterior region.

The screw thread provides more free space for increasing the vertical dimensions than for decreasing the vertical dimensions. For this reason, Ivoclar Vivadent recommends taking the impression in a lower bite position and to then raise to the correct occlusal height with the registration. Please also take into consideration that there will be a minimal increase in occlusal height when the functional impression is taken. Once the height has been adjusted, check the distance at the dorsal ends of the two impression trays with a flat instrument. The trays or excess impression material must not touch each other. The freedom of movement of the mandible must not be compromised.

2.3.4. Conducting the needle point registration

For needle point tracing, remove the upper impression from the oral cavity, remove the registration plate (8), warm slightly and apply a thin layer of wax colour (crayon, marker). Then replace the registration plate (8).

Both functional impressions are inserted in the mouth.

Now instruct the patient to execute the following movements:

It is recommended to practice these movements beforehand!


- Move the mandible forward (protrusion) – and back,
- Repeat once
- Move the mandible to the right (lateral movement) – and back,
- Repeat once
- Move the mandible to the left (lateral movement) – and back,
- Repeat once


Important information

The patient must sit upright and hold the head straight (no leaning).

Remove the registration from the mouth and check the markings (Gothic Arch) on the upper plate (8). If the Gothic Arch is unclear, repeat the procedure.


To secure the registration and for the final check, attach the fixing plate (3) with the fixation screw (9) in one of the two openings over the centre of the needle point tracing using a screwdriver.

Watch whilst the patient repeatedly closes his mouth to check whether the stylus (13) slides into the corresponding funnel hole (12) in the fixing plate (3). If the stylus (13) lands next to the funnel hole (12) again and again, adjust the fixing plate (3) to this point or repeat the registration procedure.

2.3.5. Verifying and securing the centric position

Confirm the correct position by repeating the closing movements.

Secure the registration by injecting a fast-setting impression plaster or registration silicone (CADBite), using a plastic syringe, into the space between the arches. In doing so, the registration pin (6) must be located in the lower fixing plate hole (3). Only the lips should be open.

If no extra-oral registration using a transferbow is made, an impression of the lips from the inside can be made in the anterior region using a PVS material or plaster.


In prosthetics, the individual transfer of the Bonwill triangle with the UTS 3D transferbow is preferably made according to Camper's plane. The individual transfer of the Bonwill triangle with the transferbow helps to prevent occlusal interferences by the exact orientation of the models with the articulator joints.

When securing the intraoral registration, ensure that the UTS adapter remains clean. If necessary, transfer the transferbow before securing. Mount the registration joint onto the UTS adapter up to the stop. Once it has been removed, this is the only way that the registration joint can be re-allocated correctly in the laboratory.

2.3.6. Skull-related transfer by means of the UTS 3D transferbow

Mount the universal transferbow on the patient's head according to the Operating Instructions. The connection to the Gnathometer M is established with the UTS registration joint, which can now also be attached to the Gnathometer M.

2.3.7. Disinfecting the impression and registration components

Before passing on the components, they must be cleaned and disinfected using suitable disinfectants. In this context, please also observe the general hygiene guidelines and stipulations, as well as Chapter 3.

3. Maintenance, cleaning, diagnosis

This chapter describes the user maintenance and cleaning procedures for Gnathometer M. Only those tasks are listed that may be performed by dental professionals and dental staff. All other tasks must be performed by qualified service personnel at a certified Ivoclar Vivadent Service Centre.

Monitoring and maintenance

The time for these maintenance procedures depends on the frequency of use and the working habits of the user. For that reason, the recommended times are only approximations.

Cleaning

Item	Frequency	Cleaning material
Registration components and small parts (metal)	After every use	Sterilization in the autoclave at 134 °C for 5.5 min using the prevacuum method. Please observe the valid standards and regulations.
Fixing plates (plastic)	After every use	Disinfection in the immersion bath according to the instructions of the disinfectant manufacturer
Bite rim supports (plastic)	After every use	Disinfection in the immersion bath according to the instructions of the disinfectant manufacturer

Maintenance

Avoid any contact with strong acids and solvents (e.g. MMA, acetone) to prevent damage of the surfaces. Clean, sterilize and/or disinfect the components before the first use and after every use and/or after contact with saliva or blood.

- Use only disinfectants suitable for stainless steel and light metals.
- Avoid any contact with strong acids, brines or solvents.
- Clean in an ultrasonic bath, with water, steam jets or sterilization devices.

4. What if...

This chapter will help you to recognize malfunctions and take appropriate measures or, if possible, to perform some minor repairs.

The expendable parts are as follows:

- Fixing plates
- Bite rim supports

Technical malfunction

Description	Instructions for users	Measures
The surface of the device is scratched and damaged	Was the surface cleaned or manipulated with a sharp object?	Do not scratch the surface with a sharp object
The surface of the device is damaged or partially corroded	Was the surface cleaned with acid or a solvent?	Do not use any acids or solvents for cleaning
Fixing plate is corroded and destroyed	Was the surface cleaned with an unsuitable disinfectant?	Do not use any unsuitable disinfectants. Order a new part (expendable part)
Bite rim support is corroded and destroyed	Was the surface cleaned with an unsuitable disinfectant?	Do not use any unsuitable disinfectants. Order a new part (expendable part)

Repair


Repairs may only be carried out by a certified Ivoclar Vivadent Service Centre. Please refer to the addresses of the Service Centres on the back of these Operating Instructions.

If repairs during the warranty period are not carried out by a certified Service Centre, the warranty will expire immediately. Please also refer to the corresponding warranty regulations.


Setting the clamping effect

After repeated use, the clamping effect between the basic arch and the registration plate may diminish. This can be easily corrected:

The clamping effect is only between points A and B/C.

Points E and F always show a tolerance and are not involved in the holding effect. The clamping effect can be increased by slightly tapping the centre of the basic arch. The distance D can be reduced by fractions of a millimetre.

If the clamping effect is too strong, reduce the distance D by fractions of a millimetre as shown in the picture. Make sure that the basic arch remains flat.

5. Safety information

In the case of serious incidents related to the product, please contact Ivoclar Vivadent AG, Bendererstrasse 2, 9494 Schaan/Liechtenstein, website: www.ivoclar.com and your responsible competent authority. The current Operating Instructions are available in the download section of the Ivoclar Vivadent AG website (www.ivoclar.com).

Disposal information

The product must be disposed of according to the corresponding national legal requirements.

Risks and dangers

Risk of swallowing and choking hazard

Since there are small components that may fall off if not attached correctly, there is a risk of swallowing these parts.

Risk of injury

To prevent injury by the registration pin, the patient must be instructed to slide the tongue below the registration plate.

Gnathometer M must only be used for the purpose described in Chapter 2. Further instructions to assure proper use of Gnathometer M:

- The instructions, regulations, and notes in these Operating Instructions must be observed.
- Gnathometer M must be properly maintained (see Chapter 3).

6. Shelf life and storage

This product does not require any special storage conditions.

7. Product specifications

Technical data

Dimensions:

Mandibular registration plate (into which the stylus is screwed)

Width: 24.6 mm

Length: 54.4 mm

Material thickness: 2 mm

Weight: 9.4 g

Maxillary registration plate (to which the fixing plate is attached)

Width: 36.9 mm

Length: 54.4 mm

Material thickness: 2 mm

Weight: 12.2 g

Bite rim support

Width: 32.5 mm

Length: 54.5 mm

Material thickness: 2 mm

Weight: 1.5 g

Material

Registration plates: Stainless steel

Bite rim supports: Polyamide

Fixing plates: PMMA, crystal clear

Height adjustment of the stylus: 5 mm

Tips for Gothic Arch analysis and interpretation


Please also refer to the specialized literature available on this subject matter.

8. Additional information

Keep out of the reach of children.

The device has been developed solely for use in dentistry. Start-up and operation should be carried out strictly according to the Operating Instructions. Liability cannot be accepted for damages resulting from misuse or failure to observe the Instructions or the stipulated area of use. The user is responsible for testing the apparatus for its suitability and use for any purpose not explicitly stated in the Instructions.

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