SR Ivocron®

[en] Instructions for Use

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1 Intended use

Intended purpose

Fabrication of temporary and fixed crowns and bridges by means of veneering metal alloy frameworks as well as without frameworks. Modification of resin teeth and denture base materials.

Patient target group

- Patients with permanent teeth
- Adult patients with dental implants
- Edentulous adult patients

Intended users

- Dental laboratory technicians (fabrication of restorations in the dental laboratory)
- Dentists (clinical procedure)

Special training

No special training required.

Use

For dental use only.

Description

SR Ivocron is a high-quality PMMA veneering material for the crown and bridge technique as well as for the fabrication of temporary restorations. As the product accommodates a wide range of work procedures, it is suitable for the fabrication of resin veneers (hot/press technique) and short-term, medium-term and long-term temporaries as well as for securing denture teeth on model-cast frameworks (cold technique).

The PMMA veneering material SR Ivocron allows all 20 Chromascope shades and Gingiva pink to be reproduced.

Product name	Product description
SR Ivocron Opaquer	Masking of metal frameworks and adjustment of the base shade
SR Ivocron Dentin	Veneering of metal frameworks; fabrication of single crowns and bridges; modification of denture teeth
SR Ivocron Incisal	Veneering of the occlusal and incisal areas of metal frameworks; fabrication of single crowns and bridges; modification of denture teeth
SR Ivocron Gingiva	Veneering of the gingival portion of metal frameworks; fabrication of single crowns and bridges
SR Ivocron Intensiv	Intensive powder shades for adjusting the shade of restorations made of SR Ivocron. They allow the reproduction of natural characteristics.
SR Ivocron Opaquer Liquid	Mixing liquid especially coordinated with SR Ivocron Opaquer
SR Ivocron Cold Liquid	Mixing liquid for the cold technique (pouring technique)
SR Ivocron Hot Liquid	Mixing liquid for the hot technique (layering technique)
SR Ivocron Press Liquid	Mixing liquid for the press technique (flasking technique)

Technical data

SR Ivocron is a PMMA-based crown and bridge veneering material of type 2 and class 1 (ISO 10477:2020) featuring the following mechanical properties:

Properties	Specification ¹
Flexural strength	> 50 MPa
Water absorption	≤ 40 µg/mm ³
Solubility	≤ 75 µg/mm ³

¹according to ISO 10477:2020

Indications

- Missing tooth structure in anterior and posterior teeth
- Partial edentulism in the anterior and posterior region
- Complete edentulism

Areas of application:

Cold technique (pouring technique)

- For short-term and medium-term temporary restorations
- Securing of ground denture teeth on the model-cast framework
- Repair of denture teeth
- Masking of metal frameworks and adjustment of the base shade

Hot technique (layering technique) press technique (flasking technique)

- Crown and bridge veneering technique
- Masking of metal frameworks and adjustment of the base shade

Shade modification and characterization

- Characterization and shade modifications of dentin, incisal or gingiva materials

Contraindications

The use of this product is contraindicated if the patient is known to be allergic to any of its ingredients.

Limitations of use

- The final restoration must not be reused.
- Do not apply directly in the oral cavity.

Side effects

There are no known side effects to date.

Interactions

There are no known interactions to date.

Clinical benefit

- Reconstruction of chewing function
- Restoration of esthetics

Composition

- **SR Ivocron powders (dentin, incisal, gingiva, intensive)**
PMMA, dibenzoyl peroxide
Total content of inorganic fillers: 0.02–0.6 wt%
Particle size of the inorganic fillers: 0.02–0.8 µm
- **SR Ivocron powder (opaquer)**
Copolymer, TiO₂, PMMA, iron oxide, dibenzoyl peroxide
Total content of inorganic fillers: 58–62 wt%
Particle size of the inorganic fillers: 0.05–0.8 µm
- **SR Ivocron Cold Liquid**
MMA, EGDMA, TEGDMA
- **SR Ivocron Hot Liquid**
MMA, TEGDMA
- **SR Ivocron Press Liquid**
MMA, EGDMA
- **SR Ivocron Opaquer Liquid**
MMA

2 Application

Notice:

- Only use the system components mentioned in chapter "Intended use".
- In order to apply shade adjustments or characterizations, mix incisal, dentin and gingiva materials with intensive material.
- Use SR Link as the bonding agent.
- SR Ivocron Opaquer is self-curing.

2.1 Pouring technique (cold technique)

Fabricating temporary restorations

Fabricating a silicone key

1. Fabricate the working model.
2. Wax up the restoration.
3. Check the shape and function of the wax-up.
4. Fabricate the silicone key.

Boiling out and isolating

1. Boil out the wax.
2. Block out undercuts on the model. Use pink modelling wax.
3. Isolate the model using separating fluid.
4. Allow the model to dry for 1 min.
5. Repeat steps 3 to 4.

Pouring and polymerizing the dentin and incisal materials

Notice:

- Observe the following processing parameters:

Mixing ratio		Dough time	Working time at 23 °C	Polymerization in the pressure apparatus
by volume	in g			
1 part polymer : 1 part monomer	1 g polymer : 0.83 g monomer	3–4 min	approx. 8 min	2–6 bar 15 min 40–50°C

1. Mix the material in a rubber cup according to the processing parameters indicated.
2. Cover and let it sit according to the processing parameters indicated.
3. Pour the prepared incisal material into the silicone key.
4. Let it sit.
5. Fill up the silicone key with prepared dentin material.
6. Reposition the filled silicone key on the isolated model.
7. Immobilize the silicone key and the model with a rubber band.
8. Place the immobilized model in the pressure apparatus.
9. Polymerize according to the processing parameters indicated.

Finishing and polishing

1. Finish with burs and rubber polishers.
2. Prepolish the restorations with goat's hair brushes and pumice.
3. Polish the restorations to a high gloss with cotton wool buffers and high-gloss polishing medium.

Filling and securing ground denture teeth on model-cast frameworks

Fabricating a silicone key

1. Grind the denture teeth to fit on the working model.
2. Secure them in place with wax.
3. Check occlusion.
4. Fabricate the silicone key.

Boiling out and isolating

1. Boil out the wax.
2. Block out undercuts on the model. Use pink modelling wax.
3. Isolate the model using separating fluid.
4. Allow the model to dry for 1 min.
5. Repeat steps 3 to 4.

Applying monomer

1. Roughen the basal areas of the teeth.
2. Replace the roughened teeth in the silicone key.
3. Wet the teeth with monomer.

Conditioning the framework

1. Blast the framework with Al_2O_3 (80–100 μm) at 2–3 bar pressure.
2. Tap to remove grinding dust.
3. Optional: Wet the framework with bonding agent.
4. Allow the framework to dry for 3 min.

Applying the opaquer and allowing it to polymerize

1. Mix opaquer with opaquer mixing liquid.
2. Cover the mixture and let it sit for 2–3 min.
3. Apply the opaquer.
4. Allow to polymerize for 15 min.
5. Before proceeding, check whether the opaquer has cured completely.

Pouring and polymerizing the dentin and incisal materials

Notice:

- In order to make shade adjustments in interdental areas, mix incisal and dentin material in a ratio of 1:1.
- Observe the following processing parameters:

Mixing ratio		Dough time	Working time at 23 °C	Polymerization in the pressure apparatus
by volume	in g			
1 part polymer : 1 part monomer	1 g polymer : 0.83 g monomer	3–4 min	approx. 8 min	2–6 bar 15 min 40–50°C

1. Mix the material in a rubber cup according to the processing parameters indicated.
2. Cover and let it sit according to the processing parameters indicated.
3. Pour the prepared material into the silicone key.
4. Reposition the filled silicone key on the isolated model.
5. Immobilize the silicone key and the model with a rubber band.
6. Place the immobilized model in the pressure apparatus.
7. Polymerize according to the processing parameters indicated.

Finishing and polishing

1. Finish with burs and rubber polishers.
2. Prepolish the restorations with goat's hair brushes and pumice.
3. Polish the restorations to a high gloss with cotton wool buffers and high-gloss polishing medium.

2.2 Layering technique (hot technique)

Designing the framework

1. Use metal to fabricate the palatal, incisal and occlusal areas of the bridge.
2. Add mechanical retentions to the surface to be veneered.

Blocking out and isolating

1. Immerse the model in water for 5 min.
2. Block out the undercuts on the model. Use pink modelling wax.
3. Isolate the model using separating fluid.
4. Allow the model to dry for 1 min.
5. Repeat steps 3 to 4.

Conditioning the framework

1. Blast the framework with Al_2O_3 (80–100 μm) at 2–3 bar pressure.
2. Tap to remove grinding dust.
3. Optional: Wet the framework with bonding agent.
4. Allow the framework to dry for 3 min.

Applying the opaquer and allowing it to polymerize

1. Mix opaquer with opaquer mixing liquid.
2. Cover the mixture and let it sit for 2–3 min.
3. Apply the opaquer.
4. Allow to polymerize for 15 min.
5. Before proceeding, check whether the opaquer has cured completely.

Layering and polymerizing the incisal and dentin materials

Notice:

- Observe the following processing parameters:

Mixing ratio		Dough time	Working time at 23 °C	Polymerization in the pressure apparatus	
by volume	in g			with model	without model
1–1.5 parts polymer : 1.5 parts monomer	1–1.5 g polymer : 1.25 g monomer	2–3 min	8–25 min	2–6 bar 25–30 min 100 °C	2–6 bar 10 min 120 °C

1. Mix the materials in a rubber cup according to the processing parameters indicated.
2. Cover and let them sit according to the processing parameters indicated.
3. Layer the prepared materials.
4. Polymerize according to the processing parameters indicated.

Finishing and polishing

1. Finish with burs and rubber polishers.
2. Prepolish the restorations with goat's hair brushes and pumice.
3. Polish the restorations to a high gloss with cotton wool buffers and high-gloss polishing medium.

2.3 Flasking technique (press technique)

Designing the framework

1. Use metal to fabricate the palatal, incisal and occlusal areas of the bridge.
2. Add mechanical retentions to the surface to be veneered.

Conditioning the framework

1. Blast the framework with Al_2O_3 (80–100 μm) at 2–3 bar pressure.
2. Tap to remove grinding dust.
3. Optional: Wet the framework with bonding agent.
4. Allow the framework to dry for 3 min.

Applying the opaquer and allowing it to polymerize


1. Mix opaquer with opaquer mixing liquid.
2. Cover the mixture and let it sit for 2–3 min.
3. Apply the opaquer.
4. Allow to polymerize for 15 min.
5. Before proceeding, check whether the opaquer has cured completely.

Modelling and investing

1. Build up the veneer using tooth-coloured wax.
2. Invest the waxed-up restoration in the flask using dental stone.

Notice:

- Place the restoration at a 25–30° angle. Make sure that the restoration offers proximal access.

Boiling out and isolating

1. Boil out the wax.
2. Isolate the warm plaster halves with separating fluid.
3. Allow to dry for 1 min.
4. Repeat steps 2 to 3.

Applying monomer

Wet the exposed opaquer with monomer.

Pressing the dentin material

Notice:

- Mix the dentin material to a viscous consistency.
- Observe the following processing parameters:

Mixing ratio		Dough time	Working time at 23 °C	Polymerization in the flask in the clamping frame	
by volume	in g			Single crown	Bridge
2.5 parts polymer : 1 part monomer	2.5 g polymer : 0.83 g monomer	2–3 min	8–10 min	1. Place in cold water. 2. Boil for 30 min. 3. Allow to cool slowly.	1. Preheat in water measuring 70 °C for 30 minutes. 2. Boil for 30 min. 3. Allow to cool slowly.

1. Mix the material in a rubber cup according to the processing parameters indicated.
2. Cover and let it sit according to the processing parameters indicated.
3. Preheat the upper flask half in boiling water.

Notice:

- Prevent any contact of the isolated surface with water.
4. Apply the prepared material.
 5. Cover with a plastic film.
 6. Place the preheated upper flask half on the lower half.
 7. Apply 2 t of pressure to the flask.
 8. Allow to cool under pressure.

Reducing the incisal area

1. Open the flask
2. Remove the plastic film.
3. Bevel the dental material with a sharp instrument.

Pressing and polymerizing the incisal material

Notice:

- Mix the incisal material to a viscous consistency.
- Observe the following processing parameters:

Mixing ratio		Dough time	Working time at 23 °C	Polymerization in the flask in the clamping frame	
by volume	in g			Single crown	Bridge
2 parts polymer : 1 part monomer	2 g polymer : 0.83 g monomer	3–4 min	8–10 min	1. Place in cold water. 2. Boil for 30 min. 3. Allow to cool slowly.	1. Preheat in water measuring 70 °C for 30 minutes. 2. Boil for 30 min. 3. Allow to cool slowly.

1. Mix the material in a rubber cup according to the processing parameters indicated.
2. Cover and let it sit according to the processing parameters indicated.
3. Preheat upper flask half in boiling water.

Notice:

- Prevent any contact of the isolated surface with water.
4. Apply the prepared material.
 5. Cover with a plastic film.
 6. Place the preheated upper flask half on the lower half.
 7. Apply 1.5 t of pressure to the flask.
 8. Allow to cool under pressure.
 9. Polymerize according to the processing parameters indicated.

Finishing and polishing

1. Finish with burs and rubber polishers.
2. Prepolish the restorations with goat's hair brushes and pumice.
3. Polish the restorations to a high gloss with cotton wool buffers and high-gloss polishing medium.

Seating the restoration

Notice:

- Seating is exclusively done by the dentist.
- Observe the instructions for use of the framework material when seating framework-supported restorations.

Depending on the type of restoration, seat it using either a temporary, conventional, adhesive or self-adhesive luting procedure.

3 Safety information

- In the case of serious incidents related to the product, please contact Ivoclar Vivadent AG, Bendererstrasse 2, 9494 Schaan/Liechtenstein, www.ivoclar.com and your responsible competent authority.
- The current Instructions for Use and the explanation of symbols are available on the website: www.ivoclar.com/elifu
- The Summary of Safety and Clinical Performance (SSCP) can be retrieved from the European Database on Medical Devices (EUDAMED) at <https://ec.europa.eu/tools/eudamed>.
- Basic UDI-DI: 76152082ACOMP002LE

Warnings

- Observe the Safety Data Sheet (SDS) (available at Ivoclar Vivadent AG www.ivoclar.com).
- Observe the safety notes on the individual primary packages and labels.
- SR Ivocron Opaquer Liquid, as well as the Cold, Hot, and Press Liquids contain methyl methacrylate. MMA is irritant and highly flammable (flash point: 10 °C/50 °F).
- Do not inhale vapours.
- Contact with unpolymerized material may have a slight irritating effect and may lead to a sensitization against methacrylates. The material irritates eyes, respiratory organs, and skin.
- Commercial medical gloves do not provide protection against the sensitizing effect of methacrylates.
- Do not inhale grinding dust.

Disposal information

Remaining stocks or removed restorations must be disposed of according to the corresponding national legal requirements.

Residual risks

Users should be aware that any dental intervention in the oral cavity involves certain risks.

The following known residual clinical risks exists:

- Chipping, fracture or loss of retention of the crown and veneering material may lead to accidental swallowing or inhalation and dental re-treatment.

4 Shelf life and storage

- Storage temperature for powders and mixing liquids: 2–28°C (36–82 °F).
- Keep out of sunlight.
- Expiry date: see note on bottles and packaging.
- Do not use the product after the indicated expiry date.
- Check to make sure that the packaging and the product are intact and undamaged before use. If in doubt, contact Ivoclar Vivadent AG or your local sales partner.

5 Additional information

Keep material out of the reach of children!

The material has been developed solely for use in dentistry. Processing should be carried out strictly according to the Instructions for Use. Liability cannot be accepted for damages resulting from failure to observe the Instructions or the stipulated area of use. The user is responsible for testing the product for its suitability and use for any purpose not explicitly stated in the Instructions.

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- Perskaitykite saugos duomenų lapą (SDL) (jį rasite Ivoclar Vivadent AG www.ivoclar.com).