WOHLWEND AG

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Introduction to

VISION-CLASSIC 'The metal-ceramic'

Our Vision-Classic is a high-fusing micro structured leucite glass-ceramic. It makes the veneering of metal frameworks a safe and simple job. It is characterized by its enamel-like light refraction and reflection behavior. Very individual and aesthetic results can be achieved by using fluorescent and opalescent additional powders.

Wohlwend's unique manufacturing process produces an ideal grain size and pigment distribution. The advantages of this unique process include.

- Consistently accurate color rendering-V-colors A1-D4
- Easy handling
- High modelling stability, accurate reproduction after firing
- Low shrinkage during firing
- No green discoloration in combination with silver-containing alloys

Vision-Classic chemical composition

Material designation: Amorphous leucite glass-ceramic

Essential components firmly integrated into the glass-ceramic network:

Vision-Classic: Opaker, Dentin, Incisal, Transparent, Shoulder, Gingiva, Mamelons

SiO_2 , Al_2O_3 , K_2O , Na_2O , CaO ,

Physical properties

Vision-Classic	Flexural strength (Mpa)	Solubility ($\mu\text{g}/\text{cm}^2$) 100 $\mu\text{g}/\text{cm}^2$ [± 5]	Thermal expansion CTE (25/500°C) [ppm/K]	Temperature of the glass transition TG [°C] ± 20
Powder Opaker	100	50	13.8	600
Paste Opaker	100	50	13.8	630
Intensiv Opaker	100	50	13.8	600
Schulter	80	25	13.8	590
Effect	80	25	13.3	590
Dentin	80	25	13.3	590
Incisal	80	25	13.3	590
Transparent	75	25	13.3	590
Opal-Incisal	80	25	13.3	590
Clear	80	25	13.3	590
Mamelon	80	25	13.3	590
Gingiva	80	25	13.3	590

Classification

- Dental ceramics type I/class 1
- The tested materials comply with the EN ISO 6872 standard

Indications

- Veneering of dental alloys for the fabrication of fixed dentures, for the metal-ceramic technique with a CTE of 13.8–15.2 10^{-6} - K-1 $\mu\text{m}/\text{mK}$ (25-600°C) on CoCr or NiCr alloys, precious metal alloys, cast, milled or laser-sintered.

Contraindication

- Combination with ceramic materials that do not belong to the products specified in Vision Ceramics.
- Use of non-approved framework materials.
- If the permissible framework or ceramic thickness is exceeded or undercut.
- Dental ceramic restorations are not recommended for patients with parafunction or bruxism.

Intended use

- VISION products are ceramic materials for dental treatments.

Patient target group

- No restrictions

Intended users

- Exclusively professional users. Dentist and dental technician (Rx only)
- The material is intended for professional use in dental laboratories and should only be used by trained personnel.

Reference to risks

- Serious incidents must be reported to the manufacturer and the competent authority
- The inhalation of ceramic dust released during processing must be avoided by means of suitable extraction devices and, if necessary, mouth and nose protection.

Storage

- Store the product in a dry place at 4–28°C (39–82°F)

Important!

Observe the alloy information!

- CTE of the alloy should be between 13.8–15.2 (25–600°C)
- The ideal value of the alloy is 14.1 to 14.2 (25–600°C)
- If the CTE value is higher than 14.5 (25–600°C), long-term cooling is recommended.
- If the CTE is lower than 14.1 (25–600°C), it is recommended that the object be fired should be removed from the firing table immediately after opening the firing chamber.
- For non-precious metal alloys, only new metal may be used for casting the frameworks! **The use of scrap metal can lead to a change in the CTE value of the alloy and there is a risk of contamination of the alloy with foreign substances!**
- We recommend the use of Vision-Bonding
- Please follow the instructions of your alloy manufacturer!

Color	A1	A2	A3	A3.5	A4	B1	B2	B3	B4	C1	C2	C3	C4	D2	D3	D4
Opaquere		O	O	O	O	O	O	O	O	O	O	O	O	O	O	O
Opakdentin		OD	OD	OD	OD	OD	OD	OD	OD	OD	OD	OD	OD	OD	OD	OD
Dentin		D	D	D	D	D	D	D	D	D	D	D	D	D	D	D
Incisal		57	58	59	59	60	57	59	59	59	60	59	59	60	60	59
Transpa materials		superclear, clear, neutral, opal, red, blue, grey														







Modified Shade	A1	A2	A3	A3.5	A4	B1	B2	B3	B4	C1	C2	C3	C4	D2	D3	D4	
Shoulder material		neutral + BA				neutral + BB				neutral + BC				neutral + BD			
Effect shoulder		ivory, straw, yellow-gold, brown, yellow-grey, brown-grey															
Opaque to intensify chroma		10%	20%	30%	35%	40%	10%	20%	30%	40%	10%	20%	30%	40%	20%	30%	40%
Intensive opaques		white, yellow, orange, ocher, brown, violet, gingiva															
Body to intensify chroma		10%	20%	30%	35%	40%	10%	20%	30%	40%	10%	20%	30%	40%	20%	30%	40%
Intensive dentine		ivory, straw, honeyyellow, brown, grey, blue, neon															
Intensive incisal		7, 8, 9, 10															
Mamelon		ivory, honey, sunset															
Effect masses		ivory, straw, honeyyellow, khaki															
Bleach		Opaque AB-0, Opaque-Dentin A0,B0, Dentine A0, B0, Incisal AB-0															
Cervical		straw, brown															
Gingiva		softpink, darkpink, softred, darkred															
Repair masses		Opaque, Dentine, Incisal, Transpa, Stain															

Metal processing prior to porcelain application

1. Prepare the alloys according to the manufacturer's instructions and round off sharp edges. Unless otherwise indicated, the framework should be sandblasted with min. 110 µm aluminium oxide and 2 bar pressure. Make sure that the colour of the oxide is uniform. For non precarious alloys the oxide layer should again be radiated with min. 110 µm Aluminium oxide (follow the alloy manufacturer's instructions!). Clean surface and degrease with steamer or ultrasonic bath.
2. Apply Vision-Bonding Paste evenly and thinly. The following is fired according to bonder instructions.
3. The result after firing also replaces a wash firing.
4. Alternatively Vision-Classic Opaquer can be used for the wash firing. Apply the opaquer thinly and evenly, see picture number 6.

	Start temperature	Dry time	Heat rate	Final temperature	Hold time	Vacuum
Degassing / Oxydbrand	Follow alloy manufacturer manual					
Vision-Bonding	550°C / 1022°F	6 min.	80°C / 176°F min.	980°C / 1796°F	1 min.	Yes
Vision-Classic / Powder Opaquer 1	600°C / 1112°F	3 min	80°C / 176°F min.	970°C / 1778°F	1 min.	Yes
Vision-Classic / Powder Opaquer 2	600°C / 1112°F	3 min	80°C / 176°F min.	950°C / 1742°F	1 min.	Yes
Vis-Opak / Paste Opaquer 1	600°C / 1112°F	6 min	80°C / 176°F min.	970°C / 1778°F	1 min.	Yes
Vis-Opak / Paste Opaquer 2	600°C / 1112°F	4 min	80°C / 176°F min.	950°C / 1742°F	1 min.	Yes

1st. Opaque:

Opaque powder & liquid

Mix the opaque to a creamy consistency with a brush or glass instrument than apply with 70% opacity. Firing temperature 970°C. With the use of OL-2 Vision-Classic Opaque Liquid a similar viscosity as with Vis-Opaque Paste can be achieved.

Vis Opaque / Pastes Opaque

Apply with a brush or glass instrument. Firing temperature 970°C.

If necessary, the temperature can be reduced to the oxide firing temperature of the alloy. Holding time for 1 minute. Close tightly after use to prevent dehydration. Glossy surface after firing.

Second step: 2nd opaque.

1. Apply the opaque even and covering.
2. Slightly glossy surface after firing.

Options

To achieve effects characterized from the depth with the help of opaque modifiers „IO“

For example:

Place IO-Orange in the central fossa or in the cervical area with IO-Brown.

In order to reach a higher chroma from the depth, Chroma Opaque is applied to the equator or mixed according to colour chart.

Third step: Ceramic shoulder / margin

Our shoulder masses in the basic colours A, B, C, D and Neutral will give you the possibility to create a shoulder for every colour situation. The mixing ratio can be taken from the table. Optionally, an even larger colour palette can be achieved with effect shoulder masses.

Shoulder / margin combination chart:


Vision – Classic Combination shoulder/margin chart																		
Shade	AO	BO	A1	A2	A3	A3.5	A4	B1	B2	B3	B4	C1	C2	C3	C4	D2	D3	D4
Shoulder / Margin	Ivory + NT		Ivory + BA					Ivory + BB				Ivory + BC				Ivory + BD		
Mixing ratio in %	20	80	80	70	50	40	20	80	70	50	30	80	70	50	30	70	50	20
			20	30	50	60	80	20	30	50	70	20	30	50	70	30	50	80
Effect shoulder	ivory, straw, yellow-gold, brown, yellow-grey, brown-grey																	

	Start temperature	Dry time	Heat rate	Final temperature	Hold time	Vacuum
Shoulder / Margin 1/2	450°C / 842°F	3 min.	45°C / 113°F min.	940°C / 1724°F	1 min.	Yes

Fourth step: layering


- 1. By limited space apply a thin layer of opaque-dentin (occlusal, cervical, etc.).
 - 2. Apply dentin (of the given shade) in full anatomical shape.
 - 3. Create a small cut back in the incisal and approximal region. Build up the marginal ridges and mamelons with Transpa Neutral and mamelon masses (ivory, honey, sunset).
 - 4. Complete the tooth shape of the tooth with the corresponding incisal mass. A slight over dimensioning during the construction is ensured by the light Shrinkage during firing for a perfect result.
- » The interdental space is filled with opaque-dentine.

	Start temperature	Dry time	Heat rate	Final temperature	Hold time	Vacuum
Dentin 1	450°C/842°F	6 min	45°C/113°F min.	920°C/1688°F	1 min.	Yes
Dentin 2	450°C/842°F	6 min	45°C/113°F min.	910°C/1670°F	1 min.	Yes

Important:

If Vision-Classic ML modelling liquid is used, the starting temperature of the ceramic furnace must be observed. This should not exceed 450°C. The constituents of the modelling liquid need time to evaporate. If this time is not granted due to too high starting temperature and too low pre-drying time, the dry outer layer of the ceramic may cause it to an inner build up of steam, which can cause damage to the surface of your layered crown during the firing circle!





Fifth step: Correction


1. After the first Dentin firing.
2. If necessary, add shape corrections with Transpa materials.
3. Second dentin firing
4. After finishing (picture 4) the contact points, tooth shape and surface –macro and micro structures, proceed with staining and glaze firing.
There are three glazing options available:
 - Glaze firing without glaze material and a subsequent mechanical glaze firing.
 - Polishing (rubber polishers, brushes, diamond polishing paste and so on).
 - Glaze firing with glaze material. Mix these with glaze liquid and apply thinly. If necessary, characterize with Visio-Classic Stains.
 - Glaze firing with Vision-Universal glaze and Stains.


Glaze Firing:

	Start temperature	Dry time	Heat rate	Final temperature	Hold time	Vacuum
Glaze without Glaze material	600°C/ 1112°F	2 min	60°C/ 140°F min.	920°C/ 1688°F	1 min.	No
Glaze with Glaze materiall	600°C/ 1112°F	3 min	60°C/ 140°F min.	910°C/ 1670°F	1 min.	No
Vision-Universal Glaze	480°C/ 869°F	3 min	45°C/ 113°F min.	890°C/ 1634°F	1 min.	No

Layering Guide

- **Neck**
Straw / Brown / or
D.BA / BB / BC / BD
- Dentin
- TR-NT
- MM 1 = Ivory
MM 2 = Honey
MM 3 = Sun-Set
- TR-SCL
- TR-BL
- II-10
- Incisal / Enamel
- II-7
- TR-NT


- TR-SCL
- MM 1 = Ivory
MM 2 = Honey
MM 3 = Sun-Set
- Dentin
- **Neck**
Straw / Brown / or
D.BA / BB / BC / BD
- TR-NT
- ID.Honey

- TR-NT
- **Incisal edge**
- II-7
- Dentin
- **Light blocker**
OD.White & Dentin (1:2)
- TR-SCL
- **Mamelon**
MM 1 = Ivory
MM 2 = Honey
MM 3 = Sun-Set
- Incisal / Enamel
- II-10
- ID.Honey
- **Neck**
Straw / Brown / or
D.BA / BB / BC / BD

- Incisal
- Dentin
- Opak Dentin
- Opaker

Classic Firing chart

Guidelines! Deviations are possible due to different furnace performances.

	Start temp.	Dry time	Heat rate	Final temperature	Hold	Vacuum	Appearance
Degassing	Follow alloy manufacturers manual						
Opaque 1*	1112°F 600°C	3 min.	176°F/min. 80°C	1778°F 970°C	1 min.	yes	Shiny
Opaque 2*	1112°F 600°C	3 min.	176°F/min. 80°C	1742°F 950°C	1 min.	yes	Shiny
Shoulder 1/2	842°F 450°C	3 min.	113°F/min. 45°C	1724°F 940°C	1 min.	yes	Shiny
Dentine 1	842°F 450°C	6 min.	113°F/min. 45°C	1688°F 920°C	1 min.	yes	Shiny
Dentine 2	842°F 600°C	6 min.	113°F/min. 45°C	1670°F 910°C	1 min.	yes	Shiny
Glaze	1112°F 600°C	2 min.	140°F/min. 60°C	1670°F 910°C	1 min.	–	Shiny

* With Vis-Opak paste, extend the drying time by 4 minutes to 6–7 pre drying time.
 Observe alloy instructions! For horizontal furnace systems (e.g. Programat®)
 Reduce start temperature to 400°C / 752°F
 For REPAIRKIT please use separate firing chart!

Mistakes	Reason	What to do
Bubbles, splits in the opaquer	Contamination in the metal; Wrong drytime; trapped air; paste too thick	Process alloy according to instruction; creamy consistence; Lengthen up drytime
Cracks: incisal, pontic	Framework design wrong; Check CTE-range;	Run a slow cooling cycle
Horizontal cracks	Slow cooling; Temper	Temper
Cracks; Tensile stress	Framework too thin; Check CTE-range; CTE of the ceramic is too high; Irreparable	Alloy compatibility; Minimum thickness of the frame-work (0.2 mm non precious 0.3 mm precious)
Shade too bright; Less translucency	Preheattemperature too high	Lower preheattemp.; for app. 50°C / 122°F
Surface too rough	Firing temperature too low	Increase firing temp.
Ceramic has porosities	Firing temp. too low; Vacuum too late; Vacuumlevel too low	Increase firing temp.; Lower vacuum start temp.; Check furnace and pump Comprobar horno y bomba
Ceramic has not enough glaze	Hold time without vacuum too short	Lengthen up hold time
Ceramic has too much glaze; Edges are getting round	Firing temp too high; Hold time too long	Lower firing temp.; Shorten hold time

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Wohlwend AG