IPS InLine®

[en] Instructions for Use

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1 Intended use

Intended purpose

Ceramic veneering materials for metal frameworks (CTE range of $13.8\text{--}15.0 \times 10^{-6}/\text{K}$ ($25\text{--}500^\circ\text{C}$); fabrication of layered veneers.

Patient target group

Patients with permanent teeth

Intended users / Special training

- Dentists (clinical procedure)
- Dental laboratory technicians (fabrication of restorations)

No special training required.

Use

For dental use only.

Description

IPS InLine® is a ceramic veneering material containing leucite. It is suitable for the fabrication of metal-ceramic restorations at firing temperatures higher than 900°C (1652°F). The product can be used for veneering alloys in the CTE range of $13.8\text{--}15.0 \times 10^{-6}/\text{K}$ ($25\text{--}500^\circ\text{C}$) and for fabricating veneers on refractory dies.

Product name	Product description
IPS InLine System Opaquer IPS InLine Gingiva Opaquer IPS InLine System Powder Opaquer IPS InLine System Intensive Powder Opaquer IPS InLine Gingiva Powder Opaquer	Masking of metal frameworks and adjustment of the base shade
IPS InLine Margin	Fabrication of ceramic margins on metal frameworks
IPS InLine Deep Dentin IPS InLine Dentin	Veneering of metal frameworks; fabrication of layered veneers in the anterior region
IPS InLine Incisal IPS InLine Transpa Incisal	Veneering of the occlusal and incisal areas of metal frameworks; fabrication of layered veneers in the anterior region
IPS InLine Transpa	Characterization of ceramic veneers on metal frameworks; characterization of layered veneers in the anterior region
IPS InLine Cervical Dentin IPS InLine Occlusal Dentin IPS InLine Mamelon IPS InLine Opal Effect IPS InLine Cervical Incisal	Characterization of ceramic veneers on metal frameworks; Characterization of layered veneers in anterior teeth
IPS InLine One Dentoisal	Veneering of metal frameworks
IPS InLine Add-On Margin IPS InLine Add-On IPS InLine System Add-On $690^\circ\text{C}/1274^\circ\text{F}$	Adjustment of ceramic margins and restorations on metal frameworks
IPS InLine Gingiva	Veneering of metal frameworks in the gingival area
IPS InLine Intensive Gingiva	Characterization of Gingiva materials
IPS InLine System Opaquer Liquid	For diluting IPS InLine System Opaquer
IPS InLine System Powder Opaquer Liquid	For mixing IPS InLine System Powder Opaquer
IPS InLine Margin Build-Up Liquid	For mixing IPS InLine Margin layering materials
IPS InLine System Build-Up Liquid (L, P)	For mixing IPS InLine layering materials

Technical data

Properties of IPS InLine according to ISO 6872:2015 + Amd 1:2018 and ISO 9693:2019. IPS InLine has been classified as a Type I, Class 1 dental ceramic.

Physico-chemical properties of the components of IPS InLine according to ISO 6872:2015 + Amd 1:2018 and ISO 9693:2019

Product component	Property		Radioactivity		Biaxial flexural strength σ		Chem. solubility		Bond strength τ_b	
			[Bq/g U ²³⁵]		[MPa]		[$\mu\text{g}/\text{cm}^2$]		[MPa]	
	A ¹⁾	B ²⁾	A ¹⁾	B ²⁾	A ¹⁾	B ²⁾	A ¹⁾	B ²⁾	C ³⁾	B ²⁾
IPS InLine System Opaquer IPS InLine Gingiva Opaquer	≤ 1.0	$<< 1.0$	> 50	170	< 100	28	> 25	40		
IPS InLine System Powder Opaquer IPS InLine System Intensive Powder Opaquer IPS InLine Gingiva Powder Opaquer	≤ 1.0	< 0.03	> 50	127	< 100	28	> 25	43		
IPS InLine Margin	≤ 1.0	$<< 1.0$	> 50	98	< 100	23	n. a.	n. a.		
IPS InLine Deep Dentin IPS InLine One Dentoisal IPS InLine Dentin IPS InLine Incisal IPS InLine Transpa Incisal	≤ 1.0	$<< 1.0$	> 50	89	< 100	12	> 25	40		
IPS InLine Cervical Incisal IPS InLine Cervical Dentin IPS InLine Occlusal Dentin IPS InLine Mamelon Masse IPS InLine Transpa IPS InLine Gingiva IPS InLine Intensiv Gingiva	≤ 1.0	$<< 1.0$	> 50	89	< 100	12	> 25	40		
IPS InLine Opal Effect	≤ 1.0	$<< 1.0$	> 50	90	< 100	10	n. a.	n. a.		
IPS InLine Add-On	≤ 1.0	$<< 1.0$	> 50	92	< 100	10	n. a.	n. a.		
IPS InLine Margin Add-On	≤ 1.0	$<< 1.0$	> 50	94	< 100	12	n. a.	n. a.		
IPS InLine System Add-On 690 °C ³⁾	≤ 1.0	$<< 1.0$	> 50	108	< 100	20	n. a.	n. a.		

¹⁾ A: Acceptance criterion according to ISO 6872:2015 + Amd 1:2018

²⁾ B: Typical measured values (mean values)

³⁾ C: Acceptance criterion according to ISO 9693:2019

Thermal properties of the components of IPS InLine according to ISO 6872:2015 + Amd 1:2018 and ISO 9693:2019

Property	Property	CTE 2x ¹⁾	CTE 4x ²⁾	CTE 0 ³⁾	Glass transition temperature T_g	Thermal cycling test
		[10 ⁻⁶ /K]	[10 ⁻⁶ /K]	[10 ⁻⁶ /K]	[°C]	
IPS InLine System Opaquer IPS InLine Gingiva Opaquer		13.5 ± 0.5	13.7 ± 0.5	13.6 ± 0.5	605 ± 20	compliant
IPS InLine System Powder Opaquer IPS InLine System Intensive Powder Opaquer IPS InLine Gingiva Powder Opaquer		13.2 ± 0.5	13.4 ± 0.5	13.3 ± 0.5	605 ± 20	compliant
IPS InLine Margin		13.2 ± 0.5	13.5 ± 0.5	13.4 ± 0.5	600 ± 20	compliant
IPS InLine Deep Dentin IPS InLine One Dentoisal		12.4 ± 0.5	13.2 ± 0.5	12.8 ± 0.5	580 ± 20	n. a.
IPS InLine Dentin IPS InLine Incisal IPS InLine Transpa Incisal		12.4 ± 0.5	13.2 ± 0.5	12.8 ± 0.5	580 ± 20	compliant
IPS InLine Cervical Incisal IPS InLine Cervical Dentin IPS InLine Occlusal Dentin IPS InLine Mamelon Masse IPS InLine Transpa IPS InLine Gingiva IPS InLine Intensiv Gingiva		12.4 ± 0.5	13.2 ± 0.5	12.8 ± 0.5	580 ± 20	n. a.
IPS InLine Opal Effect		12.6 ± 0.5	13.5 ± 0.5	13.1 ± 0.5	595 ± 20	n. a.
IPS InLine Add-On		12.3 ± 0.5	12.8 ± 0.5	12.6 ± 0.5	455 ± 20	n. a.
IPS InLine Margin Add-On		13.3 ± 0.5	13.7 ± 0.5	13.5 ± 0.5	585 ± 20	n. a.
IPS InLine System Add-On 690 °C ³⁾		12.6 ± 0.5	13.3 ± 0.5	13.0 ± 0.5	440 ± 20	n. a.

¹⁾ CTE after two firing cycles, measuring range 25 – 500 °C

²⁾ CTE after four firing cycles, measuring range 25 – 500 °C

³⁾ CTE mean value (after two/four firing cycles), measuring range 25 – 500 °C

⁴⁾ The CTE values apply to a temperature range of 25 to 400 °C, as the glass transition temperature is below 500 °C.

Firing table IPS InLine® – Firing parameters

IPS InLine	Firing temperature t [°C/°F]	Stand-by temperature B [°C/°F]	Closing time S [min]	Heating rate t [°C/°F/min]	Holding time H [min]	Vacuum on V1 [°C/°F]	Vacuum off V2 [°C/°F]
1 st /2 nd Powder Opaquer firing	960/1760	403/757	4:00	100/180	2:00	450/842	959/1758
1 st /2 nd Paste Opaquer firing	930/1706	403/757	6:00	100/180	2:00	450/842	929/1704
1 st /2 nd Margin firing	930/1706	403/757	4:00	60/108	1:00	450/842	929/1704
1 st Dentin / Incisal / Gingiva / Dentoisal firing	910/1670	403/757	4:00	60/108	1:00	450/842	909/1668
2 nd Dentin / Incisal / Gingiva / Dentoisal firing	900/1652	403/757	4:00	60/108	1:00	450/842	899/1650
Margin Add-On firing	900/1652	403/757	4:00	60/108	1:00	450/842	899/1650
Corrective firing after Dentin / Incisal firing, Add-On	860/1580	403/757	4:00	60/108	1:00	450/842	859/1578
Stain firing using IPS Ivocolor	830/1526	403/757	6:00	60/108	1:00	450/842	829/152
Glaze firing using IPS Ivocolor	830/1526	403/757	6:00	60/108	1:00	450/842	829/152
Add-On after glaze firing (690 °C/1274 °F)	690/1274	403/757	4:00	60/108	1:00	450/842	689/1274

IPS InLine® Veneer – Firing parameters

IPS InLine Veneer	Firing temperature t [°C/°F]	Stand-by temperature B [°C/°F]	Closing time S [min]	Heating rate t [°C/°F/min]	Holding time H [min]	Vacuum on V1 [°C/°F]	Vacuum off V2 [°C/°F]
Wash firing	830/1526	403/757	4:00	60/108	1:00	450/842	829/1524
Cervical firing	940/1724	403/757	8:00	60/108	1:00	450/842	939/1722
Dentin / Impulse firing	940/1724	403/757	8:00	60/108	1:00	450/842	939/1722
Incisal firing	930/1706	403/757	8:00	60/108	1:00	450/842	929/1704
Glaze firing using IPS Ivocolor	830/1526	403/757	8:00	60/108	1:00	450/842	829/1524

Note: All the firing programs described in these instructions for use do not involve active cooling. After the holding time, the furnace heater switches off and the furnace head cools down with a device-related cooling gradient due to the time-controlled opening of the furnace head. If long-term cooling is conducted, the closed furnace head cools down to 800 °C /1472 °F or 700 °C/1292 °F after the heater is switched off, at which point the time-controlled opening of the furnace head provides the device-related cooling gradient. These firing parameters are guidance values. They are valid for the Programat® furnaces from Ivoclar. Deviations (approx. ± 10 °C/18 °F) may occur:

- depending on the furnace generation
- with ceramic furnaces from other manufacturers
- due to regional differences in the power supply or if several electrical devices are operated on the same circuit.

Indications

Missing tooth structure in anterior and posterior teeth, partial edentulism in the anterior and posterior region

Types of restorations

- Crowns
- 3–14-unit bridges
- Veneers

Contraindications

The use of the product is contraindicated if the patient is known to be allergic to any of its ingredients.

Limitations of use

- Untreated bruxism (a splint is indicated after incorporation)
- Mixing with and processing in conjunction with other dental ceramics
- Failure to observe the necessary minimum connector and framework thicknesses
- Exceeding or falling short of the stipulated veneering layer thicknesses
- Failure to observe the layer thickness ratio between the framework and layering ceramic
- Veneering of dental alloys not within the stipulated CTE range
- Veneering of titanium and zirconium oxide frameworks
- The final restoration must not be reused.

Side effects

There are no known side effects to date.

Interactions

There are no known interactions to date.

Clinical benefit

- Reconstruction of chewing function
- Restoration of esthetics

Composition

Product components powders and pastes		Ingredient
Opaquer	IPS InLine System Opaquer IPS InLine Gingiva Opaquer	– Leucite glass-ceramic – Glass – Zirconium dioxide – Inorganic pigments – Glycols (butandiol, glycerine) – Polyvinylpyrrolidone – Highly dispersed silicic acid
	IPS InLine System Powder Opaquer IPS InLine System Intensive Powder Opaquer IPS InLine Gingiva Powder Opaquer	– Leucite glass-ceramic – Zirconium dioxide – Inorganic pigments
Margin	IPS InLine Margin	– Leucite glass-ceramic – Inorganic pigments
Deep Dentin Dentoisal Dentin Incisal Transpa Incisal	IPS InLine Deep Dentin IPS InLine One Dentoisal ^{*)} IPS InLine Dentin* IPS InLine Incisal* IPS InLine Transpa Incisal ^{*)}	– Leucite glass-ceramic – Opal glass – Inorganic pigments
Impulse	IPS InLine Cervical Incisal IPS InLine Cervical Dentin IPS InLine Occlusal Dentin IPS InLine Mamelon Masse IPS InLine Transpa IPS InLine Gingiva IPS InLine Intensiv Gingiva	
	IPS InLine Opal Effect	– Opalescent leucite glass-ceramic – Glass – Inorganic pigments
Add-On	IPS InLine Add-On	– Leucite glass-ceramic – Glass – Inorganic pigments
	IPS InLine Margin Add-On	– Leucite glass-ceramic – Inorganic pigments
	IPS InLine System Add-On 690 °C/1274 °F	– Leucite glass-ceramic – Glass – Inorganic pigments

^{*)} additionally contains small amounts of identification colour, which burns without residue when the ceramic is fired.

Product components ceramic mixing liquids	Ingredient
IPS InLine System Opaquer Liquid	– Glycols (butandiol, glycerine) – Polyvinylpyrrolidone
IPS InLine System Powder Opaquer Liquid	– Water – Glycol (propandiol) – Metallic salt – Carbonic acid
IPS InLine Margin Build-Up Liquid	– Water – Cellulose derivative
IPS InLine System Build-Up Liquid P IPS InLine System Build-Up Liquid L	– Water – Polyethyleneglycol – Metallic salt

2 Application

Preparation of the framework

Design the framework so that it reflects the shape of the tooth in a reduced form. This will ensure an even layer thickness and a uniform support of the veneering ceramic. Fabricate and oxidize the framework according to the directions of the alloy manufacturer.

Ceramic layering and processing

Application of the powder or paste opaquer – 1st + 2nd Opaquer firing

IPS InLine System paste and powder opaques evenly mask the framework in the respective shade. Two opaquer firings are recommended. The consistency of the paste opaquer can be individually adjusted using the IPS InLine System Opaquer Liquid. IPS InLine System Powder Opaquer Liquid is used to mix the powder opaquer to its desired consistency.

IPS InLine Paste Opaquer F

Paste Opaquer F can be used to reinforce the in-depth fluorescence.

- Either apply Opaquer F as a thin, third opaquer layer and fire (930 °C/1706 °F)
- Or mix up to 20% of Opaquer F with the conventional Paste Opaquer before the second layer is applied and fire at 930 °C/1706 °F.

Note: The firing tray with the opaqued metal framework should only be placed in the firing chamber and removed from it once the furnace head is completely open and the beeper has sounded.

IPS InLine® – conventional metal-ceramic

Ceramic margin – 1st + 2nd Margin firing

- When creating fired ceramic margins, make sure that the framework rather than the veneer is supported by the prepared tooth. The framework is thus reduced exactly to the inner edge of the chamfer or shoulder preparation.
- When creating ceramic margins, make sure to reduce the crown margin in such a way that it terminates approx. 0.5 to 0.8 mm above the lowest point of the chamfer or shoulder.
- Blast the inner and outer surfaces of the coping (especially the margins) and subsequently clean the coping with a steam cleaner.
- After that, apply a generous amount of IPS Margin Sealer and then, after drying, IPS Ceramic Separating Liquid.
- Supplement the reduced crown margin using Margin material.
- The gap created due to the shrinkage of the ceramic can be filled by conducting a second firing cycle.

Ceramic layering – 1st + 2nd Dentin and Incisal firing

Individual layering using the Deep Dentin and Opal Effect materials of the IPS InLine range enables you to achieve high-quality, true-to-nature restorations. Coat the framework with Deep Dentin materials to enhance the chroma and create a fluorescent base to characterize and individualize the tooth shade. Then build up the dentin portion as usual. Following the cut-back, supplement the incisal portion using Transpa Incisal materials and e.g. Opal Effect materials (OE1 – OE5 and violet) or more translucent Incisal materials in the incisal area and the Cervical Incisal materials orange and yellow in the cervical and body areas. A shade guide can be used in order to make selection of the shade of the different materials easier.

IPS InLine® One – One-Layer Metal-Ceramic

Ceramic layering – 1st + 2nd Denticisal firing

When the one-layer material (Denticisal) is used, the complete restoration is created with only one ceramic material. The restoration can be individualized using IPS InLine layering materials if required.

Adjustments using Add-On materials

The restoration can be adjusted and supplemented with the add-on materials IPS InLine Add-On Margin, IPS InLine Add-On and IPS InLine System Add-On 690 °C/1274 °F.

Individual finishing - Preparing for stain and glaze firing

Finish the restoration using diamonds and give it a true-to-nature shape and surface structure, such as growth lines and convex / concave areas.

IPS InLine® and IPS InLine® One – Completion – Stain and Glaze firing

For the individualized staining and glazing of restorations fabricated with IPS InLine and IPS InLine One, the IPS Ivocolor universal stains and glaze assortment can be used.

IPS InLine® – Veneers on refractory dies

Model fabrication

Duplicate the working model and subsequently cast with a customary refractory die material.

Wash firing

After degassing the refractory dies, apply IPS InLine Add-On mixed with the IPS Ivocolor Mixing Liquid allround in a thin layer and fire.

Cervical firing

Build up the marginal areas using a mixture of IPS InLine Dentin and, for example, Occlusal Dentin brown.

Dentin / Impulse firing

Internal layering is performed based on the natural model and consists of a dentin build-up and various effects. Individual layering with the Impulse materials enables mamelons as well as opalescent and translucent effects to be created.

Incisal firing

After applying the Incisal material, fire the veneer on the die.

Glaze firing

Apply IPS Ivocolor Glaze to the surface and fire.

Divesting of the veneer

Remove large amounts of die material using a grinding disc. Fine divestment is carried out by blasting with polishing beads at 1 bar (30 psi) pressure.

Conditioning the veneer for adhesive cementation

Etch the inner aspect of the veneer with IPS Ceramic Etching Gel for 120 seconds in preparation for adhesive cementation. Next, thoroughly rinse the restoration under running water and dry. Important: IPS InLine veneers must be seated using the adhesive technique.

Cementation

Restorations fabricated with IPS InLine can be cemented by the dentist either conventionally or adhesively.

3 Safety information

- In the case of serious incidents related to the product, please contact Ivoclar Vivadent AG, Bendererstrasse 2, 9494 Schaan/Liechtenstein, website: www.ivoclar.com, and your responsible competent authority.
- The current Instructions for Use are available in the download section of the Ivoclar Vivadent AG website (www.ivoclar.com).
- Explanation of symbols: www.ivoclar.com/eIFU.
- The Summary of Safety and Clinical Performance (SSCP) can be retrieved from the European Database on Medical Devices (EUDAMED) at <https://ec.europa.eu/tools/eudamed>. Basic UDI-DI: 76152082ACERA011ET

Warnings

- Observe the Safety Data Sheet (SDS) (available in the download section of the Ivoclar Vivadent AG website www.ivoclar.com).
- Do not inhale ceramic dust during finishing. Use suction equipment, a face mask and protective goggles.
- Hexane is highly flammable and detrimental to health. Avoid contact of the material with the skin and eyes. Do not inhale vapours. Keep material away from sources of ignition.

Disposal information

Remaining stocks or removed restorations must be disposed of according to the corresponding national legal requirements.

Residual risks

Users should be aware that any dental intervention in the oral cavity involves certain risks. Some of these risks are listed below:

- Chipping / fracture / decementation of the restorative material may lead to accidental swallowing or inhalation and dental re-treatment.
- Cement excess may lead to irritation of the soft tissue / gingiva. Progressive inflammation may lead to bone resorption or peri-implant disease.

4 Shelf life and storage

- Storage temperature 2–28 °C/36–82 °F: IPS InLine System Opaquer, IPS InLine System Opaquer Liquid, IPS InLine System Powder Opaquer Liquid, IPS InLine Margin Liquid, IPS InLine System Build-Up Liquids L/P
- Expiry date: See note on the packages.
- Ceramic powder: No special storage conditions are required.
- Do not use the product after the indicated date of expiration.
- Store in a vibration-free environment.
- Keep dry.
- Keep out of sunlight
- Before use, visually inspect the packaging and the product for damage. In case of any doubt, consult Ivoclar Vivadent AG or your local trade partner.

5 Additional information

Keep material out of the reach of children!

Not all products are available in all countries.

The product has been developed solely for use in dentistry. Processing should be carried out strictly according to the Instructions for Use. Liability cannot be accepted for damages resulting from failure to observe the instructions or the stipulated area of use. The user is responsible for testing the products for their suitability and use for any purpose not explicitly stated in the Instructions.

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Propriétés	Radioactivité		Résistance à la flexion biaxiale σ		Solubilité chimique		Valeurs d'adhésion τ_b	
	[Bq/g U ²³⁸]		[MPa]		[$\mu\text{g}/\text{cm}^2$]		[MPa]	
Composant du produit	A ¹⁾	B ²⁾	A ¹⁾	B ²⁾	A ¹⁾	B ²⁾	C ³⁾	B ²⁾
IPS InLine System Opaquer IPS InLine Gingiva Opaquer	≤ 1,0	<< 1,0	> 50	170	< 100	28	> 25	40
IPS InLine System Powder Opaquer IPS InLine System Intensive Powder Opaquer IPS InLine Gingiva Powder Opaquer	≤ 1,0	< 0,03	> 50	127	< 100	28	> 25	43
IPS InLine Margin	≤ 1,0	<< 1,0	> 50	98	< 100	23	n/d	n/d
IPS InLine Deep Dentin IPS InLine One Dentoisal IPS InLine Dentin IPS InLine Incisal IPS InLine Transpa Incisal	≤ 1,0	<< 1,0	> 50	89	< 100	12	> 25	40
IPS InLine Cervical Incisal IPS InLine Cervical Dentin IPS InLine Occlusal Dentin IPS InLine Masse Mamelon IPS InLine Transpa IPS InLine Gingiva IPS InLine Intensive Gingiva	≤ 1,0	<< 1,0	> 50	89	< 100	12	> 25	40
IPS InLine Opal Effect	≤ 1,0	<< 1,0	> 50	90	< 100	10	n/d	n/d
IPS InLine Add-On	≤ 1,0	<< 1,0	> 50	92	< 100	10	n/d	n/d
IPS InLine Margin Add-On	≤ 1,0	<< 1,0	> 50	94	< 100	12	n/d	n/d
IPS InLine System Add-On 690 °C ³⁾	≤ 1,0	<< 1,0	> 50	108	< 100	20	n/d	n/d

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Componente del producto	Propiedad		Radiactividad		Resistencia a la flexión biaxial σ		Solubilidad química		Resistencia de adhesión tb	
			[Bq/g U ²³⁵]		[MPa]		[μg/cm ²]		[MPa]	
	A ¹⁾	B ²⁾	A ¹⁾	B ²⁾	A ¹⁾	B ²⁾	A ¹⁾	B ²⁾	C ³⁾	B ²⁾
IPS InLine System Opaquer IPS InLine Gingiva Opaquer	≤1,0	<<1,0	>50	170	<100	28	>25	40		
IPS InLine System Powder Opaquer IPS InLine System Intensive Powder Opaquer IPS InLine Gingiva Powder Opaquer	≤1,0	<0,03	>50	127	<100	28	>25	43		
IPS InLine Margin	≤1,0	<<1,0	>50	98	<100	23	n.d.	n.d.		
IPS InLine Deep Dentin IPS InLine One Dentoisal IPS InLine Dentin IPS InLine Incisal IPS InLine Transpa Incisal	≤1,0	<<1,0	>50	89	<100	12	>25	40		
IPS InLine Cervical Incisal IPS InLine Cervical Dentin IPS InLine Occlusal Dentin IPS InLine Mamelon Masse IPS InLine Transpa IPS InLine Gingiva IPS InLine Intensive Gingiva	≤1,0	<<1,0	>50	89	<100	12	>25	40		
IPS InLine Opal Effect	≤1,0	<<1,0	>50	90	<100	10	n.d.	n.d.		
IPS InLine Add-On	≤1,0	<<1,0	>50	92	<100	10	n.d.	n.d.		
IPS InLine Margin Add-On	≤1,0	<<1,0	>50	94	<100	12	n.d.	n.d.		
IPS InLine System Add-On 690 °C*)	≤1,0	<<1,0	>50	108	<100	20	n.d.	n.d.		

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Komponenti i produktit	Radioaktiviteti		Rezistenca në përkuqe biaksiale σ		Tretshmëria kim.		Rezistenca e lidhjes tb	
	[Bq/g U ²³⁵]		[MPa]		[μg/cm ³]		[MPa]	
	A ¹⁾	B ²⁾	A ¹⁾	B ²⁾	A ¹⁾	B ²⁾	C ³⁾	B ²⁾
IPS InLine System Opaquer IPS InLine Gingiva Opaquer	≤ 1,0	<< 1,0	> 50	170	< 100	28	> 25	40
IPS InLine System Powder Opaquer IPS InLine System Intensive Powder Opaquer IPS InLine Gingiva Powder Opaquer	≤ 1,0	< 0,03	> 50	127	< 100	28	> 25	43
IPS InLine Margin	≤ 1,0	<< 1,0	> 50	98	< 100	23	n.d.	n.d.
IPS InLine Deep Dentin IPS InLine One Dentoisal IPS InLine Dentin IPS InLine Incisal IPS InLine Transpa Incisal	≤ 1,0	<< 1,0	> 50	89	< 100	12	> 25	40
IPS InLine Cervical Incisal IPS InLine Cervical Dentin IPS InLine Occlusal Dentin IPS InLine Mamelon Masse IPS InLine Transpa IPS InLine Gingiva IPS InLine Intensiv Gingiva	≤ 1,0	<< 1,0	> 50	89	< 100	12	> 25	40
IPS InLine Opal Effect	≤ 1,0	<< 1,0	> 50	90	< 100	10	n.d.	n.d.
IPS InLine Add-On	≤ 1,0	<< 1,0	> 50	92	< 100	10	n.d.	n.d.
IPS InLine Margin Add-On	≤ 1,0	<< 1,0	> 50	94	< 100	12	n.d.	n.d.
IPS InLine System Add-On 690°C ¹⁾	≤ 1,0	<< 1,0	> 50	108	< 100	20	n.d.	n.d.

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