Press to use, carefully read the instructions for use.

EN

GC Initial™ CAST NP Cobalt-Chrome (CoCr) based casting alloy for use in ceramic veneering (type 4)

For use only by a dental professional in the indications for use.

INDICATIONS FOR USE

GC Initial™ CAST NP is a cobalt chrome based casting alloy for metal veneering. Initial CAST NP is used for the fabrication of crowns & bridges with a suitable metal veneering ceramic or acrylic veneering material.

Thanks to its thermal expansion of 14,1 x 10⁻⁶ K⁻¹ it is suitable for the most common metal-bonded veneering ceramics on the market. It has outstanding properties as an optimal casting performance, an extra fine molecular structure resulting in smooth, compact surfaces offering an easy and straightforward processing. The alloy offers a reduced oxide formation even after repeated ceramic firings.

GC Initial™ CAST NP is produced conform the ISO 90693-1 and ISO 22674 standards. It has passed the cytotoxicity test according to the norm ISO 10993-5. The alloy is highly corrosion resistant and it is free of beryllium, indium and gallium.

COMPOSITION

Cobalt	Co	63%
Chromium	Cr	28%
Tungsten	W	3%
Niobium	Nb	1%
Silicon	Si	4%
Fe / Mo	Fe / Mo	1%

ALLOY CHARACTERISTICS

Yield strength: 475 MPa
Tensile strength: 1253 - 1304°C
Casting temperature: 1450 °C
Highest firing temperature: 950°C
Thermal expansion coefficient: 14,1 x 10⁻⁶ K⁻¹ (25 - 600°C)
Density: 8,4 g/cm³
Vickers Hardness: 324 HV10
Percentage elongation at fracture: 3,4%
Yield strength at 1000 °C: 475 MPa
Modulus of elasticity: 194 GPa
Colour: White
Recommended ceramics: GC Initial MC

INSTRUCTIONS FOR USE

1. WAX-UP

Wax should be 0.5 mm thick to ensure that after manufacturing and process the metal wall thickness is at least 0.3 mm.
Avoid sharp edges and undercuts. Frameworks need to be anatomically reduced using a uniform thickness of the ceramic layer. Contractors should be modelled as strong and as high as possible.

Single crowns: cast on directly spined. Suggested round casting sprues: Ø 3.0 mm to 4.0 mm.
Bridges: a crossbar spanning length is advised. Connections to crossbar Ø 3 to 2.5 mm.
Distance from modulation to crossbar 2.5mm. Thickness of the crossbars 0.5 mm.
Thickness of sprues to crossbar Ø 3 to 4 mm.

The higher the volume of the casting the thicker the sprue should be. The casting sprue depends on the volume dimensions.

II. INVESTING

Use only phosphate bonded investment material for casting. GC Initial™ CAST NP is optimized for producing high precision fit castings in combination with smooth cast surfaces. Guidelines for the use of GC's phosphate bonded investments for Crown & Bridge can be found in the "GC Focus Edition - Phosphate Bonded Investments for C&B techniques". A pre-heating temperature of 820±5°C is advised. Follow the manufacturer's Instruction for Use related to the pre-heating times before casting.

III. CASTING

We recommend to only use fresh alloys in own crucibles. Do not re-use GC Initial CAST NP. Optimal results of the alloy can be guaranteed only if new material is used.
Melt GC Initial CAST NP in pre-heated crucibles used for this alloy only.

Amount of casting alloy to use = wax framework weight with sprues / wax density x metal density.

High-frequency melting equipment

When last ingots have melted in the pre-heated ceramic crucible and the incandescence shadow has disappeared immediately start the casting process. Do not use any flux.

Open flame melting

Use the pre-heated ceramic crucible and heat them evenly with circular movements. When ingots have melted, start centrifugal unit. Use multi-flame welding torches only. Do not use any flux.

Overheating the melted material could cause shrinkling hollows, micro-porosity and building of a coarse-grained structure and cause bridges to break or ceramic veneering to crack.

Follow the Instructions of the respective manufacturer of the casting devices for parameters and casting procedures.

After casting bench-cool the casting ring. Remove the investment material with care and sandblast using a sandblasting sandblaster with 100 to 125 µm aluminium oxide at a pressure of 3-4 bar.

IV. PREPARATION FOR VENEERING

Use tungsten to grind and finish the prepared framework. Use always the same grinding instruments to avoid contamination.
Do not use diamond burs. Grind the surface with slow speed and low pressure.
Grind the surface in one direction only to avoid overlying layers to form air bubbles after the first investment.
Sandblast the prepared framework using a non-recycling sandblaster with 110 to 150 µm aluminium oxide at a maximum pressure of 3-4 bar and then steam clean.

V. OXIDE FIRING

Fire for 5 min. under vacuum at 950-980°C (10°C more than oxide melting firing temperature). After firing, carefully sandblast the oxide layer using a non-recycling sandblaster with 110 to 150 µm aluminium oxide at a pressure of 2-5.3 bar, then steam clean. The framework should have a homogeneous grey surface. Proceed with firing process using a compatible veneering ceramic following the manufacturer's instruction.

VI. VENEERING

- Veneering with GC Initial MC: please refer to GC Initial MC technical Manual
- Veneering with other ceramics: please refer to the manufacturer's instruction for use
- Veneering with acrylics: please refer to the manufacturer's instruction for use

PACKAGING

- GC Initial CAST NP 250g
- GC Initial CAST NP 500g
- GC Initial CAST NP 1000g

CAUTION

- Cobalt-Chrome based alloys can seldom cause dermatitis on sensitive subjects.
A Patch-Test is thus advisable.
- Metal dusts and smoke are dangerous for health. Please wear safety glasses and a mask when using the alloy.
We recommend to keep batch records and link it to patient's file to allow it's complete traceability.
- This dust does not need any special preservation precautions. Nonetheless it has to be employed within the expiry date on the packaging.
- Disposal of all waste according to local regulations.
- Personal protective equipment (PPE) such as gloves, face masks and safety eyewear should always be worn.
- GC Initial CAST NP contains Cobalt, substance classified as CMR.
- Do not re-use GC Initial CAST NP. Optimal results of the alloy can be guaranteed only if new material is used.

GUARANTEE

All recommendations are based on GC findings and GC Intents and are therefore only to be used as guidelines. Any claim for compensation GC receives in any event, shall not exceed the value of the individual order of GC Initial CAST NP.

Some products referenced in the present IPU may be classified as hazardous according to GHS. Always familiarize yourself with the safety data sheets available at https://www.gc.europa.eu/sectors/medical-devices/contacts_en.

For the Summary of Safety and Clinical Performance (SSCP) and EUDAMED database (https://ec.europa.eu/eudamed/eudamed) or contact us at Regulatory@gc.europa.eu.

Undesired effects - Reporting
If you become aware of any kind of undesired effect, reaction or similar events experienced by use of this product, including those not listed in this instruction for use, please report them directly through the relevant authority system, by selecting the proper authority of your country accessible through the following link: https://ec.europa.eu/growth/sectors/medical-devices/contacts_en

In this way you will contribute to improve the safety of this product.

Last revised: 09/2024

For the Summary of Safety and Clinical Performance (SSCP) and EUDAMED database (https://ec.europa.eu/eudamed/eudamed) or contact us at Regulatory@gc.europa.eu.

Undesired effects - Reporting
If you become aware of any kind of undesired effect, reaction or similar events experienced by use of this product, including those not listed in this instruction for use, please report them directly through the relevant authority system, by selecting the proper authority of your country accessible through the following link: https://ec.europa.eu/growth/sectors/medical-devices/contacts_en

In this way you will contribute to improve the safety of this product.

Last revised: 09/2024

For the Summary of Safety and Clinical Performance (SSCP) and EUDAMED database (https://ec.europa.eu/eudamed/eudamed) or contact us at Regulatory@gc.europa.eu.

Undesired effects - Reporting
If you become aware of any kind of undesired effect, reaction or similar events experienced by use of this product, including those not listed in this instruction for use, please report them directly through the relevant authority system, by selecting the proper authority of your country accessible through the following link: https://ec.europa.eu/growth/sectors/medical-devices/contacts_en

In this way you will contribute to improve the safety of this product.

Last revised: 09/2024

For the Summary of Safety and Clinical Performance (SSCP) and EUDAMED database (https://ec.europa.eu/eudamed/eudamed) or contact us at Regulatory@gc.europa.eu.

Press to use, carefully read the instructions for use.

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