



Guided &
conventional
SURGERY



Developed By P-I Brånemark




DIGITAL



Fully Guided Surgery

novel stepwise approach



Precision • Sensitivity • Versatility



Guided &
conventional



DIGITAL
prosthetics

The P-I Brånemark • Guided & Conventional novel stepwise approach was developed to simplify clinical procedures and reduce the potential risks associated with guided surgery providing flexibility for conventional procedures by using the same Instruments and Kit.

Designed to increase tridimensional accuracy of Implant positioning, the Instruments overcome the inherent tactile feedback reduction, visualization and access limitations of surgical site preparation and Implant installation when using surgical guides.

Combined with the MT-F Implant System, the novel Guided & Conventional concept is a stepwise approach to the modern, immediate and digitalized implantology which keeps clinicians focused on patients and in control of procedures, even when plans change during surgery.

The P-I Guided & Conventional provides a predictable, minimally invasive and versatile solution for precise tridimensional positioning of Implants with an optimal esthetic outcome when used with integrated digital diagnostic, planning and prosthetic software.



 MT-F

 **D** DIGITAL
prosthetics





Contents

Guided & Conventional

Novel Stepwise Approach Technology 5

Combinations 10

Dimensions 11

Alignments 13

Kit • Guided & Conventional Surgery 14

Instruments 15

MT-F Guided Surgical Sequence 19



 Guided & conventional

Wash Tray
Surgical + Prosthetic





State-of-the-art

Less friction

Only titanium Stop is in contact with Sleeve

Less trauma

Constant apical angle • Diamond Like Carbon

High performance

P-I custom technology




Exceptional cutting performance

P-I Conical Drills' performance in dense bone, at the highest recommended rotation, without gradual diameter increments and applying constant feeding, present a very low friction coefficient therefore lower temperature transmission to bone tissue. Data on file.




Drill frictionless through Sleeve

 Screw-retained Stops are the only contact with Sleeves • Titanium to titanium
 Initial and Conical Drills do not rotate against surgical guide or Sleeves
 Potential reduction of both temperature and debris release on surgical site



Fully Guided • Precision in all steps

Crestal Drills are continuously guided through Sleeves in Offset [9] and [10.5]
 Initial and Conical Drills are oriented by previous drilling steps and Stops
 Screw-retained Stops provide accuracy and stability



Same Instruments for Guided & Conventional

Guided Surgery - GS identification

For conventional procedure

-  Guided • All
-  Narrow - N
-  Regular - R • Long - L
-  Regular - R • Short - S



Short and Long
Drills GS



Insertion
Drivers GS

Unmatched versatility


Less Instruments • Keyless • Mountless

Shorter Drills and Stops can be used for incremental preparation depth of osteotomy

Stops can be pre-assembled according to planning


▼ Incremental depth option

MT-F 3.75 x 10 • Offset [9]




Full depth option

MT-F 3.75 x 10 • Offset [9]



Surgical Sequence according to bone density

Implant installation



Effective Insertion Torque Value ITV ↗

Reduced contact area • Diamond Like Carbon • Less friction on Sleeve
Safe pick-up, installation and removal • Avoids undesired movement
Possible further Implant submersion through Sleeve

Handpiece • Manual • Torque Wrench

Combinations


Narrow


Regular
Short


Regular
Long

Drill • GS



Stop • GS 



Sleeve • GS





[10.5] Offset [10.5] and Biological Width oriented positioning should consider drilling depth compensation. Offset [10.5] cannot be used with MT-F h = 15. 

h=18 For MT-F h=18 use Drill Stop - GS - 15, and the final drilling depth is determined conventionally without surgical guide.

4.8 For MT-F Ø4.8 use Conical Drill - GS - 3.8 for soft bone. The following drilling steps are performed conventionally without surgical guide. Implant insertion can be performed through Sleeve - GS - R.

 Consult Instructions for Use and Guided Surgical Sequence. Measurements in millimeters. Images for illustrative purposes only. Verify available libraries and updates before planning.

Dimensions



Pin

h 35



Narrow

40



Regular Short

35



Regular Long

40


Drill • GS



Stop • GS



Sleeve • GS



d	1.8	3.5	5
D	4.2	4	6
h	0.7	1	1
H	7	4	4
f	5	4.5	6
F	5	4.5	7



Sleeve - GS - N can also be used in limited spaces and for marking the initial drilling steps preceding conventional (free hand) surgical site finalization for the sleeveless installation of MT-F Ø3.75, 4.1 and 4.8.



Consult Instructions for Use, Implant Dimensions and Guided Surgical Sequence. Measurements in millimeters. Images for illustrative purposes only. Verify available libraries and updates before planning.




MT-F Dimensions



One Interface for all Implants



Alignments





- The flat part of the Sleeve - R and the Sleeves' slots can be positioned buccally to align with the prosthetic orientation of the Implant Insertion Driver (flat areas, upper hexagon and precision dot) which are aligned to Implant and Prosthetic Components' indexation.

- For angled Conical Abutments, and in case of mesiodistal space limitation, the flat part of the Sleeve - R can be mesially oriented leading to a Sleeve corner to buccal positioning.
- When the Sleeve - R is positioned corner to buccal, the prosthetic orientation of the Implant Insertion Driver can be directed to the corner of the Sleeve - R, and between the Sleeves' slots, for a buccal indexation of Implant and Prosthetic Components (except with angled Conical Abutments).

⚠ A minimum distance of 3mm should be ensured between 2 Implants

Depth verification

Offset [9] example



Wash Tray • GS

Surgical + Prosthetic



Guided & Conventional

181040

All Instruments



Guided +

181041

Guided Surgery Instruments only
includes Prosthetic Instruments



Tray • GS

131162

Tray only

width	190
height	61
depth	138



Stop Organizer • GS

141399

Organizer [9] only

Instruments • GS

Drill



Crestal Drill

Narrow

Regular

141330

141335



Ø N • Long

R • Short

R • Long



Initial

2.2 141378

141336

141340



Conical

2.8 141381

141337

141341

3.4 141381

141338

141342

3.8

141339

141343



Stop • GS

Offset [9]

Set • 11 Stops - GS

181043

without organizer

h



6

141365

7

141366

8.5

141367

10

141360

141368

11.5 141361

141369

13 141362

141370

15 141363

141371


Sleeve



141375

141376

Implant Insertion Driver



131153

131152



Insertion Driver must be completely connected into the Implant internal hexagonal indexation. When used manually or with Torque Wrench, the upper hexagonal portion of the Insertion Driver must be entirely connected to the hexagon of the Driver Adapter.

Guide Fixation



Sleeve - Pin

141382



Pin Drill Ø1.7

141377



Fixation Pin

131167



Surgical Accessories

Narrow

Regular

Punch



Soft Tissue Punch

141331

141332

Dense Drill



Dense Drill

141344

141345

Torque Wrench Kit



C14569

HIGH TORQUE VERSION
originally developed for OSPOL • P-I

Prosthetic Instruments



Hexagonal Ø 1.2

Short	131010
Medium	131011
Long	131012



Conical Abutment Ø 2.0

Short	131016
Medium	131017

Handling Tool, Angled

102964



Retriever MT

Short	141564
Medium	131131

Surgical Accessories • GS



In the event Stop cannot be easily removed by hand, use handpiece in reverse or the Surgical Adapter, with or without Torque Wrench, to stabilize the Drill and apply manual counter torque on the Stop using the Stop Removal Tool.



Depth Probe • GS
[9] Offset tip + Round tip

141439



Stop Removal Tool

131172

■ Manufactured by Elos MedTech Pinol A/S. Torque Wrench Kit includes Surgical & Prosthetic Adapters.
● Stop Removal Tool should be used when Surgical Adapter is not sufficient for the removal of Stop from Drill.
Depth Probe is not part of Kit • GS and should be ordered separately.


Stop • Conventional



STOP
conventional

FREE
HAND

Short • Long
Drills GS



R • Short

R • Long

h	R • Short	R • Long
6	141592	
7	141593	
8.5	141594	
10	141590	
11.5	141591	141584
13	141600	141586
15	141601	141586
18		141602

MT-F Guided Surgical Sequence



rpm 600 - 1,200
Lowest possible rpm



In-Out Δ Coordinated in-and-out movement of Drills for better cooling



ITV ≤ 70 Ncm
Insertion Torque Value



Full Length Prepare at planned full length of Implant position combining the Drills GS and Stops GS



Irrigation Constant irrigation to the insertion margin of Conical Drills

Δ Drills are less than 1 mm longer than Drill marks



3.3

Crestal	2.2	2.8	3.4
✓	✓	S ↓	D ↓
		N	



3.75

Crestal	2.2	2.8	3.4
✓	✓	S	N ↓
			D



4.1

Crestal	2.2	2.8	3.4	3.8
✓	✓	S	N	D
	 Optional for Dense>		



4.8

Crestal	2.2	2.8	3.8	4.6
✓	✓	✓	S	N ↓
				D

▼ Incremental depth option

Bone Density
S Soft N Normal D Dense

↓ 6 mm Drill tip



The subsequent Conical Drill, in terms of diameter, should be considered with a drilling depth of 6 mm, in order to not exceed 70 Ncm of Insertion Torque Value. The use of Dense Drills (15 - 50 rpm) can also be considered to lower the Insertion Torque Value. The use of Soft Tissue Punch (15 - 50 rpm) precedes Crestal Drill - GS use and can also be used in conventional surgery procedures.

Consult Instructions for Use and Guided Surgical Sequence in detail.



Guided Surgical Sequence

in detail

MT-F		Guide		Guided Surgical Sequence						
Ø	h	Offset	Sleeve 	Punch optional 15-50 rpm	1	2	3 Bone Density			Implant Driver ≤ 70 Ncm
					Crestal 600 - 1,200 rpm	Initial 600 - 1,200 rpm	SOFT	NORMAL	DENSE	
3.3	8.5	[9]	R	R	R	Ø 2.2 + 8.5	Ø 2.8 + 6 ▲	Ø 2.8 + 8.5	Ø 3.4 + 6	R
3.3	10.0	[9]	R	R	R	Ø 2.2 + 10	Ø 2.8 + 6 ▲	Ø 2.8 + 10	Ø 3.4 + 6	R
3.3	11.5	[9]	R	R	R	Ø 2.2 + 11.5	Ø 2.8 + 6 ▲	Ø 2.8 + 11.5	Ø 3.4 + 6	R
3.3	13.0	[9]	R	R	R	Ø 2.2 + 13	Ø 2.8 + 6 ▲	Ø 2.8 + 13	Ø 3.4 + 6	R
3.3	15.0	[9]	R	R	R	Ø 2.2 + 15	Ø 2.8 + 6 ▲	Ø 2.8 + 15	Ø 3.4 + 6	R
3.75	6.0	[9]	R	R	R	Ø 2.2 + 6	Ø 2.8 + 6	Ø 3.4 + 6 ▲	Ø 3.4 + 6	R
3.75	7.0	[9]	R	R	R	Ø 2.2 + 7	Ø 2.8 + 7	Ø 3.4 + 6 ▲	Ø 3.4 + 7	R
3.75	8.5	[9]	R	R	R	Ø 2.2 + 8.5	Ø 2.8 + 8.5	Ø 3.4 + 6 ▲	Ø 3.4 + 8.5	R
3.75	10.0	[9]	R	R	R	Ø 2.2 + 10	Ø 2.8 + 10	Ø 3.4 + 6 ▲	Ø 3.4 + 10	R
3.75	11.5	[9]	R	R	R	Ø 2.2 + 11.5	Ø 2.8 + 11.5	Ø 3.4 + 6 ▲	Ø 3.4 + 11.5	R
3.75	13.0	[9]	R	R	R	Ø 2.2 + 13	Ø 2.8 + 13	Ø 3.4 + 6 ▲	Ø 3.4 + 13	R
3.75	15.0	[9]	R	R	R	Ø 2.2 + 15	Ø 2.8 + 15	Ø 3.4 + 6 ▲	Ø 3.4 + 15	R
3.75	18.0	[9]	R	R	R	Ø 2.2 + 15	Ø 2.8 18 mark	Ø 3.4 6 mark ▲	Ø 3.4 18 mark	R
4.1	6.0	[9]	R	R	R	Ø 2.2 + 6	Ø 2.8 + 6	Ø 3.4 + 6 ▲	Ø 3.8 + 6	R
4.1	7.0	[9]	R	R	R	Ø 2.2 + 7	Ø 2.8 + 7	Ø 3.4 + 7 ▲	Ø 3.8 + 7	R
4.1	8.5	[9]	R	R	R	Ø 2.2 + 8.5	Ø 2.8 + 8.5	Ø 3.4 + 8.5 ▲	Ø 3.8 + 8.5	R
4.1	10.0	[9]	R	R	R	Ø 2.2 + 10	Ø 2.8 + 10	Ø 3.4 + 10 ▲	Ø 3.8 + 10	R
4.1	11.5	[9]	R	R	R	Ø 2.2 + 11.5	Ø 2.8 + 11.5	Ø 3.4 + 11.5 ▲	Ø 3.8 + 11.5	R
4.1	13.0	[9]	R	R	R	Ø 2.2 + 13	Ø 2.8 + 13	Ø 3.4 + 13 ▲	Ø 3.8 + 13	R
4.1	15.0	[9]	R	R	R	Ø 2.2 + 15	Ø 2.8 + 15	Ø 3.4 + 15 ▲	Ø 3.8 + 15	R

▼ Incremental depth option

Possible use of shorter Stops - GS preceding final depth for increased angular accuracy through Sleeve

R = Regular

● Short ● Long

▲ Optional step for dense bone



Guided Surgical Sequence

in detail

MT-F		Guide		Guided Surgical Sequence						
Ø	h	Offset	Sleeve 	Punch optional 15-50 rpm	1	2	3 Bone Density			Implant Driver ≤ 70 Ncm
					Crestal 600 - 1,200 rpm	Initial 600 - 1,200 rpm	SOFT	NORMAL	DENSE	
4.8	6.0	[9]	R	R	R	Ø 2.2 + 6	Ø 2.8 + 3.8 + 6	Ø 4.6 6 mark ▲	Ø 4.6 6 mark	R
4.8	7.0	[9]	R	R	R	Ø 2.2 + 7	Ø 2.8 + 3.8 + 7	Ø 4.6 6 mark ▲	Ø 4.6 7 mark	R
4.8	8.5	[9]	R	R	R	Ø 2.2 + 8.5	Ø 2.8 + 3.8 + 8.5	Ø 4.6 6 mark ▲	Ø 4.6 8.5 mark	R
4.8	10.0	[9]	R	R	R	Ø 2.2 + 10	Ø 2.8 + 3.8 + 10	Ø 4.6 6 mark ▲	Ø 4.6 10 mark	R
4.8	11.5	[9]	R	R	R	Ø 2.2 + 10	Ø 2.8 + 3.8 11.5 mark	Ø 4.6 6 mark ▲	Ø 4.6 11.5 mark	R
4.8	13.0	[9]	R	R	R	Ø 2.2 + 10	Ø 2.8 + 3.8 13 mark	Ø 4.6 6 mark ▲	Ø 4.6 13 mark	R
3.3	10.0	[9]		N	N	Ø 2.2 + 10	Dense N ▲	Ø 2.8 + 10	Dense N	N
3.3	11.5	[9]	N	N	N	Ø 2.2 + 11.5	Dense N ▲	Ø 2.8 + 11.5	Dense N	N
3.3	13.0	[9]	N	N	N	Ø 2.2 + 13	Dense N ▲	Ø 2.8 + 13	Dense N	N
3.3	15.0	[9]	N	N	N	Ø 2.2 + 15	Dense N ▲	Ø 2.9 + 15	Dense N	N

▼ Incremental depth option

Possible use of shorter Stops - GS preceding final depth for increased angular accuracy through Sleeve

R = Regular N = Narrow

● Short ● Long ● Narrow

▲ Optional step for dense bone



Full Prepare at planned full length of Implant position combining the Drills - GS and Stops - GS



In-Out ▲ Coordinated in-and-out movement of Drill for better cooling



Irrigation Constant cooled irrigation to the insertion margin of Drills



The subsequent Conical Drill, in terms of diameter, should be considered with a drilling depth of 6 mm, in order to not exceed 70 Ncm of Insertion Torque Value. The use of Dense Drills (15 - 50 rpm) can also be considered to lower the Insertion Torque Value.



[10.5] Offset [10.5] and Biological Width oriented positioning should consider drilling depth compensation

 Consult Instructions for Use

MT-F Guided Surgical Sequence

in detail



⚠ When using Stops, only use them with designated Drills.
pibranemark.com



Guided &
conventional
SURGERY




Developed By P-I Brånemark

P-I GUIDED-&-CONVENTIONAL-00HENG-1710-2024
P-I Brånemark. All rights reserved.

pibranemark.com

 SIC invent Sweden AB
Stora Åvägen 21
43634 Askim, Sweden



files for open
source software
pibranemark.com

exocad DentalCAD® exoplan

