



Developed By P-I Brånemark



○ MT-F
THE NEXT GENERATION



SMART GUIDE

MT-F

THE NEXT GENERATION



SIMPLE EXPERIENCE FOR EXCEPTIONAL OUTCOMES

Enhanced Biological Metrics to unlock immediate replacement potential

The P-I Implant Systems were developed by Professor Per-Ingvar Brånemark, the Osseointegration pioneer, jointly with scientists from renowned universities and the P-I Research & Development team to meet the modern implant dentistry demands.

In 2012, Ospol AB Sweden was acquired, and key technologies were integrated in the P-I solutions.

With the human biology, long-term expertise, clinical and scientific evidences as a foundation, our main objective is to support you in patient-focused treatments by providing Implant Systems that represent: Simplification • High Performance • Safety and Longevity

MT-F is the Next Generation System, a result of the P-I Brånemark fundamentals evolutionized by outstanding Biological Metrics and Simplicity.



Developed By P-I Brånemark

Contents

MT-F UNIQUENESS	06
------------------------	----

P-I Technologies

 MT Interface	10
 OSPOL Surface	14
 ENHANCED Ceramic Surface	18
Custom Instruments	19

MT-F Implants	26
----------------------	----

Soft Tissue Healing	29
---------------------	----

Conical Abutment	31
------------------	----

Abutment Cemented Cylinder	33
----------------------------	----

Contour & Esthetic Abutments	35
------------------------------	----

Cylinder over Implant	37
-----------------------	----

 DIGITAL	39
--	----

Kit & Instruments	44
-------------------	----

 Guided & Conventional SURGERY	47
---	----

Resonance Frequency	53
---------------------	----

Locator, Overdenture	53
----------------------	----

Surgical Sequence	54
-------------------	----


Torques, Implant Packaging	57
----------------------------	----

LifeTime Guarantee	59
--------------------	----



Adaptive bone contact

Multiple transitions • Interpolated core



Less bone displacement

In all bone densities



● MT-F

Less Trauma • Site Engagement

Cutting threads • Pronounced depth in all sections

● MT-F UNIQUENESS

The multiplicity of interpolated core transitions associated with the P-I Conical Drills site preparation, and the gradual evolution of the pronounced depth cutting threads, result in a gentle implant-to-osteotomy engagement in all sections independently.

These unique geometrical combinations provide greater initial contact area with significantly less bone displacement and compression enhancing the Biological Metrics.


↑ ISQ

ITV ↗

RTQ %

Less compression • Progressive torque • Greater area

MT-F displaces significantly less bone volume and achieves similar or higher Insertion Torque Value in all bone densities, exhibiting greater area in comparison to the leading competitive tapered-active implants of similar dimensions. Data on file.





Enhanced Biological Metrics

↑ ISQ

ITV ↗

RTQ %

Biological Metrics

High initial and secondary Implant Stability Quotient • ISQ measurements by Resonance Frequency Analysis • RFA in association with sufficient Insertion Torque Value • ITV and low rotational micro-mobility, indicated by the proportional Removal Torque • RTQ% to the obtained ITV, are relevant Biological Metrics and critical success factors for the prosthetic rehabilitation of patients with implants in post extraction, healed and compromised sites, low density bone and in combination with tissue regeneration techniques.

The P-I expertise

Our expertise related to ISQ using RFA micromovement measurements to clinically monitor Osseointegration and to determine when to load implants, originates from the acquisition of Ospol AB in 2012.

Ospol AB and Osstell AB were sister companies established in Sweden and developers of an Implant System and RFA measurement technologies, respectively. The Ospol AB developments of the last 20 years are comprised within the P-I Implant Systems and the newest technologies are present in The Next Generation • MT-F System.



sister companies



Increased coronal space

Slightly inward flange

Cortical stability

Micro Patterns

Adaptive bone contact

Interpolated core transitions



Gradual thread evolution


Pronounced depth in all sections

Early engagement


Gentle cutting • Dual thread

Axial insertion control

Biological Width positioning



○ MT-F



Platform switching

One interface
for all Implant dimensions

Concave emergence
Soft tissue volume

Peri-implant tissue preservation

● **MT Interface**

Biological
Width



Double
Sealing

System stability

Interface in clinical use for 15+ years • Superior biomechanics • Double Sealing

The P-I Morse Taper is an original technology. Highlighting 8.5° x 2 conical indexed, 3mm long, the P-I MT Interface offers a high torsional yield and fatigue strength when compared to other leading systems and was adopted by a global leader in 2015. The MT-F Ø 3.3 Implant can withstand static load of approximately 600N. Data on file.

The high-preload Double Sealing mechanism has easy prosthetic reversibility, seals the Abutment on the MT Interface and the MT Screw on the Abutment, stabilizing the system, minimizing micromovement and microleakage in comparison to certain leading systems under simulated occlusal stress. The Double Sealing is an important hypothesis for the clinical consideration of MT-F Implant placement observing Biological Width principles.



One Interface


Various prosthetic platforms



High Pre-Load

Effective sealing

P-I Coating +



Easy reversibility

Low stress to peri-implant tissue



Sealing starts at provisionalization

P-I Coating+ is a biocompatible layer that reduces friction and, combined with the MT Screw optimized geometry, provides a substantially higher and homogeneous pre-load, clamping, in comparison to titanium screws at the same tightening torque of 25 Ncm. Data on file. MT Retriever is used to cancel the morse sealing and safely remove Abutments.

Prosthetic Overview

Screw-Retained

Cement-Retained

Screw • Cemented

D DIGITAL
prosthetics



D



N **R** **W**



D One MT Interface for all MT-F Implant dimensions

 Locator® Abutments available

 Same MT Screw and Prosthetic Driver Ø 1.2 for all Abutments, except straight Conical Abutment and Locator®. All P-I Components are supplied with the respective screw, except Milling Blank.

D See DIGITAL Prosthetics. Libraries are available at piبرانemark.com and/or from CAD/CAM software system. Direct Geometries, 3D model Analogs and Geometry Contours available. Link C post has Dentsply Sirona, CEREC® dimensions. Always check latest library version and availability.


D files for open source software piبرانemark.com

3shape  exocad

Strong Osseointegration

REDUCTION OF BIOFILM INFECTIONS

BIOACTIVE



 **OSPOL**
SURFACE

Improved bone response

In comparison to rougher oxidized and blasted surfaces

Less bacterial adhesion

Equivalent to turned surfaces • Minimally rough

Chemically enhanced

Anodized • Bioactive ions

CP Ti
grade 4

Anodized
+
Chemistry

ZERO %
Aluminum
or Acid




Widely documented
Evolution of moderately rough surfaces




Direct and early
response

Oxide, micropores and crystal structures
greatly influence bone response



Biochemical bond,
bone in-growth
and mechanical
interlocking



Courtesy of : YT Sul, A. Wennerberg,
T. Albrektsson

Surface chemistry, anodic oxidation and ion
incorporation, enhance Osseointegration and
compensate for minimal roughness

OSPOL Surface was developed in the Gothenburg University, Sweden, and is documented in several publications. In continual evolution since 2000 and in clinical use for over 15 years, the OSPOL Surface is a modern technology for a rapid and strong bone response and its modification method was adopted by a global leader in 2019. Less prone to bacterial adhesion, it is a pioneer technology for chemical modification of thin anodized, oxidized, ion incorporation of smoother implants surfaces.



Higher ISQ for chemically- modified Surface

OSPOL Surface modification method achieves faster secondary Implant Stability Quotient • ISQ measured by Resonance Frequency Analysis RFA indicating potential for shorter healing periods.

Less bacterial
adhesion and
biofilm formation

BIOACTIVE

⚠ OSPOL Surface is easier to clean than rougher surfaces and its bioactivity reduces biofilm formation. The bacterial adhesion is similar to turned, machined surfaces.



Blasted, SLA



Turned



OSPOL

(!) Some conditions, whether combined or not, represent contraindications, limitations and risks, relative and absolute, for the treatment of patients with implants. There are several risk factors in Osseointegration widely described in literature. ISQ is a critical factor to clinically monitor Osseointegration. Data from pre-clinical studies.



ENHANCED Ceramic Surface



D DIGITAL
prosthetics

The P-I ENHANCED Ceramic Surface for Prosthetic Components is a smooth, tailored anodized surface technology that can potentially offer additional benefits for soft tissue in comparison to untreated titanium.

The acid-free P-I ENHANCED Ceramic Surface also exhibits a golden color for improved esthetics.



Kit

Surgical & Prosthetic



Stainless Steel

BIOSAFETY & DURABILITY

Conventional

Conical Drills



Less friction. Less trauma

Constant apical conical angle • 3 cutting areas

Corrosion protection

Wear resistance • Diamond Like Carbon

State-of-the-art

Special P-I design

Exceptional cutting performance

P-I Conical Drills' performance in dense bone, at the highest recommended rotation, without gradual diameter increments and applying constant feeding, present a very low friction therefore lower temperature transmission to bone tissue. Data on file.



Easy, simplified installation

Maximum of 3
low speed steps




(1) Except for MT-F Ø 4.8 Implant, 3 or 4 low speed steps are used. See Surgical Sequence. Spade and Round burr are optional. No pilot drill, counter sink or screw-tap required.



Insertion Driver


Handpiece • Manual • Torque Wrench





Biological Width

Driver reference



(!) The horizontal Implant Insertion Driver's mark is at approximately 3 mm and serves as a minimum Biological Width vertical reference for Implant platform positioning when completely covered by the lowest point of the soft tissue, the gingival margin. For further submersion, verify available prosthetic Component dimensions. Implant Insertion Driver dots, upper hexagon and flat areas are indexed to the Implant's hexagonal index and Sleeve • GS flat and slots. When used with Torque Wrench, the upper hexagonal portion of the Implant Insertion Driver should be entirely connected to the hexagon of the Driver Adapter.


 Guided &
conventional
SURGERY



Precision • Sensitivity • Versatility



Fully Guided Surgery novel stepwise approach



Drill frictionless
through sleeve

Potential reduction of both
temperature and debris
release on surgical site

Initial and Conical Drills do not rotate
against surgical guide or Sleeves



files for open
source software
piبرانemark.com

exocad DentalCAD® *exoplan*
BIEMME
BIOMAGING TECHNOLOGIES



● MT-F
Implants

Platform Ø	3.3	3.5	3.9	4.6
h				
18		172319		
15	172297	172302	172384	
13	172296	172301	172383	172306
11.5	172295	172300	172382	172305
10	172294	172299	172381	172304
8.5	172293	172298	172380	172303
7		172318	172379	172321
6		172317	172378	172320



Implant Ø

3.3

3.75

4.1

4.8



Short
Implants

(!) Same Interface in all Platform Ø including narrow and short implants.

D DIGITAL
prosthetics



◎ Components

Soft Tissue Healing

Cover Screw • Healing Abutment



Biological Width

- Concave or Parallel emergence Healing • Soft Tissue contouring
- Potential for more soft tissue volume
- Minimized cortical bone removal for sub-crestal Implants



One Cover Screw

- For all Implants and Platforms • MT Interface



(!) Healing Abutment Parallel h = 4.5 mm features height marks serving as a reference for Abutment selection.

Conical Abutment

Single or multiple, screw-retained prosthesis



Biological Width

Concave emergence • Increased soft tissue volume
Minimizes cortical bone removal to install Abutment

Multiple and Single prosthesis


Universal Ø 4.8 Multi Unit platform
Straight Conical Abutment has double indexation
Select engaging components for single unit

D DIGITAL

Provisional and CAD/CAM prosthesis
Adjustable post heights: 10, 9, 8, 7, 5.5 and 4.5mm



Direct Geometry for single and multiple unit • bridge



Cylinder



	Non-Engaging	
	CoCrMo	101141
	Castable	101143
	Titanium	101142
	Engaging	
	CoCrMo	171249
	Castable	171250
	Titanium	171248

Analog



Cast 171247

Impression Coping



	Non-Engaging	
	Open Tray	102385
	Closed Tray	101113
	Engaging	
	Open Tray	171245
	Closed Tray	171246

Healing Cap



101155

Abutment
Multiple



	30°	h			
		3	●	171129	●
		1.5	●	172500	●



	17°	4.5	●	172561	●
		3	●	171128	●
		1.5	●	171127	●

Abutment
Multiple and Single



	Straight	4.5	●	171126	●
		3	●	171125	●
		1.5	●	171124	●
		0.8	●	171123	●

 **Conical Abutment**

D See DIGITAL Prosthetics.

● Regular Abutment is used.

(!) Conical Abutment prosthetic Platform has Ø 4.8 mm.

Cylinder • Conical Abutment ENG (Engaging) option is available for straight Conical Abutment only.

● Abutment Cemented Cylinder

Single or multiple, cement-retained prosthesis



Biological Width

Concave emergence • Potential for more soft tissue volume
Minimizes cortical bone removal to install Abutment

Anterior and posterior

Indexed • 6 and 4mm post heights

Single and multiple Castable Cemented Cylinders ○ ●



One-time one-abutment option

Prosthetic procedures over Abutment or Implant Platform

Zero margin Abutment

For limited interproximal spaces “0”



● Contour & Esthetic Abutments

Single or multiple, cement-retained prosthesis



Increased Biological Width

Concave emergence • Potential for more soft tissue volume

Minimizes cortical bone removal to install Abutment

Robust design

Adjustable • Straight and 17°

Delicate slim profile

Adjustable • Straight and 15°

Impression at Implant Platform

Short and Long Impression Coping

Open and Closed Tray

Contour



Esthetic



Impression



N R W

Contour 17°



h	N	R	W
4.5	171116	171119	171122
3	171115	171118	171121
1.5	171114	171117	171120

Contour Straight



4.5	171107	171110	171113
3	171106	171109	171112
1.5	171105	171108	171111

Contour Abutment

Esthetic 15°



h	N	R	W
4.5	171178	171181	●
3	171177	171180	●
1.5	171176	171179	●

Esthetic Straight



4.5	171171	171175	●
3	171170	171174	●
1.5	171169	171173	●
0.8	171168	171172	●

Esthetic Abutment

Implant Analog



●	171212	●
---	--------	---

Impression Coping
Implant



Open Tray	●	171206	●
Closed Tray	●	171209	●
Open Tray, Long	●	172418	●
Closed Tray, Long	●	172417	●

Implant Impression

Cylinders over Implant

Single or multiple, cement or screw-retained prosthesis



E ENHANCED
Ceramic Surface



D DIGITAL
prosthetics

Provisional, conventional and CAD/CAM
Platform diameter and margin height options
Adjustable post heights: 6, 5 and 4mm

Increased Biological Width

Concave emergence • Potential for more soft tissue volume
Minimizes cortical bone removal to install Cylinder

Definitive • Overcasting

Main body CoCrMo and waxing sleeve POM with retentions



Provisional and Custom Healing • Titanium

Diameter options • Flat areas and trapezoidal retentions



Impression at Implant Platform

Short and Long, Open and Closed Tray, Impression Coping



		N	R	W	
Cylinder DIGITAL and Conventional		Non-Engaging			
		h			
		4.5	172534	172536	172541
		3	172533	172539	172544
		1.5	172532	172538	172543
		0.8	172531	172537	172542
		Engaging			
		4.5	172401	172405	172409
		3	172404	172408	172412
		1.5	172403	172407	172411
	0.8	172402	172406	172410	
Cylinder Conventional		Non-Engaging			
		Titanium	172530	172535	172540
		Engaging			
	CoCrMo	171183	171185	171187	
	Titanium	171182	171184	171186	

● Cylinder over Implant

Implant Analog			●	171212	●
Implant Impression Coping		Open Tray	●	171206	●
		Closed Tray	●	171209	●
		Open Tray, Long	●	172418	●
		Closed Tray Long	●	172417	●

Implant Impression

D See DIGITAL Prosthetics.

● Regular Components are used.

(!) Cylinder over Implant • Ti with height options can also be used to support prosthesis manufactured by CAD/CAM.



D DIGITAL prosthetics



Manufacturing
Prosthetics • 3D Model



Components
and Direct Geometries



Design
Libraries



Scan
Intraoral • Desk



Scan Body
Implant • Conical Abutment



D files for open
source software
pibranemark.com

3shape

exocad

(!) Libraries are available for upload at pibranemark.com and/or from CAD/CAM software system, except Dentsply Sirona. Always check latest library version and availability.

D Links

ENHANCED
Ceramic Surface

N **R** **W**

Link CAD | CAM



h	N	R	W
4.5	●	171137	●
3	●	171136	●
1.5	●	171135	●
0.8	●	171134	●

Link C



	Small	Large	
4.5	172422	172426	●
3	172421	172425	●
1.5	172420	172424	●
0.8	172419	172423	●

Link

Implant • MT
for cast or 3D models



●	172416	●
---	--------	---

Analog + Digital

Implant • MT



PEEK	●	161469	●
Titanium	●	161563	●

Scan Body

● MT interface is the same for all platforms. Scan Body and Analog + Digital are used.
(!) Links C 3.5 • Small and 4.1 • Large have similar post geometries to Dentsply Sirona, CEREC®.

D Cylinder over Implant

ENHANCED
Ceramic Surface

N

R

W

Cylinder
DIGITAL and
Conventional



	h		
Non-Engaging			
4.5	172534	172536	172541
3	172533	172539	172544
1.5	172532	172538	172543
0.8	172531	172537	172542



	h		
Engaging			
4.5	172401	172405	172409
3	172404	172408	172412
1.5	172403	172407	172411
0.8	172402	172406	172410

Cylinder over Implant

Implant • MT
for cast or 3D models



	●	172416	●
--	---	--------	---

Analog + Digital

Implant • MT



PEEK	●	161469	●
Titanium	●	161563	●

Scan Body

● MT interface is the same in all platforms. Scan Body and Analog + Digital are used.

D Cylinder • Conical Abutment

N R W

Titanium
DIGITAL and
Conventional



○ Non-Engaging

101142

● Engaging

171248

Cylinder

Conical Abutment
for cast or 3D models



172453

Analog + Digital

Conical Abutment



Non-Engaging



PEEK

161471

Titanium

161565

Engaging



PEEK

161556

Titanium

161566

Scan Body

Milling Blank



A - Arum



172443



Milling Blank

Analog Fixation



Conical Abutment Tool

70170003

Implant MT Tool

70170002

Screw

70160502

Digital Analog
with fixation tool



Conical abutment, NENG

70150007

Implant MT

70150005

Scan Body
Titanium



Conical abutment, NENG

80620277

Implant MT Interface



80620276



Third Party Accessories


● MT interface is the same in all platforms. Scan Body, Analog + Digital and Milling Blank are used.

■ Manufactured by CoreMedTech. Check latest CoreMedTech libraries version and availability.

(!) Cylinder • Conical Abutment ENG (Engaging) option is available for straight Conical Abutment only.

D Dimensions and References

Link



Cylinder over Implant



Cylinder • Ti • Conical Abutment



Direct Geometry



(!) Direct Geometries are used to manufacture direct screw-retained CAD/CAM prosthesis without P-I titanium base (Links or Cylinders), using spare P-I Screws. Always check latest library version and availability.

Kit

Conventional



width	254 mm
height	40
depth	130

Advanced
Stainless Steel

181036

(!) Surgical and prosthetic
(!) Reference number is for ordering purposes only. Instruments and Tray delivered separately.

Kit

Tray options



w	202
h	67
d	158

**Advanced
Polymer**
181022



w	120
h	40
d	80

**Advanced
Compact**
181023



w	120
h	40
d	80

**Prosthetic
Compact**
181029

- (!) Prosthetic Kit does not include Torque Wrench and Adapters.
- (!) Conventional Kit Composition at pibranemark.com
- (!) Reference number for ordering purposes only. Instruments and Tray delivered separately.

Instruments • Conventional

Drills

		Ø	
	Initial	2.2	141138
	Conical	2.8	141146
		3.4	141148
		3.8	141314
		4.6	141152
		4.8	141315
	Dense	3.3	141213
		3.75	141316
		4.0	141215
		4.8 5.0	141317

Implant Insertion Driver

		Medium	Long
	All Systems	131139	131140

Tools

	Guide Pin	2.2 2.8	131114
		2.8 3.8	131115
	Guide Pin	2.2 2.8 C	141535
		2.8 3.8 C	141536
	Drill Extension		131028
	Spade	1.5	141319 ▲
	Round Burr	1.8	141001 ▲

 Guided & conventional
SURGERY



Wash Tray • GS
Surgical + Prosthetic



Guided & Conventional

181040

All Instruments



Guided +

181041

Guided Surgery Instruments only
includes Prosthetic Instruments



Tray • GS

131162

Tray only

width	190
height	61
depth	138



Stop Organizer • GS

141399

Organizer [9] only

Instruments • Guided & Conventional


Drills		Narrow	Regular	
	Crestal Drill	141330		141335
← 27				
		∅ N • Long	R • Short	R • Long
	Initial	2.2 141378	141336	141340
		2.8 141381	141337	141341
	Conical	3.4	141338	141342
← 35 S 40 L		3.8	141339	141343
Stop • GS Offset [9]				
		h		
		6	141365	
		7	141366	
		8.5	141367	
		10	141360	141368
		11.5	141361	141369
		13	141362	141370
		15	141363	141371
				
		141375	141376	
Sleeve				
Implant Insertion Driver				
		131153	131152	
← 31				
	Insertion Driver must be completely connected into the Implant internal hexagonal indexation. When used manually or with Torque Wrench, the upper hexagonal portion of the Insertion Driver must be entirely connected to the hexagon of the Driver Adapter.			

Guide Fixation



Sleeve - Pin

141382



Pin Drill Ø1.7

141377



Fixation Pin

131167



Surgical Accessories

		Narrow	Regular
Punch			
	Soft Tissue Punch	141331	141332
Dense Drill			
	Dense Drill	141344	141345

Torque Wrench Kit



Prosthetic Instruments

	Hexagonal Ø 1.2	Short	131010
		Medium	131011
		Long	131012
	Conical Abutment Ø 2.0	Short	131016
		Medium	131017
	Handling Tool, Angled		102964
	Retriever MT	Short	141564
		Medium	131131

Accessories



	Depth Probe • GS [9] Offset tip + Round tip	141439	●
	Stop Removal Tool	131172	●

■ Manufactured by Elos MedTech Pinol A/S . Torque Wrench Kit includes Surgical & Prosthetic Adapters.
● Stop Removal Tool should be used when Surgical Adapter is not sufficient for the removal of Stop from Drill.
Depth Probe is not part of Kit • GS and should be ordered separately.

Stop • Conventional



+



R • Short



R • Long

h	R • Short	R • Long
6	141592	
7	141593	
8.5	141594	
10	141590	
11.5	141591	141584
13	141600	141586
15	141601	141586
18		141602

Trays



Polymer
131134



Stainless Steel
131117



Prosthetic
131135

Compact
131138

Lab & Replacement Conical Abutment



Screw Cylinder 101984



Screw Impression Coping
Open Tray 101737
Closed Tray 171260

Implant



Coated Screw MT 171239



Screw Impression Coping
Open Tray 171207

Resonance Frequency Analysis

Penguin^{RFA} ■ ⁽¹⁾



ISQ

Kit
55002-S ▲



MultiPeg

MT-F • Type 38

55046

Conical Abutment • Type 72

55080

Locator® Overdenture

Abutment ■ ⁽²⁾



N
R
W

h	●	●	●
4	●	2203	●
3	●	2202	●
2	●	2201	●
1	●	2200	●

■ (1) Manufactured by Integration Diagnostics Sweden AB .
▲ Penguin RFA Kit includes instrument, charger, MultiPeg driver and user's manual.
■ (2) Manufactured by Zest Dental . Components and instruments not included in the P-I Catalog.
● MT interface is the same in all platforms. Regular Abutments are used.
 P-I pegs are available from Osstell. MT-F Implant SmartPeg, type 21.

○ MT-F

SURGICAL SEQUENCE



- 
rpm
600 - 1,200
Lowest possible rpm
- 
In-Out
△ Coordinated in-and-out movement of Drills for better cooling
- 
ITV
≤ 70 Ncm
Insertion Torque Value
- 
Full Length
Prepare at planned full length of Implant position combining the Drills GS and Stops GS in accordance to planning
- 
Irrigation
Constant irrigation to the insertion margin of Conical Drills



△ Drills are less than 1 mm longer than Drill marks

SURGICAL SEQUENCE

 3.3	<table border="0" style="width: 100%; text-align: center;"> <tr><td>2.2</td><td>2.8</td><td>3.4</td></tr> <tr><td>✓</td><td>S ↓</td><td>D ↓</td></tr> <tr><td></td><td>N</td><td></td></tr> </table>	2.2	2.8	3.4	✓	S ↓	D ↓		N		 3.75	<table border="0" style="width: 100%; text-align: center;"> <tr><td>2.2</td><td>2.8</td><td>3.4</td></tr> <tr><td>✓</td><td>S</td><td>N ↓</td></tr> <tr><td></td><td></td><td>D</td></tr> </table>	2.2	2.8	3.4	✓	S	N ↓			D										
2.2	2.8	3.4																													
✓	S ↓	D ↓																													
	N																														
2.2	2.8	3.4																													
✓	S	N ↓																													
		D																													
 4.1	<table border="0" style="width: 100%; text-align: center;"> <tr><td>2.2</td><td>2.8</td><td>3.4</td><td>3.8</td></tr> <tr><td>✓</td><td>S</td><td>N</td><td>D</td></tr> <tr><td></td><td colspan="3">..... Optional></td></tr> <tr><td></td><td colspan="3">..... for Dense</td></tr> </table>	2.2	2.8	3.4	3.8	✓	S	N	D	 Optional>			 for Dense			 4.8	<table border="0" style="width: 100%; text-align: center;"> <tr><td>2.2</td><td>2.8</td><td>3.8</td><td>4.6</td></tr> <tr><td>✓</td><td>✓</td><td>S</td><td>N ↓</td></tr> <tr><td></td><td></td><td></td><td>D</td></tr> </table>	2.2	2.8	3.8	4.6	✓	✓	S	N ↓				D
2.2	2.8	3.4	3.8																												
✓	S	N	D																												
 Optional>																														
 for Dense																														
2.2	2.8	3.8	4.6																												
✓	✓	S	N ↓																												
			D																												

Bone Density

S Soft **N** Normal **D** Dense




6 mm
Drill tip

(!) The subsequent Conical Drill, in terms of diameter, should be considered with a drilling depth of 6 mm, in order to not exceed 70 Ncm of insertion torque value. The use of Dense Drills (15 – 50 rpm) can also be considered to lower the insertion torque value. Soft Tissue Punch (15 – 50 rpm).

📄 See Guided Surgery Surgical Sequence at pibranemark.com. Guided Surgery Drills consider a [9mm] offset and, when used with Drill stops, allow for limiting the total length of osteotomy with the objective of providing predetermined Drill length and orientation through the surgical guide. Height repositioning for Sleeve and Stop selection required for [10.5] Offset.

One interface

for all Implant dimensions



Dimensions



(!) Images are for illustrative purposes only. Measurements in millimeters.

One Driver

Manual • Torque Wrench



1.2



Torques

Ncm

MT-F Implants

≤ 70

Abutments

Cylinders over Implant

25

Links, Milling Blank

Cylinders • Conical Abutment

15

Cover Screw


Healing Abutments

Impression Copings

Manual

Scan Bodies

(!) Recommended Torques. Abutment and Components torques should not exceed the torque obtained at Implant installation.
(!) One Prosthetic Driver, except Locator® and straight Conical Abutment.



LifeTime GUARANTEE





Developed By P-I Brånemark

MT-F SMART GUIDE-009-ENG-2011.2024
P-I Brånemark. All rights reserved.

pibranemark.com

 SIC invent Sweden AB
Stora Åvägen 21
43634 Askim, Sweden

 2460 