Aquasil® Ultra+

Aquasil® Ultra+ XLV • Aquasil® Ultra+ LV • Aquasil® Ultra+ Medium • Aquasil® Ultra+ Heavy • Aquasil® Ultra+ Rigid

Smart Wetting® Impression Material

INSTRUCTIONS FOR USE ENGLISH 2

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Distributed by
Dentsply Canada
525 Rowntree Dairy Rd
Woodbridge, Ontario
L4L 7K6 Canada

Aquasil® Ultra+

Smart Wetting® Impression Material

Aquasil® Ultra+ XLV • Aquasil® Ultra+ LV • Aquasil® Ultra+ Medium • Aquasil® Ultra+ Heavy • Aquasil® Ultra+ Rigid

INSTRUCTIONS FOR USE – ENGLISH

Caution: This is a medical device. For dental use only by dental professionals.

USA: Rx only.

1. PRODUCT DESCRIPTION

Aquasil® Ultra+ Smart Wetting® Impression Material is a quadrafunctional hydrophilic addition reaction silicone elastomeric dental impression material with excellent hydrophilic properties, dimensional accuracy, high tear strength, and resistance to permanent deformation.

Aquasil® Ultra+ Material complies with ISO 4823:2021 Elastomeric impression materials.

1.1 Indications For Use

Aquasil® Ultra+ Impression Material is indicated for all dental impression techniques.

1.2 Contraindications

Aquasil® Ultra+ Impression Material is contraindicated for use with patients who have a known hypersensitivity to peppermint oil.

1.3 Delivery Forms (Some delivery forms may not be available in all countries)

Aquasil® Ultra+ Material is available in:

- Light-bodied consistency, Medium-bodied consistency, and Heavy-bodied consistency
- digit® Targeted Delivery System Unit Dose Dispensed delivery, digit Power® Dispensed delivery, 50 ml cartridge delivery, and 5:1 DECA™ delivery for Dynamic Mixing Machines
- Fast Set, Regular Set and Extended Set

1.4 Composition

Base Pastes: Vinyl Terminated Polydimethylsiloxane Polymer; Diatomaceous Earth (Silicon Dioxide – Crystalline and Amorphous); Polymethylhydrogen Siloxane; Silicone Surfactant; Sodium Aluminosilicate; Various pigments: (F D & C Blue #1, Day Glo Pigments, Ultramarine Pigments, Iron Oxides & Titanium Dioxide); Peppermint Oil

Catalyst Pastes: Vinyl Terminated Polydimethylsiloxane Polymer; Diatomaceous Earth (Silicon Dioxide – Crystalline and Amorphous); Titanium Dioxide; Organo Platinum Complex; Sodium Aluminosilicate; Peppermint Oil

1.5 Technical Specifications

- Classification: ISO 4823

- Type 3: Light-bodied consistency: Aquasil® Ultra+ XLV, Aquasil® Ultra+ LV
- Type 2: Medium-bodied consistency: Aquasil® Ultra+ Medium, Aquasil® Ultra+ Heavy
- Type 1: Heavy-bodied consistency: Aquasil® Ultra+ Rigid

- Working Time (22 °C/72 °F)..... **Fast Set:** 1'00" minimum
Regular Set: 2'15" minimum
Extended Set: 2'15" minimum
(room temperature mixing and preparation, maximum insertion time)
- Working Time (37 °C/98 °F)..... **Fast Set:** 35" Maximum
Regular Set: 1'10" Maximum
Extended Set: 1'45" Maximum
(intraoral syringing prior to tray insertion)
- Setting Time..... **Fast Set:** 2'30" from start of mix
Regular Set: 5'00" from start of mix
Extended Set: 5'30" from start of mix
(Mouth Removal Time)
- % Linear Dimensional Change..... < 0,50
- % Recovery from deformation ≥ 98
- Detail Reproduction..... 20 microns

1.6 Compatible Tray Adhesive

Aquasil® Ultra+ Impression Material is compatible with conventional tray adhesives designed for use with vinyl polysiloxane materials (see complete Instructions for Use).

1.7 Compatible Die Stones

- The material is compatible with a range of dental plasters and die stones including Glastone® Type III dental stone and Glastone® 3000 Type V high strength die stone.
- Following impression casting, reprocess reusable impression trays according to tray manufacturer's instructions.

2. GENERAL SAFETY NOTES

Be aware of the following general safety notes and the special safety notes in other chapters of these Instructions for Use.

Safety alert symbol.

This is the safety alert symbol. It is used to alert you to potential personal injury hazards. Obey all safety messages that follow this symbol to avoid possible injury.

2.1 Warnings

- Aquasil® Ultra+ Material contains polymerizable monomers and surfactants which may be irritating to skin, eyes and oral mucosa and may cause allergic contact dermatitis in susceptible persons.
 - **Avoid Eye contact** to prevent irritation and possible corneal damage. In case of contact with eyes, rinse immediately with plenty of water and seek medical attention.
 - **Avoid Skin contact** to prevent irritation and possible allergic response. In case of contact, reddish rashes may be seen on the skin. If contact with skin occurs immediately remove material with cotton and wash thoroughly with soap and water. In case of skin sensitization or rash, discontinue use and seek medical attention.
 - **Avoid ingestion/swallowing of material** to prevent irritation or obstruction. If accidental swallowing occurs, drink lots of water. This material is not hazardous when small quantities are ingested. Larger quantities may cause bowel obstruction. Seek medical attention in the event of digestive irregularities.

- The product is intended for transient contact. Do not use in applications other than the intended use (e.g., do not use as a temporary reliner) because it may cause irritation.
- Avoid multiple impressions at one sitting (more than three of the same arch/area) to prevent irritation to the oral mucosa. If symptoms occur discontinue use immediately. If irritation does not subside consult a physician.
- Aquasil® Ultra+ Material should not be used with patients who have a known hypersensitivity to any of the components.
- Do not use Aquasil® Ultra+ Rigid in combination with full arch impression trays for the double mix technique. The low flexibility of the set material may interfere with tray removal.
- Due to the stiffness of Aquasil® Ultra+ Rigid Material, severe undercuts should be blocked out prior to tray insertion. For full arch dentate impressions, select other Aquasil® Ultra+ Material viscosities for tray material.
- Avoid use on periodontally involved teeth or teeth with compromised attachments, severe undercuts or open furcations to prevent injury to periodontal tissues or tooth extraction. Do not use excessive force to remove impression. If necessary, section tray and material for removal.
- To avoid periodontal tissue harm, check the integrity of the impression and check that there is no residue of material to be removed from the patient's mouth. Verify that there are no residues of material in the gingival area.
- Use digit cartridges and mixing/intraoral tips are contaminated. Follow proper exposure control plans. In the event of an accidental exposure, follow recommended post-exposure procedures. Following intraoral use of the digit cartridge/mixing tip, proper handling and disposal techniques are required. Used intraoral tips are contaminated and sharp. When handling or disposing, use precautions as when handling or disposing of other contaminated sharps.

2.2 Precautions

- This product is intended to be used only as specifically outlined in the Instructions for Use. Any use of this product inconsistent with the Instructions for Use is at the discretion and sole responsibility of the practitioner.
- Insufficient data exist to support the use of Aquasil® Ultra+ Material for impressions on exposed bone. Aquasil® Ultra+ Material should only be used on or around intact mucosa.
- Wear suitable protective eyewear, mask, clothing and gloves. Protective eyewear is recommended for patients.
- Devices marked "single use" on the labeling are intended for single use only. Discard after use. Do not reuse in other patients in order to prevent cross-contamination.
- Use only clean and properly disinfected and sterilized, or ready to use accessories. Do not use dispensers which have not been properly pre-use processed or reprocessed. Before use and each re-use dispensers must be properly cleaned and autoclaved. Before storage after use/re-use, cartridges must be properly cleaned and disinfected to reduce the risk of cross contamination (see Hygiene, Processing and Disposal below).
 - As additional precautionary measure, dispensers may be protected from gross debris but not from all contamination by applying a protective barrier.
- digit Cartridges, Mix Tips and Intraoral Tips are provided ready for use and should not be processed. Use only new, clean unused tips, stored in covered storage, such as drawers and cabinets, away from potential contamination, and removed from their original packaging using either clean gloved hands or a clean suitable instrument at time of use. To reduce risk of infection, do not use tips if packaging is damaged or compromised.
- All products should be used at room temperature. Higher temperatures reduce work times and laboratory bench set times (faster) lower temperatures increase them (slower). Clinically, intraoral set time is unaffected within the range of allowable storage and initial mixing temperatures. The material sets faster in the mouth than on the bench. Check impression in the mouth (not on the bench) to be sure it is resilient before removal.
- Do not use if the mix is not homogeneous. To obtain optimum physical properties, the mix should be completely homogeneous (streak free). Base and Catalyst are contrasting colors. When dispensing, observe both components entering the mixing tip, and exiting as a homogeneous mix. If not, remove the dispenser from the operating field, discard the mixing tip and ensure flow from both ports. Replace mixing tip or cartridge and tip and return to dispensing.
- Aquasil® Ultra+ Material should be used with properly fitted, rigid trays with sides for most applications. Trays with retentive features such as Caulk Rim-Lock® Trays or perforations, pretreated with tray adhesive, are recommended. To reduce gagging, choking or swallowing, do not overfill tray.
- Material should extrude easily: DO NOT USE EXCESSIVE FORCE. Excessive pressure may result in unanticipated extrusion of the material or cause cartridge rupture. If the material doesn't extrude easily remove from the operating field and discard the mixing tip, level component (bleed cartridge) and install a new mixing tip.
- Safety and efficacy of Aquasil® Ultra+ Heavy, Aquasil® Ultra+ Rigid Material when used as an intraorally syringed (wash) material has not been demonstrated.
- Tip (hubs) and cartridge caps and collars are color coded. Use only blue tips on high viscosity cartridge materials. Use of red (for low viscosity wash) tips on blue (higher viscosity) cartridges may result in excessive extrusion force, slowed extrusion speed, incomplete mix and work/set time irregularities.
- Check the integrity of the impression and check that there is no residue of material to be removed from the patient's mouth.
- If used with a syringe refer to the manufacturer instructions to check compatibility of the syringe with VPS before. Do not use a syringe that has been used with another chemistry.
- Do not use the product if package is broken.
- Suitable cartridges for use with digit Power® Dispenser will be marked PWR on tab. Do not use digit Power® Dispenser with materials in digit® Targeted Delivery System cartridges designed for manual dispensing. Device and cartridge malfunction may occur.
- Do not reuse the original cartridge cap. Store the cartridge with the disinfected used mixing tip.
- Two-step technique is NOT recommended for double arch trays.
- The components of the device (base and catalyst) must be used in equal parts to obtain the right setting times and not to alter the performance of the product; a higher dose of catalyst does not accelerate setting times.

Interactions:

- Use of latex or "sulfur-containing" polymeric gloves, some hand soaps and lotions can lead to interference with setting reaction. Contamination may accumulate on the surface of the skin and may not easily wash off. Non-latex gloves worn during use will avoid this possibility.
- Astringents based on aluminum or iron salts, or Hydrogen peroxide-containing dental materials may interfere with the setting reaction of vinyl polysiloxanes. Rinse carefully with water before applying the material. Please consult Instructions for Use of astringent to be used.
- Handle cords with college pliers. Do not touch prepared tooth with gloves. Rinse preparation with copious amounts of water before making the impression.
- Aquasil® Ultra+ Material should not be intermixed with, or used in conjunction with, other manufacturers' vinyl polysiloxanes or any polyether, polysulfide or conventional (condensation cured) silicones.
- The dual phase, two-step ("putty/wash") crown and bridge impression technique is only recommended with Aquasil® Putty (all versions) or Reprosil® impression materials (available separately) for primary tray impressions.
- When used as a primary impression, material should be covered with a polyethylene spacer sheet, and protected from direct contact with oral tissue.

2.3 Adverse Reactions

- 1. Product may irritate the eyes and skin.
 - **Eye contact:** irritation and possible corneal damage.
 - **Skin contact:** irritation or possible allergic response. Reddish rashes may be seen on the skin.
 - **Mucous Membranes:** inflammation
- 2. Bowel obstruction or other digestive distress may result from ingestion of mixed impression material (see Warnings).
- 3. Residual material may lead to soft tissue injury or irritation (see Warnings and Step-by-Step instructions).
- 4. Long term or repeated contact may result in tissue irritation (see Warnings).

2.4 Storage Conditions

Inadequate storage conditions may shorten the shelf life and may lead to malfunction of the product.

- Store at temperatures between 10 °C-24 °C/50 °F-75 °F.
- Use the product at room temperature.
- Do not store used cartridges with new (unused) mix tips attached.
- Keep out of sunlight.
- Do not freeze.
- Protect from moisture.
- Do not use after expiration date.
- Do not use if the lot number and/or expiration date are missing or illegible.
- Work time may be affected by temperature. See Step-by-Step Instructions for Use.

3. STEP-BY-STEP INSTRUCTIONS

3.1 Cartridge Dispensing Gun

- 3.1.1 Insert cartridge into dispensing gun.

Technique Tip: Cartridge is compatible with a standard 50 ml cartridge dispenser. No special additional dispenser is required.
- 3.1.2 Break safety seal by unscrewing the cap. Remove the colored protective cap from the cartridge and discard.
- 3.1.3 Seat the plunger snugly against the cartridge.
- 3.1.4 For the first use, squeeze the trigger on the dispenser to bleed a small amount of material from the cartridge, ensuring that both sides are flowing evenly. Discard the bled material.

Technique Tip: Bleeding is not necessary prior to installing a new mixing tip for next use.
- 3.1.5 Affix the appropriate mixing tip onto the cartridge. Align the alignment key with the keyhole then screw collar tightly down about half a turn until it locks.

Technique Tip: Use only mixing tips supplied by Dentsply Sirona for this cartridge. Mixing tips from other manufacturers will not fit this advanced cartridge design. Attach the intraoral tip to the RED TIP. Intraoral tips are not recommended for material in blue-capped cartridge/blue mixing tips.
- 3.1.6 Before use on each new patient, always apply a disposable cover/barrier system to the cartridge, to be removed after each use. The disposable cover/ barrier system does not replace the recommended cleaning and disinfection operations.

3.2 Dispensing Material

- 3.2.1 Squeeze dispenser trigger with moderate and even pressure to begin mixing material. Material will stop flowing after trigger is released.
 - 3.2.2 Preparation of syringe/wash material:
 - Dispense directly to the mouth with the intraoral tip; OR
 - An intraoral syringe may be filled with mixed material from the cartridge; OR
 - Prepare digit Power® Dispenser or digit® Targeted Delivery System Unit Dose wash/syringe material according to the Instructions for Use.
- NOTE:** Do not inject Aquasil® Ultra+ Heavy or Aquasil® Ultra+ Rigid Impression Material intraorally.
- 3.2.3 Allow material to harden in the mixing tip to serve as a cap on the cartridge until the next use.
 - 3.2.4 Sterilize cartridge dispensing gun, disinfect cartridge and attached used mixing tip for storage according to instructions (see Hygiene section).
 - 3.2.5 Immediately prior to subsequent usage, remove the used mixing tip, wipe end of cartridge, and install a new mixing tip. No cartridge bleeding is necessary.

To reduce risk of excessive extrusion force or incomplete mix.
Do not dispense higher viscosity tray material through an attached red (small) mixing tip. Excessive force, cartridge rupture and/or mixing and discrepancies may occur.

- 3.2.6 To remove the cartridge, hold the release lever vertically while simultaneously withdrawing plunger in a backward direction. Dispose of cartridge when empty. Do not attempt to refill or consolidate used cartridges.

3.3 digit® Targeted Delivery System Unit Dose

- 3.3.1 Insert plunger into open end of syringe barrel.
- 3.3.2 Snap circular disk off of end of cartridge.
- 3.3.3 Insert unit dose cartridge into mix tip. Ensure both cartridge spurs fully lock into mix tip slots. Attach intraoral tip.
- 3.3.4 Insert assembled cartridge/mix tip (with intraoral tip) into the syringe side opening. Slide mix tip through barrel end with cross notches. Pull on mix tip to ensure complete seating.
- 3.3.5 Dispense the material by applying slow, steady pressure on plunger. Do not use excessive force.

Material extrusion – To reduce risk of injury due to excessive force.

- Apply slow and steady pressure on the plunger.
- Do not use excessive force – cartridge rupture may result.

- 3.3.6 Bleed and discard a small amount of mixed material, then proceed immediately to clinical application.

To reduce risk of incomplete mix – Tissue irritation or compromised properties.

- Observe both base and catalyst entering mix tip.
- If not, remove cartridge/mix tip assembly and replace.

- 3.3.7 To disassemble, withdraw the plunger, tap mix tip against the counter to dislodge, and remove through side opening.
- 3.3.8 Properly discard used cartridge/mix tip and prepare digit Syringe Dispenser for subsequent reuse by following instructions (see Hygiene below) for cleaning and sterilization.

3.4 digit Power® Dispenser (see complete Instructions for Use)

- 3.4.1 Power handpiece on.
- 3.4.2 Select pre-set timer.
- 3.4.3 Assemble digit Power® cartridge. Ensure selected cartridge is marked PWR on tab.
- 3.4.4 Assemble cartridge cap.
- 3.4.5 Attach cartridge cap assembly to handpiece. Ensure cap is attached securely.
- 3.4.6 Place a protective barrier sleeve over Dispenser.
- 3.4.7 Bleed cartridge, then dispense material.
- 3.4.8 Properly discard used cartridge/mix tip and reprocess according to the Instructions for Use.

3.5 DECA™ Cartridge Loading

NOTE: Heavy Cartridge – Dropping may result in damage to the cartridge rendering it unusable. Never use a damaged cartridge.
The DECA™ cartridge fits duomix™, duomix II, Pentamix®, Dynamix and Volume Mix dynamic mixers*. Use with MixStar dynamic mixer is not recommended. Cartridges are not compatible with 5:1 Pentamix Lite. Check compatibility with the manufacturer instructions.

Carefully remove dispensing cartridge from the package. Also refer to separate dynamic mixer manufacturer's instructions.

- 3.5.1 Remove the cap from the cartridge by gently turning it counterclockwise and lifting the cap from the cartridge.
- 3.5.2 Affix the dynamic mixing tip by turning it clockwise to the cartridge. Ensure that the alignment key of the dynamic mixing tip slides into the port of the cartridge. Turn the mix tip until tip is locked in place. No locking ring (or bayonet ring) is required to retain the dynamic mixing tip. Use only Dentsply Sirona mix tips. Other tips may not seal properly.
- 3.5.3 Insert the cartridge into the dynamic mixer.
- 3.5.4 Engage machine plunger until the machine's drive shaft is engaged in the mixer drive socket (mixing element is turning). Do not extrude material until ready for use.
- 3.5.5 When ready to dispense, activate the machine, dispense, and discard initial mixed material to ensure homogeneous mix.

NOTE: For first use, always dispense and discard approximately 75 mm (approximately 3 inches) of material before filling tray. For all subsequent uses, always dispense and discard approximately 30 mm of material before filling tray.

To reduce risk of incomplete mix, work time or set time discrepancies.

- For first use, always dispense approximately 75 mm of material through the mixing tip and discard.
- For subsequent uses, always dispense approximately 30 mm of material through the mixing tip and discard.
- Material must be bled through the mixing tip (not the cartridge before tip attachment) before loading tray.
- Once dispensing and bleeding through the mixing tip is completed, immediately load tray. Observe material work time starting from cartridge bleeding.

- 3.5.6 Immediately after activation and bleeding material through the tip, proceed to dispensing material to tray as desired.
- 3.5.7 Disinfect cartridge with attached used mixing tip. It will work as cap for storage.
- 3.5.8 Immediately prior to subsequent usage, remove the used mixing tip, wipe end of cartridge, install a new mixing tip, bleed and discard material through the tip and proceed to dispensing.

3.6 Impression Technique

Handling

- Use only clean and properly disinfected and sterilized, or ready to use accessories.
- Single or Dual Phase one-step ("double-mix") technique with a rigid tray or with double arch technique

NOTE: Tray and wash material mix timing – It is highly recommended to utilize an assistant to assure simultaneous syringing and tray loading.

- 3.6.1 Assure hemostasis and clean the field with air water spray.
- 3.6.2 Syringe material: Backfill the impression syringe or attach an intraoral tip to the end of the mix tip or prepare digit® cartridge syringe/Dispenser.
- 3.6.3(a) Tray material: Dispense material (assistant) into impression tray. Do not overfill tray.
- 3.6.3(b) Simultaneous with tray loading, syringe (dentist) completely around the preparation/aboutment.

Do not use if the mix is not homogeneous. To obtain optimum physical properties, the mix should be completely homogeneous (streak free). Base and Catalyst are contrasting colors. When dispensing, observe both components entering the mix tip, and exiting as a homogeneous mix. If not, remove the dispenser from the operating field, discard the mixing tip and ensure flow from both ports. Replace mix tip or cartridge and tip and return to dispensing.

To reduce risk of incomplete mix – Tissue irritation or compromised properties.

- Observe both base and catalyst entering mix tip.
- If not, remove tip, bleed cartridge and replace tip or replace cartridge and tip.

NOTE: Tray and wash material insertion timing.

- **Fast Set** wash and tray material: Loaded tray must be inserted within 35 seconds from time of first syringing of wash material intraorally.
- **Regular Set** wash and tray material: Loaded tray must be inserted within 1 minute 10 seconds from time of first syringing of wash material intraorally.
- **Extended Set** wash and tray material: Loaded tray must be inserted within 1 minute 45 seconds from time of first syringing of wash material intraorally.

- 3.6.4 Insert filled tray into the mouth.
- 3.6.5 Retain impression in position until firmly set. The material sets faster in the mouth than on the bench. Check impression in the mouth (not on the bench) to be sure it is firm, resilient and non-tacky before removal.

Minimum Removal Time (from the start of mix)

Fast Set: 2 minutes 30 seconds

Regular Set: 5 minutes

Extended Set: 5 minutes 30 seconds

NOTE: If Regular Set Wash is used with Fast Set Tray material, Minimum Removal Time is 4 minutes from start of mix.

Timing Guide	Fast Set: up to 0:35 (35 sec)	Fast Set: 2:30 (2 min 30 sec)
0:00 (0 sec)	Regular Set: up to 1:10 (1 min 10 sec)	Regular Set: 5:00 (5 min)
	Extended Set: up to 1:45 (1 min 45 sec)	Extended Set: 5:30 (5 min 30 sec)
Begin simultaneous syringing and tray loading	Complete syringing and seat loaded tray	Remove set impression

* Pentamix®, Dynamix and Volume Mix dynamic mixers are not registered trademarks of Dentsply Sirona.

- 3.6.6 Carefully release the impression and remove.
- 3.6.7 Check for residual material in undercut areas.

- To reduce risk of residual material – Periodontal irritation.**
- undercuts should be blocked out prior to use.
 - check the integrity of the impression.
 - Verify that there is no material residue.

3.6.8 Rinse impression thoroughly under running water. Follow Hygiene processing procedures below before shipping and/or use.

3.7 Dual Phase Two-Step (“Putty-Wash”) Technique With A Primary Tray Impression

Complete primary impression with Aquasil® Putty or Reprosil® Tray material.

- 3.7.1 Place one or more polyethylene plastic sheets over the impression material surface. This sheet creates the uniform spacing required for this technique. It also protects the surface from contamination.
- 3.7.2 Insert filled tray into the mouth. Move tray slightly from side to side and back and forth. This creates the spacing for the wash material.
- 3.7.3 Remove tray from the mouth when material has reached a rubbery consistency and set aside. Do not disturb impression until fully set.
- 3.7.4 Remove polyethylene sheet, rinse impression with water and air dry.
- 3.7.5 It is recommended to utilize an assistant and two cartridge guns, or one cartridge gun and digit® Targeted Delivery System Unit Dose syringe or digit Power® Dispenser to assure simultaneous syringing and tray loading.
- 3.7.6 Follow impression steps above, using wash material for both intraoral syringing and loading into primary tray impression.

3.8 Impression Technique: Edentulous/Overdenture Impression Technique

Aquasil® Ultra+ Medium, Aquasil® Ultra+ Heavy, and Aquasil® Ultra+ Rigid are acceptable for border molding.

- 3.8.1 Border Molding Technique.
- Dispense material on to the borders of the tray.
 - Insert the tray and perform conventional border molding techniques.
 - When material has reached a firm consistency, remove tray and examine borders. The addition of a small mix can be used if voids or deficient areas are present. Over-extensions can be easily removed with a scalpel blade.
- 3.8.2 The final impression can be made with Aquasil® Ultra+ XLV, LV, Medium, or Heavy Material.
- Insert filled tray into the mouth.
 - Observe the Minimum Removal Time of the selected Aquasil® Ultra+ impression material.
 - Rinse impression thoroughly under running water. Follow Hygiene processing procedures below before shipping and/or use.

3.9 Impression Technique: Duplication Of Models

- 3.9.1 Load the impression tray with Aquasil® Ultra+ tray material.
- 3.9.2 It is suggested to simultaneously syringe, using Aquasil® Ultra+ wash material, completely around the model teeth.
- 3.9.3 Seat model into the loaded tray. Retain impression in position until firmly set.
- 3.9.4 Remove model from impression by pulling slowly to break seal.

4. HYGIENE, PROCESSING AND DISPOSAL

Cross-contamination

- Do not reuse single use products. Dispose of in accordance with local regulations.
- Protective barriers are designed for single use and must be disposed of after each use in accordance with local regulations. The barrier is not a replacement for applicable cleaning, disinfection, and sterilization.
- Do not attempt to refill or consolidate used cartridges. Dispose of empty cartridges in accordance with local regulations.
- Reprocess reusable products as described below.
- Used cartridge and mix tips are contaminated and must be processed as described below before storage. Keep used, processed mix tip on cartridge as a cap.
- Process Impression as described below before shipping and/or use.
- Refer to dynamic mixing machine manufacturer’s Instructions for Use for dispenser cleaning and disinfection reprocessing procedures.
- When not in use, all supplies and instruments should be in covered storage, such as drawers and cabinets, and away from potential contamination.
- Dispensing materials with clean/disinfected gloves or hands or suitable transfer forceps in a separate room, bringing only what will be used into the operator room is strongly recommended.

4.1 50 ml Cartridge Dispenser

Process		✓ Acceptable
Cleaning	Mechanical (Ultrasonic)	
	Manual	✓
	Automated (Washer/Disinfector)	
Disinfection	Automated (Washer/Disinfector)	
	Manual	✓
	Immersion	
Sterilization	Steam Autoclave	✓
	Immersion	
Single Use	Dispose of after Use	
Reusable/Not Reprocessable	Dispose of if contaminated	

4.2 50 ml Cartridge/Mix Tip

Process		✓ Acceptable
Cleaning	Mechanical (Ultrasonic)	
	Manual	✓
	Automated (Washer/Disinfector)	
Disinfection	Automated (Washer/Disinfector)	
	Manual	✓
	Immersion	
Sterilization	Steam Autoclave	
	Immersion	
Single Use	Dispose of after Use	✓
Reusable/Not Reprocessable	Dispose of if contaminated	

4.3 digit Power® Cartridge Cap, and digit® Syringe Dispenser

Process		✓ Acceptable
Cleaning	Mechanical (Ultrasonic)	✓
	Manual	✓
	Automated (Washer/Disinfector)	✓
Disinfection	Automated (Washer/Disinfector)	✓
	Manual	✓
	Immersion	
Sterilization	Steam Autoclave	✓
	Immersion	
Single Use	Dispose of after Use	✓
Reusable/Not Reprocessable	Dispose of if contaminated	

4.4 VPS Bite Registration Impression

Process		✓ Acceptable
Cleaning	Mechanical (Ultrasonic)	
	Manual	✓
	Automated (Washer/Disinfector)	
Disinfection	Automated (Washer/Disinfector)	
	Manual	
	Immersion	✓
Sterilization	Steam Autoclave	
	Immersion	
Single Use	Dispose of after Use	
Reusable/Not Reprocessable	Dispose of if contaminated	

For complete processing instructions and validated parameters, please refer to the “50 ml Dispenser-Cartridge-VPS Bite Registration/Impression Hygiene” Instructions, “digit® Targeted Delivery Syringe Dispenser” Instructions and/or “digit Power® Dispenser” Instructions, which are available on our webpage at www.dentsplysirona.com/en/service-contact/ifu.html. If requested, we will send you a free printed copy of the Instructions for Use in the language you require within 7 days. In the USA, call 1-800-532-2855. Outside the USA, use the order form provided on www.dentsplysirona.com/en/service-contact/ifu.html#orderform-ifu for this purpose.

Technique Tip: Unset Material can be removed from clothing with dry cleaning solvent.

4.5 Casting

- The impression should be removed from the disinfectant, rinsed with water, and lightly air dried.
- The impression may be poured in plaster or stone 30 minutes after disinfection or delayed up to fourteen days. For epoxy dies, delay pour 60 minutes.
- The material is compatible with a range of dental plasters and die stones including Glastone® Type III dental stone and Glastone® 3000 Type V high strength die stone. Impressions may be silver or copper-plated.
- Store the impressions at ambient temperature, in a dry place away from direct sunlight.
- If the impression is to be shipped, use suitable packaging to prevent distortion.
- Impression scanning can be performed without the use of opaque sprays using specific scanners. Follow the instructions of the scanner manufacturer

4.6 Disposal

- Based on the information in the SDS, the waste generated by the disposal of the product does not present characteristics of physical danger, for the person or the environment; it can therefore be managed, according to local regulations, as a non-hazardous waste.
- Contaminated containers, accessories and finished impressions that have been in contact with the patient must be cleaned and disinfected as outlined in Hygiene above prior to disposal or disposed of as hazardous waste at risk of biological contamination.
- Always wear gloves when handling packaging, accessories and finished impressions.
- Used digit cartridge mixing tip/intraoral tip is contaminated and sharp. When handling or disposing, use precautions as when handling or disposing of other contaminated sharps.

5. LOT NUMBER, EXPIRATION DATE AND CORRESPONDENCE

1. Do not use after expiration date. Do not use if the lot number and/or expiration date are missing or illegible. ISO standard uses: “YYYY-MM-DD.”
2. The following numbers should be quoted in all correspondences:
 - Reorder number
 - Lot number
 - Expiration date
3. Any serious incident in relation to the product should be reported to the manufacturer and the competent authority according to local regulations.
4. Device Identification (Basic UDI-DI): ++D002IMPRESSIONMATSGC

© 2023 Dentsply Sirona Inc. All Rights Reserved. **578014WEB** (R 4/20/23)

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