Neo Spectra® ST

Neo Spectra® ST Effects

Universal Composite Restorative

Composite universel de restauration

Composite restaurador universal

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Neo Spectra® ST Effects

Universal Composite Restorative

CAUTION: For dental use by dental professionals only.

USA: Rx only.

Content	Page
1 Product description.....	2
2 Safety notes	4
3 Step-by-step instructions	5
4 Hygiene.....	10
5 Lot number, expiration date and correspondence	13

1 Product description

Neo Spectra® ST and **Neo Spectra® ST Effects** are nano-ceramic, light-cured, radiopaque, universal composites based on the novel **SphereTEC® Filler Technology** indicated for both, direct and indirect restorations.

The spherical, granulated filler technology called SphereTEC® leads to a ball-bearing effect that takes the handling properties of Neo Spectra® ST and Neo Spectra® ST Effects to the next level:

- **Neo Spectra® ST and Neo Spectra® ST Effects easily adapt to cavity surfaces but do not stick to hand instruments.**
- **Neo Spectra® ST and Neo Spectra® ST Effects are easy to sculpt precisely and have a high slump resistance.**
- **Neo Spectra® ST and Neo Spectra® ST Effects provide a firm consistency but are at the same time easy to extrude from Compules® Tips.**

Neo Spectra® ST comprises five universal CLOUD™ shades of intermediate translucency named A1 to A4.

Each universal CLOUD™ shade is designed to achieve shade match with a range of VITA®¹ shades. As a whole, the five universal CLOUD™ shades cover all 16 VITA® classic shades and allow for esthetic clinical results. Additionally, one bleach shade (BW) is available for bleached teeth.

Neo Spectra® ST Effects offers 2 opaque dentin shades (D1 and D3) and 1 translucent enamel shade (E1) that work with the **Neo Spectra® ST** universal CLOUD™ shades in a simplified layering system to provide added opacity and/or translucency as needed in esthetically demanding cases.

The spherical, granulated fillers consisting of submicron primary glass fillers, provide a fast and easy polishing procedure resulting in an outstanding gloss.

¹ Not a registered trademark of Dentsply Sirona, Inc.

1.1 Indications

- Direct restorations of all cavity classes in anterior and posterior teeth.
- Indirect fabrication of inlays, onlays and veneers.

1.2 Contraindications

- Use in patients with a history of severe allergic reaction to methacrylate resins or any of the other components.

1.3 Delivery forms

Neo Spectra® ST Universal Composite Restorative is available as HV (High Viscosity) and LV (Low Viscosity) for preferred handling. Neo Spectra® ST Effects is available as HV (High-Viscosity) only. Both Neo Spectra® ST and Neo Spectra® ST Effects offer two delivery options:

- Predosed Compules® Tips
- Syringes

1.4 Composition

- Barium-aluminum-borosilicate glass
- Pre-polymerized filler
- TEGDMA dimethacrylate multifunctional polymethacrylate
- Ytterbium fluoride
- Methacrylate modified polysiloxane (organically modified ceramic)
- Ethyl-4(dimethylamino)benzoate (photoaccelerator)
- Initiator
- Butylated hydroxy toluene (BHT)
- UV stabilizer
- Camphorquinone (photoinitiator)
- Titanium dioxide
- Iron oxide pigments
- Fluorescent agent

Filler system based on SphereTEC®.

The filler system comprises barium glass, pre-polymerized filler and ytterbium fluoride:

HV: 78-80 weight-% or 60-62 volume-%.

LV: 76-78 weight-% or 57-60 volume-%.

Additional ingredients are catalysts, additives, stabilizers and pigments

(~1 weight-%, according to shade). The total amount of inorganic fillers is:

HV: 72-74 weight-% or 50-52 volume-%.

LV: 71-72 weight-% or 48-50 volume-%.

The particle size of the inorganic fillers ranges from 0.1 to 3.0 µm.

1.5 Compatible adhesives

Neo Spectra® ST and Neo Spectra® ST Effects are used following application of a suitable dentin/enamel adhesive and are chemically compatible with conventional (meth)acrylate-based dentin/enamel adhesives including Dentsply Sirona adhesives designed for use with visible light-cured composite restoratives (see complete Instructions for Use of selected adhesive).

2 Safety notes

Be aware of the following general safety notes and the special safety notes in other chapters of these Instructions for Use.

Safety alert symbol.

- This is the safety alert symbol. It is used to alert you to potential personal injury hazards.
- Obey all safety messages that follow this symbol to avoid possible injury.

Information symbol.

This information symbol is used to highlight important information for the correct use.

2.1 Warnings

The material contains methacrylates and polymerizable monomers which may be irritating to skin, eyes and oral mucosa and may cause sensitization by skin contact and allergic contact dermatitis in susceptible persons.

- **Avoid eye contact** to prevent irritation and possible corneal damage. In case of contact with eyes rinse with plenty of water and seek medical attention.
- **Avoid skin contact** to prevent irritation and possible allergic response. In case of contact, red-dish rashes may be seen on the skin. If contact with skin occurs, remove material with cotton and alcohol and wash thoroughly with soap and water. In case of skin sensitization or rash, discontinue use and seek medical attention.
- **Avoid contact with oral soft tissues/mucosa** to prevent inflammation. If accidental contact occurs, remove material from the tissues. Flush mucosa with plenty of water and expectorate/evacuate the water. If inflammation of mucosa persists, seek medical attention.

2.2 Precautions

This product is intended to be used only as specifically outlined in these Instructions for Use.

Any use of this product inconsistent with these Instructions for Use is at the discretion and sole responsibility of the dental practitioner.

- If refrigerated, allow material to reach room temperature prior to use.
- Contact with saliva, blood or sulcus fluid during application may cause failure of the restoration. Use adequate isolation such as rubber dam.
- Use protective measures for the dental team and patients such as glasses and rubber dam in accordance with local best practice.
- Devices marked "single use" on the labeling are intended for single use only. Discard after use. Do not reuse in other patients in order to prevent cross-contamination.
- The syringes cannot be reprocessed. To prevent syringes from exposure to spatter or spray of body fluids or contaminated hands, it is mandatory that the syringes are handled offside the dental unit with clean/disinfected gloves. Discard syringes if contaminated.
- As additional precautionary measure, syringes may be protected from gross debris but not from all contamination by applying a protective barrier.
- Use of Compules® Tips with the Compules® Tips Gun is recommended.
- The Compules® Tips Gun is only intended to be used with Dentsply Sirona Compules® or Compula® Tips.
- For further information please refer to the Compules® Tips Gun Instructions for Use.
- The shade guide and the individual shade guide tabs cannot be reprocessed. To prevent the shade guide and the individual shade guide tabs from exposure to spatter or spray of body fluids or contaminated hands, it is mandatory that the shade guide and the individual shade guide tabs are handled with clean/disinfected gloves.

- The material should extrude easily. DO NOT USE EXCESSIVE FORCE. Excessive pressure may cause the Compules® Tip to rupture or to eject from the Compules® Tips Gun.
- Tightly close syringes immediately after use.
- Interactions:
 - Do not use eugenol- or hydrogen peroxide-containing materials in conjunction with this product since they may interfere with hardening of the product.
 - Neo Spectra® ST and Neo Spectra® ST Effects restoratives are light-cured materials. Therefore, they should not be exposed to high intensity light sources (e.g., operatory light or LED lights of loupes). Proceed immediately once the material has been placed.
 - If mineral-impregnated (e.g., ferric compounds) retraction cords and/or hemostatic solutions are used in conjunction with adhesive procedures, the marginal seal may be adversely affected, allowing microleakage, subsurface staining and/or restoration failure. If gingival retraction is necessary, use of plain, non-impregnated cord is recommended.

2.3 Adverse reactions

- Eye contact: Irritation and possible corneal damage.
- Skin contact: Irritation or possible allergic response. Reddish rashes may be seen on the skin.
- Contact with mucous membranes: Inflammation (see [2.1 Warnings](#)).

2.4 Storage conditions

- Inadequate storage conditions may shorten the shelf life and may lead to malfunction of the product.
- Store at temperatures between 2 °C and 28 °C (36 °F and 82 °F).
 - Keep out of direct sunlight and protect from moisture.
 - Do not freeze.
 - Do not use after expiration date.

3 Step-by-step instructions

3.1 Direct restoration

Such as: Cavity restoration, veneering, cosmetic and functional reshaping, core build-up, splinting, amalgam, composite, ceramic and porcelain repairs.

3.1.1 Shade selection

1. Prior to shade selection, clean the tooth surface with prophylaxis paste to remove plaque or surface stain.
2. Select the shade while teeth are hydrated.

Neo Spectra® ST restorative has been scientifically designed to allow the tooth structure enamel and dentin to blend with the composite with a life-like translucent result. The 5 A shades of Neo Spectra® ST restorative adequately reproduce the traditional 16 VITA® shades.

For esthetically demanding cases requiring special effects to achieve the desired esthetic result, Neo Spectra® ST Effects offers 2 dentin shades and 1 enamel shade that work with the Neo Spectra® ST universal CLOUD™ shades in a simplified layering system to provide added opacity and/or translucency as needed.

For shade selection, there are two tools that may be used:

- The Neo Spectra® ST and Neo Spectra® ST Effects shade guide. The Neo Spectra® ST and Neo Spectra® ST Effects shade guide consists of tabs, which are manufactured from composite with optical properties close to those of Neo Spectra® ST and Neo Spectra® ST Effects.
- The VITA® classical shade guide in combination with the i-shade label. The Neo Spectra® ST and Neo Spectra® ST Effects i-shade label facilitates shade selection by assigning one of the five universal CLOUD™ shades to each of the 16 VITA® shades.

3.1.1.1 Shade selection with the Neo Spectra® ST and Neo Spectra® ST Effects shade guide

1. Use the shade tabs of the Neo Spectra® ST and Neo Spectra® ST Effects shade guide to select the shade closest to the shade of the area to be restored.

Avoid contamination of the shade guide.

To prevent the shade guide from exposure to spatter or spray of body fluids or contaminated hands, it is mandatory that the shade guide is handled offside the dental unit with clean/disinfected gloves. For shade selection, remove individual tabs from shade guide holder.

3.1.1.2 Shade selection with VITA® classical shade guide and i-shade label

Have the i-shade label ready.

To make sure that the i-shade label is available chairside, we recommend sticking it on the rear side of the VITA® classical shade guide.

1. Select the VITA® classical reference tooth whose central part best matches the area of the tooth to be restored.
2. Use the i-shade label (see table below) to identify the matching Neo Spectra® ST universal CLOUD™ shade, and appropriate Neo Spectra® ST Effects shade(s) if enhanced esthetic effects are desired.

Tooth shade according to the VITA® Classical Shade Guide	A1	B1	C1	A2	B2	D2	A3	C2	D3	D4	A3,5	B3	B4	C3	A4	C4
Recommended Neo Spectra® ST universal CLOUD™ shade	A1			A2			A3			A3,5			A4			
Recommended Neo Spectra® ST Effects shades	D1						D3									
	E1															

3.1.2 Cavity preparation

1. Clean the tooth surface with a fluoride-free prophylaxis paste (e.g., Nupro® prophylaxis paste).
2. Prepare cavity (unless not necessary, e.g., cervical lesion).
3. Use adequate isolation such as rubber dam.
4. Rinse surface with water spray and carefully remove rinsing water. Do not desiccate the tooth structure.

3.1.3 Matrix placement and wedging (Class II restorations)

1. Place a matrix (e.g., AutoMatrix® or Palodent® matrix system) and wedge. Burnishing of the matrix band will improve contact and contour. Pre-wedging or Palodent® V3 ring placement is recommended.

3.1.4 Pulp protection, tooth conditioning/dentin pre-treatment, adhesive application

Refer to adhesive manufacturer's Instructions for Use for pulp protection, tooth conditioning and/or adhesive application. Once the surfaces have been properly treated, they must be kept uncontaminated. Proceed immediately to placement of the material.

3.1.5 Application

Apply material in increments immediately after the application of the adhesive. The material is designed to be cured in increments up to a 2 mm depth/thickness.

Using Compules® Tips

1. Insert Compules® Tip into the notched opening of the Compules® Tips Gun. Be certain that the collar on the Compules® Tip is inserted first.
2. Remove the colored cap from the Compules® Tip. The Compules® Tip may be rotated 360° to gain the proper angle of entrance into the cavity or to the mixing pad.
DO NOT USE EXCESSIVE FORCE.
Indirect use: Dispense the necessary amount of material onto a mixing pad and protect against light. Apply material into the cavity with a plastic instrument.
- or -
Direct intraoral use: Dispense the material directly into the cavity preparation using a slow, steady pressure.
3. To remove the used Compules® Tip, be sure that the Compules® Tips Gun plunger is pulled back completely by allowing the handle to open to its widest position. Apply a downward motion to the front end of the Compules® Tip and remove.

Excessive force.

Injury.

1. Apply slow and steady pressure on the Compules® Tips Gun.
2. Do not use excessive force. Compules® Tip rupture or ejection from Compules® Tips Gun may result.

Using syringes

Avoid contamination of the syringe.

To prevent the syringe from exposure to spatter or spray of body fluids or contaminated hands, it is mandatory that the syringe is handled offside the dental unit with clean/disinfected gloves.

1. Remove the cap.
2. Turn the handle of the syringe slowly in a clockwise direction and dispense the necessary amount of the material onto a mixing pad.
3. Point the front tip of the syringe upwards and turn the handle anti-clockwise to prevent oozing of the material.
4. Immediately close the syringe with the cap.
5. Protect the restorative material on the mixing pad against light.

Optional flowable liner: Placement of a compatible flowable liner (available separately) prior to placement of Neo Spectra® ST and Neo Spectra® ST Effects restoratives is optional. Follow manufacturer's Instructions for Use.

3.1.5.1 Applying restorative

1. Dispense the composite material directly into the cavity preparation/tooth surface from the Compules® Tip using slow, steady pressure. Excessive force is not necessary. Alternatively, material may be expressed onto a clean pad from the Compules® Tip or syringe and carried to the preparation with suitable placement instrument.
2. Adapt, contour and shape with appropriate composite instruments. Material may be placed and light cured in increments up to 2 mm (see [3.1.6 Light-curing](#)).

Technique Tips:

- Use a single opacity for most common cases: Because the Neo Spectra® ST universal CLOUD™ shades are designed to replace both dentin and enamel, a single universal CLOUD™ shade is sufficient for most restorations. If esthetic blending of more than one universal CLOUD™ shade is desired, shade layers may be individually cured to form a foundation for subsequent shade placement.
- Use multiple opacities for esthetically demanding cases: For esthetically demanding cases, such as a large Class IV or a through-and-through Class III, Neo Spectra® ST Effects dentin and enamel shades can be used with Neo Spectra® ST universal CLOUD™ shades in a simplified layering system to provide added opacity and/or translucency.
 - Create a thin palatal (lingual) shell (approx. 1 mm thick) using the selected universal shade to provide a guide for achieving accurate restoration shape.
 - Use dentin shade to recreate dentin core, not fully extending to the proximal and incisal edges. Leave at least 1 mm of space on the vestibular (facial) surface for the universal shade.
 - Complete restoration using the universal shade for the vestibular (facial) surface. For a highly translucent incisal edge, use the universal shade to restore the vestibular (facial) area, leaving space on the incisal third for the enamel shade.

3.1.6 Light-curing

1. Light-cure each increment with a suitable curing light designed to cure materials containing camphorquinone (CQ) initiator, i.e., peak of spectrum in the range of 440-480 nm and irradiance of 800-2000 mW/cm². Cure for at least 20 seconds if light output is between 550 and 800 mW/cm².
2. Cure each increment for the time indicated in the table below, considering increment thickness, shade and curing light output.

Curing time for 2 mm increments

Neo Spectra® ST and Neo Spectra® ST Effects Curing Recommendations		
Shades	Output mW/cm ²	Curing Time
A1-A4, BW	≥ 550	20 seconds
	≥ 800	10 seconds
D1, D3	≥ 550	40 seconds
	≥ 800	30 seconds
E1	≥ 550	10 seconds

Curing timetable.

The curing timetable above also appears on all outer packages.

Insufficient curing.

Inadequate polymerization.

1. Check compatibility of curing light.
2. Check curing cycle.
3. Check minimum irradiance.
4. Cure each area of each increment for the recommended curing time.
5. Check distance to surface to be cured.

3.1.7 Finishing and polishing

1. Contour the restoration using finishing burs or diamonds.
2. Use Enhance® finishing devices for additional finishing.
3. Prior to polishing, check that surface is smooth and free of defects.
4. For achieving a very high luster of the restoration, we recommend Enhance® and Enhance® PoGo®, finishing and polishing instruments, Prisma® Gloss™ Regular and Extrafine polishing paste.

For contouring, finishing, and/or polishing, follow the manufacturer's Instructions for Use.

3.2 Indirect restoration: Inlay, onlay or veneer

3.2.1 Shade selection

Final desired shade selection should be accomplished prior to tooth preparation. Refer to [3.1.1 Shade Selection](#).

3.2.2 Preparation

Design requirements are essentially a conventional preparation. Rounding of internal angles and refinement of the cavo-surface margin for enhancement of enamel bonding and finishing procedures is recommended. Refer to adhesive and/or luting cement manufacturer's instructions for base/liner/pulp protection requirements.

3.2.3 Impression, master cast fabrication

1. Make an accurate impression of the preparation.
2. Prepare master cast and articulate opposing cast per usual technique.

3.2.4 Restoration fabrication

1. Block out any preparation undercuts present. Apply separating medium and/or die spacer to separated die. In most cases, the restoration may be fabricated in no more than 3 increments, each up to 2 mm depth.
2. Apply first increment/shade, creating restoration body. Adapt to die, staying short of margins. Recommended light-curing is accomplished by placing in Triad® Light Curing Unit (Dentsply Sirona) for 2 minutes. Apply second, body layer, allowing cut-back for occlusal anatomy and final proximal contour. Repeat Triad® light-curing. Prior to final "enamel" shade placement, custom staining with appropriate compatible stains may be accomplished per manufacturer's instructions. Apply final layer, slightly overfilling and covering all margins. Replace die into articulated model. Establish all external, proximal and occlusal contacts and anatomy. Slight lubrication of adjacent and opposing model teeth is recommended. Cure briefly (10 seconds) with hand-held visible curing light unit to fix contours. Remove die with restoration, place in Triad® Light Curing Unit for final 2 minute curing.
3. When removing restoration from die, it may be necessary to scrape die stone away from restoration margins, to prevent accidental chipping of restoration. Clean any residual die stone from restoration. Carefully trim away visible flash beyond preparation margins with acrylic bur.
4. Gently sandblast the internal surfaces with 50 µm alumina abrasive. Remove visibly undercut material.
5. Seat restoration onto preparation of uncut master cast, adjusting as needed. Check for marginal integrity and overall fit and contour. Additional increments may be added if needed by roughening surface and applying compatible adhesive per manufacturer's instructions, followed by placement and light curing as outlined above.

3.2.5 Finishing and polishing (laboratory)

Complete as outlined in [3.1.7 Finishing and polishing](#).

3.2.6 Cementation

Recommended cementation technique is bonding using a resin cement. Follow adhesive and/or cement manufacturer's instructions for prepared tooth and restoration surface pretreatments.

3.2.7 Adjusting, finishing and polishing (clinically)

Following cementation, make all necessary occlusal adjustments and polish any clinically adjusted surfaces as outlined in [3.1.7 Finishing and polishing](#).

4 Hygiene

Cross-contamination.

Infection.

- Do not reuse single-use products. Dispose of in accordance with local regulations.
- Syringes cannot be reprocessed. Dispose of contaminated syringes in accordance with local regulations.
- Reprocess reusable products according to the instructions.

4.1 Instructions for cleaning and disinfecting indirectly fabricated restorations

Process		✓ Acceptable
Cleaning	Automated (Washer/Disinfector)	
	Mechanical	✓
	Manual	✓
Disinfection	Automated (Washer/Disinfector)	
	Mechanical	✓
	Manual	✓
Sterilization	Steam Autoclave	
	Immersion	

Instructions for cleaning and disinfecting indirectly fabricated restorations

Warnings	• These instructions are for use ONLY for an indirectly fabricated restoration, prior to intraoral delivery. • High-level disinfection procedures have not been validated. Intermediate-level disinfection is appropriate for indirect restorations. • Do not autoclave in steam autoclave. • Do not clean/disinfect in automated washer/disinfector. • Follow proper infection prevention activities, like proper hand washing and donning new puncture and chemical resistant gloves at pertinent steps. • Do not use hard wire brushes to clean the restoration, as damage can occur. • Use only a disinfecting solution that has been approved for its efficacy, is EPA registered (and/or Health Canada approved) and use in accordance with the Instructions for Use of the disinfecting solution's manufacturer.
-----------------	---

Limitations on reprocessing	- Repeated reprocessing has minimum effect on restorations. End of life is normally determined by clinical factors. Restorations are typically processed once prior to patient delivery. - Steam-autoclaving, automated washer/disinfection, chemical-vapor-sterilization, and dry-heat-sterilization methods have not been tested or validated for efficacy and are not recommended for use.
Initial treatment at the point of use	- Do not allow residue or any form of contamination to dry on the device. - Use a new, clean pair of examination gloves. - Remove excess soil/gross debris by wiping with a disposable cloth/cotton-tipped applicator saturated with 70% isopropyl alcohol (IPA). - Remove all visible soil, ensuring fluid penetrates all crevices. Use fresh wipes to rub fluid into the crevices. Discard used wipes. Additional wipes may be used. - Transport to processing area and follow validated processing instructions below.
Preparation before cleaning	Polish external surfaces and prepare internal surfaces for cementation prior to processing.
Cleaning and Disinfection: Automated	Do not use automated washer/disinfectors for processing indirectly fabricated composite restorations. Material may become distorted or compromised.
Cleaning and Disinfection: Mechanical	- The restoration may be concurrently cleaned and disinfected if using the recommended solution 70% isopropyl alcohol (IPA). - Using a cloth and a cotton-tipped applicator saturated with 70% IPA, wipe all the surfaces of the restoration, ensuring to wipe all grooves and edges. - Immerse the restoration in an ultrasonic bath in 70% IPA solution for at least 1 minute, unless otherwise specified by the solution's manufacturer. - Rinse under tap water of at least drinking quality for a minimum of 30 seconds. - Allow to air-dry until completely dry. - Visually inspect for visible soil and repeat the cleaning and disinfection procedure, if necessary.
Cleaning: Manual	- Using a cloth and a cotton-tipped applicator saturated with 70% IPA, wipe all the surfaces of the restoration, ensuring to wipe all grooves and edges. - Immerse the restoration in 70% IPA solution for at least 1 minute, unless otherwise specified by the solution's manufacturer. - Using a cotton-tipped applicator and/or soft bristled brush saturated with 70% IPA, wipe all the surfaces of the restoration, ensuring to wipe all grooves and edges. Continue to scrub until all gross soil is removed. - Rinse under tap water of at least drinking quality for a minimum of 30 seconds. - Allow to air-dry until completely dry. - Visually inspect for visible soil and repeat the cleaning and disinfection procedure, if necessary.

Disinfection: Manual	- Immerse the restoration in 70% IPA solution for at least 1 minute, unless otherwise specified by the solution's manufacturer. - Rinse under tap water of at least drinking quality for a minimum of 30 seconds. - Allow to air-dry until completely dry. - Visually inspect for visible soil and repeat the cleaning and disinfection procedure, if necessary.
Drying	Wipe the devices dry with a clean, lint-free cloth. Allow the components to fully air-dry before delivery.
Maintenance, Inspection and Testing	Immediately before delivery, visually inspect to ensure that all contamination has been removed. If the device is discolored, damaged, or distorted it should be discarded.
Packaging	No particular requirements.
Sterilization	Sterilization is not allowed. No methods have been validated. Do not subject restorations to steam-autoclaving or liquid-chemical-sterilant-immersion. Restoration damage may occur.
Storage	- Store processed restorations at room temperature, away from moisture or excessive humidity. - To prevent contamination, store processed restorations in covered storage, like a drawer or cabinet, until use.

4.2 Compules® Tips Gun

For reprocessing instructions please refer to the Instructions for Use of the Compules® Tips Gun, which is available on our webpage at <https://www.dentsplysirona.com/ifu>. If requested, we will send you a free printed copy of the Instructions for Use in the language you require within 7 days. You can use the order form provided on <https://www.dentsplysirona.com/ifu> for this purpose.

4.3 Compules® Tips – cross-contamination

Cross-contamination.

Infection.

- Do not reuse Compules® Tips. Dispose of Compules® Tips in accordance with local regulations.

4.4 Syringes – cross-contamination

Cross-contamination.

Infection.

- The syringes cannot be reprocessed.
- To prevent the syringes from exposure to spatter or spray of body fluids or contaminated hands, it is mandatory that the syringes are handled offside the dental unit with clean/disinfected gloves. Do not reuse syringes if contaminated.
- Dispose of contaminated syringes in accordance with local regulations.

To prevent syringes from exposure to spatter or spray of body fluids or contaminated hands, or oral tissues, use of a protective barrier is recommended. The use of protective barriers is an additional precautionary measure against gross debris but not against all contamination.

Incidental contact of the syringe with water, soap or a water-based hospital-level disinfection solution will not damage syringe body. Do not allow any solution contact with contained material. Discard composite material that has been in contact with any fluid or non-sterile instrument. Dry the syringe with a lint-free single-use cloth. Wipe syringe gently.

4.5 Shade guide and individual shade guide tabs – cross-contamination

Cross-contamination.

Infection.

- The shade guide and the individual shade guide tabs cannot be reprocessed.
- To prevent the shade guide and the individual shade guide tabs from exposure to spatter or spray of body fluids or contaminated hands, it is mandatory that the shade guide and the individual shade guide tabs are handled with clean/disinfected gloves. Do not reuse shade guide and the individual shade guide tabs if contaminated.
- Dispose of contaminated shade guide and the individual shade guide tabs in accordance with local regulations.

4.6 Disposal

Dispose of in accordance with local regulations.

5 Lot number (), expiration date () and correspondence

1. Do not use after expiration date.
ISO standard is used: "YYYY-MM" or "YYYY-MM-DD"
2. The following numbers should be quoted in all correspondence:
 - Reference number (REF)
 - Lot number
 - Expiration date
3. Any serious incident in relation to the product should be reported to the manufacturer and the competent authority according to local regulations.
4. Device Identification (Basic UDI-DI): ++D010FIM03R3

A summary of the safety and clinical performance (SSCP) for this product can be found (upon activation) at <https://ec.europa.eu/tools/eudamed> by searching using the Basic UDI-DI number listed above and at <https://www.dentsplysirona.com/ifu> using the reference number (REF).

© Dentsply Sirona 2023-06-29

[These Instructions for Use are based on Master Version 02]

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Dentsply DeTrey GmbH