Zhermack
Dental

Zeta 3 Soft

EN - Instructions for use 2
FR - Instructions d'utilisation 5
IT - Istruzioni per l'uso 9
DE - Gebrauchsanweisung 13
ES - Instrucciones de uso 17
PT - Instruções de utilização 21
NL - Gebruiksinstructies 25
SV - Bruksanvisning 29
DA - Brugsanvisning 32
NO - Brukerinstruksjoner 35
FI - Käyttöohjeet 38
LT - Naudojimo instrukcija 41
LV - Lietošanas norādījumi 44
ET - Kasutusjuhised 47
CS - Návod k použití 50
SK - Návod na použitie 53
HU - Használati utasítás 56
RO - Instrucţiuni de utilizare 60
HR - Upute za uporabu 63
BG - Инструкции за употреба 66
EL - Οδηγίες χρήσης 70
TR - Kullanım Kılavuzu 74

Zhermack
Dental

Zeta 3 Soft

Disinfectant and cleaner for non-invasive medical devices surfaces

750 ml

MD

CE 0123

dentsplysirona.com/zhermack

Zhermack S.p.A. | Via Bovazecchino, 100 | 45021 Badia Polesine (RO) Italy - T +39 0425 597611 | F +39 0425 53596 | info@zhermack.com | www.zhermack.com

2100-EIFU-011-001 Zeta3_Soft_750ml

Last update: 2026-05

EN - Zeta 3 Soft 750 ml

Instructions for use

1. INTENDED USE

Disinfectant and cleaner for non-invasive medical devices surfaces.

2. PRODUCT DESCRIPTION

Zeta 3 Soft is a ready-to-use alcohol-based disinfectant for the cleaning and disinfection of the surfaces of non-invasive medical devices.

3. PACKAGING

750 ml bottle.

4. COMPOSITION

Zeta 3 Soft Lemon fragrance: 100 g of Zeta 3 Soft contain 34.4 g of ethanol, 14 g of isopropanol, nonionic surfactants, additives, auxiliaries, fragrance (limonene) and water q.s. for 100 g.

Zeta 3 Soft Classic fragrance: 100 g of Zeta 3 Soft contain 34.4 g of ethanol, 14 g of isopropanol, nonionic surfactants, additives, auxiliaries, fragrance (limonene, citral, Citrus limon peel oil) and water q.s. for 100 g.

5. INDICATIONS FOR USE

The device is intended for use in the dental sector, in professional dental surgeries (public or private health facilities with the necessary legal authorisations), by professionals trained and qualified for the disinfection and cleaning of the surfaces of non-invasive medical devices.

6. CONTRAINDICATIONS

Use of the product is contraindicated for operators with known hypersensitivity to limonene, citral and Citrus limon peel oil.

7. SIDE EFFECTS

Allergic reactions may occur following contact with Zeta 3 Soft. Direct contact with eyes causes severe eye irritation, and with skin can cause skin irritation. Inhaling its vapours can cause respiratory tract irritation.

8. CLINICAL BENEFITS

Zeta 3 Soft helps to decrease the microbial load of the surfaces of non-invasive medical devices in order to reduce the risk of transmitting infection if used in strict compliance with the instructions for use and the product safety data sheet.

9. WARNINGS

- This product contains limonene, citral and Citrus limon peel oil, which may cause an allergic skin reaction. Use by operators sensitive to limonene, citral and Citrus limon peel oil is contraindicated. If irritation, redness or other signs of hypersensitivity occur, stop using the product and consult a doctor.
- Causes serious eye irritation. Avoid direct contact with eyes. Protect eyes with airtight goggles. In case of contact with eyes, rinse immediately with plenty of water for several minutes and seek medical advice if irritation persists.
- May cause respiratory irritation. Do not inhale vapours. Wear a face mask (see section 8 of the product safety data sheet), adequately ventilate places in which the product is handled and ensure that the bottle is stored tightly closed.
- May cause skin irritation. Avoid direct contact with skin and face. Protect hands with gloves, wear a face mask, protective clothing and safety footwear for professional use (see section 8 of the product safety data sheet). In the event of accidental contact with the skin, wash thoroughly with soap and water and seek medical advice if necessary.
- Flammable liquid and vapour. Keep away from heat, sparks, naked flames and hot surfaces and other sources of ignition. Ensure the bottle is stored tightly closed and in a well-ventilated place.
- Only use the product in the surgical area when there are no patients and/or individuals not wearing appropriate PPE.
- To reduce the risk of cross-contamination, always use new, uncontaminated gloves when handling the products (bottle and spray dispenser).

10. PRECAUTIONS

- Keep the product's instructions for use for the entire duration of its use.
- Zeta 3 Soft is compatible with hard surfaces, in particular those made of: Brass; Stainless steel; Polycarbonate; Polyolefins, such as (HD)PE; Nylon; Polystyrene; POM; PMMA; ABS. Do not allow the product to come into contact with aluminium surfaces and surfaces not compatible with alcohol. For further limitations of use, read the instructions for use of the devices to be disinfected.
- Avoid contact of Zeta 3 Soft liquid with other disinfectants or cleaners.
- Using the wipe impregnated with the product, wet the entire area to be disinfected with the product and ensure that the area remains wet for the duration of the action, even if the product is reapplied several times. Failure to follow these steps may compromise the disinfection process.
- Zeta 3 Soft must be used at room temperature (23°C/73°F).
- Always comply with the contact times indicated in Table 1. Exposure for longer times may result in a deterioration of the treated surfaces in the long term. Prolonged contact with plastic surfaces containing high concentrations of plasticisers (such as PVC) may alter their properties and reduce their durability.
- Do not use the product and dispose of it in accordance with the instructions and the product safety data sheet:
 - if the batch number and expiry date on the outer packaging are missing or illegible;

- if the primary packaging (bottle, spray dispenser) is damaged or does not work properly.

• Do not use after the expiry date indicated on the package and indicated on the front of the label in year-month-day format.

• Once the contents inside the bottle have been used up, dispose of the packaging according to section 14.

• Do not use the bottle for purposes other than those indicated in these instructions.

11. STEP-BY-STEP INSTRUCTIONS

11.1 PRELIMINARY OPERATIONS

1. Read the instructions for use and the product safety data sheet before using.

2. Adequately ventilate the places where the product is used.

3. Protect hands with a pair of gloves, protect eyes and face with a pair of airtight goggles and wear a face mask, work clothes and safety footwear for professional use (see section 8 of the product safety data sheet).

4. In order to verify compatibility, it is advisable always to test the product on a non-visible part of the surface before proceeding with disinfection.

11.2 USING THE DEVICE

1. Spray Zeta 3 Soft onto a disposable wipe until it is completely wet* and wipe the surfaces and areas to be cleaned until all visible residues have been removed. Dispose of the wipe as described in section 14.

2. Repeat by spraying Zeta 3 Soft onto a new disposable wipe until it is completely wet* and wipe the previously cleaned surfaces and areas to be disinfected until they are completely wet. Dispose of the wipe as described in section 14.

3. Wait for the 5-minute action time.

*A 20x25 cm wipe onto which sufficient Zeta 3 Soft has been applied to completely wet the fabric (for low absorbency fabrics such as polypropylene or PET apply 3 to 7 sprays, for high absorbency fabrics such as viscose apply 5 to 14 sprays) covers an area $\leq 0.25\ m^2$.

12. INFORMATION ON THE PROPER MANAGEMENT OF PACKAGING

Product packaging can be reused safely without being cleaned and disinfected, provided the following guidelines for packaging conservation and storage are observed:

• Only use the product when there are no patients in the surgical area;

• Always use new, uncontaminated gloves when handling packaging. If gloves become contaminated, change them immediately before handling the product;

• Do not bring packaging into contact with contaminated surfaces and/or devices;

• Store packaging tightly sealed and outside the surgical area in a clean, enclosed and dry place away from splashes, sprays or aerosols of body fluids.

In case of biological contamination, dispose of packaging and its content as indicated in paragraph 14 dedicated to disposal.

13. SPECTRUM OF ACTION

Table 1

Effectiveness (micro-organisms)	**Reference standard** *Tests carried out in dirty conditions*	**Contact time** *(minutes)*
Bactericidal (S. aureus, P. aeruginosa, E. hirae)	EN 13727, EN 16615	*5**
Yeasticidal (C. albicans)	EN 13624, EN 16615	5
Tuberculocidal (M. terrae)	EN 14348	5
Limited Virucidal (norovirus, rotavirus, adenovirus and all enveloped viruses, including HIV, HBV and HCV)	EN 14476	3

*tests carried out in 5 minutes in dirty conditions except EN 16615 passed in 5 minutes in clean conditions and in 10 minutes in dirty conditions.

Activity tested in compliance with EN 13727, EN 16615, EN 13624, EN 14348, EN 14476.

14. STORAGE, STABILITY AND DISPOSAL

Store the product tightly closed outside the operating area in a clean, closed, cool and dry environment protected from splashes, sprays or aerosols, at a maximum temperature of 27°C/80°F, away from heat sources, sparks, open flames and hot surfaces and other sources of ignition and direct sunlight. Do not use after the use-by date indicated on the packaging.

The following disposal instructions must be considered valid only if the product has been stored and used according to these instructions for use:

• Waste with only chemical residues (residues of disinfectant solution) are a physical hazard for people and the environment and must be handled as **hazardous waste** in accordance with local regulations.

• Waste with biological residues must be disposed of as **special waste** at risk of biological contamination.

15. IMPORTANT REMARKS

Information provided in any way, even during demos, does not invalidate the instructions for use. Operators are required to check that the product is suitable for the application envisaged. The manufacturer cannot be held liable for damage, including to third parties, deriving from failure to follow instructions or from unsuitability for an application. The manufacturer's liability is, in any case, limited to the value of the products supplied. Report any serious incident involving the medical device to the manufacturer and to the relevant authorities.

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