CHARISMA® E4SY Flow

CHARISMA® E4SY Flow

GB Instructions for use
ES Instrucciones de uso
GR Οδηγίες Χρήσης
CZ Návod k použití
HU Használati utasítás
LV Lietošanas instrukcija
LT Naudojimo instrukcija
PL Instrukcja obsługi
BR Instruções de uso
HR Upute za uporabu
RO Instrucțiuni de utilizare
TR Kullanma talimatı
RU Инструкция по применению
AE تعليمات الاستخدام
CN 使用说明
KR 사용 설명서
SK Návod na použitie
BG Инструкции за употреба
UA Інструкція із застосування
EE Kasutusjuhend
SI Navodila za uporabo

Manufacturer:
Kulzer GmbH
Leipziger Straße 2
63450 Hanau (Germany)

Made in Germany

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Dentistry – Polymer-based restorative materials
ISO 4049 (type 1, class 2, group 1)

CE 0197

22 PAP

Explanation of symbols on labelling

- MD Medical device
- Manufacturer
- LOT Batch code
- Use-by date
- Date of manufacture
- Consult instructions for use
- Storage temperature
- Keep dry
- Keep away from sunlight / radioactive sources
- Contains hazardous substances
- Contains nano materials
- UDI Unique Device Identification

VITA® is a registered trademark of VITA Zahnfabrik H. Rauter GmbH & Co. KG

66101326/01

CHARISMA® E4SY Flow — Instructions for use (GB)

Medical Device – for use by dental health care professionals only.

Intended purpose
CHARISMA E4SY Flow is a flowable, light-curing, radiopaque nano-hybrid composite. It is indicated for adhesively bonded, tooth colored restorations of anterior and posterior teeth including occlusal surfaces. Its radiopacity is equivalent to 3.0 – 3.5 mm of aluminium (dentine = 1 mm, enamel = 2 mm).

Indications
• Enlarged fissure sealing • Cavity lining – as the first layer for Class I and II cavities • Class-I, -II, -III, -IV and -V direct restorations • Small repairs of direct and indirect restorations combined with a suitable bonding agent • Splinting of teeth • Composite attachments for aligner treatment • Fixation of retainers

Contraindications
Application to dentine close to the pulp (less than 1 mm). In case of known or suspected allergy to components of the product, its use is contraindicated. Do not use this product in case of known or suspected allergies against (meth)acrylate compounds.

Patient target group
This dental material is suitable for patients requiring dental treatment for the mentioned indications with consideration of the contraindications. For pregnant and nursing women, the treatment risks must be weighed carefully against the benefits taking into consideration the unborn child or infant.

Intended user
For use by dental health care professionals only.

Clinical benefit
Kulzer medical devices ensure the rehabilitation of oral functions such as chewing, speaking and aesthetics.

Potential side effects
Hypersensitivities to the product or its components cannot be excluded in individual cases.

Composition
Barium Aluminium Boro Fluor Silicate Glass, E4BADMA, UDMA, Ytterbiumfluorid, Silica, Titanium Dioxid, BHT, Aminobenzoicacidester, Camphorquinone, Fluorescent Pigments, Metallic Oxide Pigments, Organic Pigments.
CHARISMA E4SY Flow contains approximately 41% filler by volume, its filler particle size is 20 nm - 5 µm.

Shades

Shades for monochromatic layering

Corresponding VITA® shades

LIGHT	A1, B1, B2, C1
MEDIUM	A2, A3, B3, B4, C2, D2, D3
DARK	A3.5, A4, C3, C4, D4
BLEACH	BXL, BL

Product application
Shade selection should be done using the CHARISMA E4SY Flow shade guide before inserting the dental dam as the teeth change their colour noticeably once dry.
1. For ideal results, use Ivory rubber dam (or similar) to control contamination from moisture, blood and saliva.
2. Remove any caries and protect any hard dental structures according to the general guidelines for preparation.
3. A quality adhesive system (e.g., GLUMA® Bondings by Kulzer) should be applied to the enamel and dentine according to the manufacturer's instructions before inserting the filling material. Adhesive-treated cavities must not be contaminated with moisture, blood or saliva.
4. Apply CHARISMA E4SY Flow in thin layers (max. 2 mm, baseliner max. 1 mm) to the cavity. When applying the CHARISMA E4SY Flow, ensure that it flows slowly and evenly to prevent the formation of air bubbles. If necessary, remove any air bubbles using a probe before polymerisation.
5. Polymerise each layer according to the table below using a suitable dental light-curing unit (max. wavelength between 450 – 480 nm; light power: 1550 – 600 mW/cm²; e.g. Kulzer Translux®). If the tooth surface being treated is not to be instrumented by a diamond stone, carbide bur, drill or other abrasive instrument, the surface should be cleaned with pumice or a fluoride free polishing paste before restoration is placed. A non-sticking placement instrument is recommended for use in clinically placing, sculpting and shaping the material.

Recommended CHARISMA E4SY Flow-curing times and maximum layer thicknesses for Translux®-light-curing units or other suitable dental light-curing units:

Shade	Time	Curing depth
LIGHT, MEDIUM, DARK, BLEACH	20 s	2 mm

6. During polymerisation a dispersion layer forms on the surface. This layer should not be touched or removed if further composite layers are to be applied
7. CHARISMA E4SY Flow can be prepared and polished immediately after polymerisation. Finishing diamonds, flexible discs, silicone polishers (e.g., CHARISMA EasyShine) and polishing brushes are suitable for preparation.

Injection technique
When using the injection technique, prepare a transparent silicone key (e.g. Memosil® 2) on the intact row of teeth or a mock-up enclosing the tooth/teeth to be treated before injecting the flowable composite. Remove undercuts. The silicone key must fit perfectly and without gaps. Apply flowable composite slowly and with even, gentle pressure via previously inserted injection channels. The cannula should be pulled out of the canal slowly and carefully while simultaneously and continuously dispensing the composite. In the case of dissolved proximal contacts, ensure adequate interproximal separation (e.g. with Teflon). Expose the composite according to the table using a suitable dental light-curing unit (wavelength 450 - 480 nm, light output: 1550 - 600 mW/cm², e.g. Kulzer Translux®), from each surface. After removing the silicone key, polymerise the filling once again, finish and polish the composite surface as usual.

Application of attachments for aligner treatment
1. Fill a composite into the attachment's recesses of a compatible template. Prepare the enamel to be treated as described under Product application. Place the template on the prepared row of teeth and cure the attachments according to the curing times table. Remove any excess composite if necessary.
2. After completion of the aligner treatment, remove the composite attachments with rotating instruments and polish the tooth surface, e.g. with CHARISMA EasyShine.
3. UV-A -light may be used for detection in order to differentiate between composite and enamel.

Instructions
During long restorative procedures the surgery light should be temporarily directed away from the treatment area to prevent premature initiation of the polymerisation reaction. The polymerisation times specified above are minimum times and should not be reduced so that an optimal polymerisation rate, and the associated physical material properties, can be reached. Afterward, apply and cure the flow composite. After ejection of the material from the syringe, release the plunger to prevent continued flow. At the end of treatment, remove the hygienic protective cover from the syringe and discard it. Discard the used applica-

Storage conditions
Store at temperatures of -10 to 25°C (14 to 77°F). Do not expose to direct sunlight. Protect from ionising radiation. Store in a dry place. Do not use the material after the expiration date.

Warnings and safety instructions
May cause skin and eye irritation. Avoid contact with skin. May cause an allergic skin reaction. Wear suitable protective gloves. In case of contact with skin, immediately wash with plenty of water and soap. Avoid contact with eyes. In case of contact with eyes, rinse with water for several minutes. If symptoms persist, consult an eye specialist. May damage fertility or the unborn child.

General information
Please pay attention to the following safety instructions and all processing instructions in other sections. This product must be used according to the instruction for use only. Any use not complying with these instructions is at the discretion and sole responsibility of the dental professional itself.

Application notes
Visually inspect for damage before use. Damaged products must not be used. Composites should only be used at room temperature. Areas near the pulp should be covered with a suitable pulp protection material according to the manufacturer's instructions. Do not use eugenol containing materials in or around the cavity preparation being restored, because these may interfere with the proper curing of the composite. When using syringes in mouth, use common hygienic protective sleeves to prevent contamination. Do not reuse after contamination. Note: Light-curing units can have an immense variation in intensity and emitted spectral output. User should check with the light-curing unit manufacturer and verify its efficacy non-clinically before using in a patient.

Dated: 2025-12

Disposal information
Recommendation: Dispose of in accordance with official regulations. Do not dispose of contents or partially emptied packaging in the household waste or allow them to enter the sewage system. European Waste Catalogue: 180106 Chemicals consisting of or containing hazardous substances.

Complaint notice
For SSCP please see EUDAMED database (https://ec.europa.eu/tools/eudamed) or contact sscp-service@kulzer-dental.com
To find the SSCP in the EUDAMED use the UDI code attached at the packaging of the product. Please quote batch number and article number or UDI in all correspondence about the product. Article number REF and batch number LOT or UDI: See note on syringe/packaging.
According to EU Medical Device Regulation users/ patients are obliged to report serious incidents with a medical device to the manufacturer and to the competent authority of the country, the incident occurred.

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