



VALO™

BROADBAND **LED** CURING LIGHT

10+ YEARS OF INNOVATION.





BROADBAND LED CURING LIGHT

Suggested Curing Times



	DENTAL	STANDARD	HIGH POWER	XTRA POWER
POWER BUTTON				
TIMING INTERVAL LIGHTS				
MODE/STATUS LIGHT				
TIME/MODE CHANGE BUTTON				
		5 sec. 10 sec. 15 sec. 20 sec.	1 sec. 2 sec. 3 sec. 4 sec.	3 sec.
PER 2 mm LAYER		1 X 10 SECONDS	2 X 4 SECONDS	1 X 3 SECONDS
FINAL CURE		1 X 20 SECONDS	3 X 4 SECONDS	2 X 3 SECONDS
ORTHO				
METAL & CERAMIC BRACKETS		1 X 10 SECONDS	2 X 4 SECONDS	2 X 3 SECONDS

Applicable for all current and future Ultradent Products light-cured materials. Please refer to the light cured materials' Manufacturer Instructions for Use before curing with the VALO light.

ULTRADENT.COM/EU

© 2020 Ultradent Products, Inc. All Rights Reserved.

