




SOOOO STABLE
IS **VISALYS[®]** CORE

CORE BUILD-UP & POST CEMENTATION



KETTENBACHDENTAL

Simply intelligent

VISALYS® CORE – EASY TREATMENT

Dual-curing core build-up & root post cementation material

Visalys® Core – especially for difficult and freehand core build-ups.



EXCELLENT STABILITY

- with very good flow properties at the same time



UNIQUE ACTIVE-CONNECT- TECHNOLOGY

- for optimized adhesive bonding
- for use with common adhesives

DENTIN-LIKE GRINDABILITY




FREE OF BISPHENOL A



Distinguished by the
experts "The Dental
Advisor"

VISALYS® CORE – STABLE RESULTS


Adhesive bonding between Visalys® Core and dentin with different adhesives



In-house measurements by Kettenbach, shear adhesive bond, test based on DIN EN ISO 29022.

Visalys® Core achieves the clinically relevant adhesive bond of approx. 10 MPa.

Average removal depth when grinding various core build-up materials



Tangku Ibrahim Mautana. Grindability analysis of dental core build-up materials under clinical relevant test procedure. Journal of Medical Materials and Technologies (Thesis Series Edition), Vol 1 No 1, 1-21 (2017)


Visalys® Core can be precisely processed similar to dentin.

VISALYS® CORE – THE RANGE

Visalys® Core		
Intro pack	White Dentin	5 mL automix syringe, 10 brown mixing tips, 5 intraoral tips, 5 endo tips
Normal pack	White Dentin	2 x 5 mL automix syringe, 20 brown mixing tips, 10 intraoral tips, 10 endo tips
Normal pack	White Dentin	25 mL cartridge, 20 yellow mixing tips, 20 intraoral tips



Kettenbach GmbH & Co. KG
Im Heerfeld 7
35713 Eschenburg, Germany
Phone +49 (0)2774 7050
info@kettenbach.com
www.kettenbach-dental.com


KETTENBACHDENTAL
Simply intelligent

022740_0919

