



MPF

BRUSH CO



Dental Surgical Hand Tools

MPF Direct

Introducing the MPF Direct composite applicators, these tools are designed with ergonomics, comfort and balance in mind. The unique design of each and every tip allow for intraoral application and

also for the dental technicians lab composite applications as well. The handles are made from anodised aluminium and the tips are stainless steel that have a special BLACK DIAMOND coating.

These tools are all autoclavable and can withstand any sterilisation process. They are NOT to be cleaned by using wire brush techniques as this can destroy the coating process. Care must be taken by using a soft brush to remove any reminiscence of composite material.



Item number: 300-7001



Item number: 300-7002



Item number: 300-7003



Item number: 300-8001



Item number: 300-8002



Item number: 300-8003

KITS



Item number: 300-7000



Item number: 300-8000

MPF Composite Brushes

The MPF COMPOSITE BRUSH series have been created with a special fiber optic filament that is very fine in texture but also firm to assist in application.

The unique design of the MPF Revolution handle brings a completely different ergonomic feel, while it's perfectly balanced and lightweight to the hand.

The tips are to be removed from the handles when sterilization of the handle is needed. All tips of the brushes are to be sterilized utilizing liquid sterilisation methods only. No heating methods should be used for sterilization of the brush tips.



Item number: 300-3001



Item number: 300-3002



Item number: 300-3003



Item number: 300-4001



Item number: 300-4002



Item number: 300-4003

KITS



Item number: 300-3000



Item number: 300-4000

M P F Composite Brushes

The MPF COMPOSITE BRUSH series have been created with a special fiber optic filament that is very fine in texture but also firm to assist in application.

The unique design of the MPF Revolution handle brings a completely different ergonomic feel, while it's perfectly balanced and lightweight to the hand.

The tips are to be removed from the handles when sterilization of the handle is needed. All tips of the brushes are to be sterilized utilizing liquid sterilisation methods only. No heating methods should be used for sterilization of the brush tips.

FLAT



Item number: 300-5001



Item number: 300-5002



Item number: 300-5003



Item number: 300-6001



Item number: 300-6002



Item number: 300-6003

KITS



Item number: 300-5000



Item number: 300-6000

MPF COMPOSITE Replacement Tips

Replaceable composite tips for all of your Composite brushes.
Each tip comes in a pack of 10 pcs.

STRAIGHT



Item number: 300-3001A



Item number: 300-3002A



Item number: 300-3003A



Item number: 300-5001A



Item number: 300-5002A



Item number: 300-5003A




ANGLED



Item number: 300-4001A



Item number: 300-4002A



Item number: 300-4003A



Item number: 300-6001A



Item number: 300-6002A



Item number: 300-6003A



MPF Cube Stand

The MPF Cube stand is a solid piece of anodised aluminium and can be sterilised in the autoclave.



Item number: 300-2000



Item number: 300-1000