

THE DENTAL ATTACHMENTS COMPANY





Live courses

The surgical-prosthetic process in all its phases, an excellent tool to facilitate learning new protocols.

Research and Development

Innovation and research is developed with the collaboration of numerous universities, both Italian and foreign.



Academy Rhein83









FIXED PROSTHESIS OT BRIDGE	
OT EQUATOR abutment, Cuff height measurer, KIT 42, All on Four Set	Pag. 4
Titanium 4 T-base extragrade set, Titanium extragrade T-base	Pag. 5
Castable t-base, Titanium t-base with lateral holes, Healing abutments, Impression copings	Pag. 6
Laboratory analogs, Scan Abutment, Screws, Seeger tool burs, OT EQUATOR screwing tools	Pag. 7
Screwing tools, Universal handle, Seeger inserter, T-Bar Components and accessories	Pag.8
Spare parts, Gingival connectors, Seeger seat cleaning cutter, Seeger tool burs	Pag. 9
REMOVABLE PROSTHESIS	
OT EQUATOR for implants, OT EQUATOR for bars, OT EQUATOR caps kit, OT EQUATOR castable	Pag. 10
Spare parts OT EQUATOR	Pag. 11
OT EQUATOR transfer and analogs, OT EQUATOR tools and accessories	Pag. 12
Single threaded spheres, Spare parts for single threaded spheres	Pag. 13
Sphero block and Sphero flex pfor implants, Cuff height measurer	Pag. 14
Tools Sphero flex and Sphero block, Anti-unscrewing system	Pag. 15
Laboratory set	Pag. 16
Basic laboratory set	Pag. 17
OT Cap / OT Box attachments, OT Strategy attachments, OT Bar attachments	Pag. 18
OT Vertical attachments, OT Unilateral attachments, OT Cap Tecno attachments, OT Lock attachments	Pag. 19
Pivots set for clinic, OT EQUATOR pivots set	Pag. 20
Pivot OT EQUATOR caps set, Pivot Block / Flex	Pag. 21
Spare parts OT Cap normal	Pag. 22
Spare parts OT Cap micro	Pag. 23
OT Cap assortment sets	Pag. 24
Tools and accessories OT Cap normal/micro	Pag. 25
OT Strategy spare parts	Pag. 26
OT Bar multiuse spare parts	Pag. 27
OT Vertical spare parts, OT Reverse spare parts	Pag. 28
OT Lock spare parts, Sphere repair sets	Pag. 29
OT Lever, OT Cem, Broken screw extractor kit	Pag. 30
Rhein83 direct services	Pag. 31

### OT EQUATOR ABUTMENT •



**Ref. 030** 

• 1 OT EQUATOR TITANIUM ABUTMENT WITH TIN COATING

### **CUFF HEIGHT MEASURER**



Ref. 008MBG

### **SET CONTENT**

CUFF HEIGHT SLIDER GAUGE MEASURER

### **KIT 42**



Ref. BE42 + n° 42 OT **EQUATOR** (Ref. 030 to be chosen) • 1 ASSORTMENT SET (autoclavable) + 42 OT EQUATOR +TIN ABUTMENTS THE SET42 CAN BE CUSTOMIZED ON THE BASIS OF THE IMPLANTS BRANDS, IMPLANTS DIAMETERS AND CUFF HEIGHTS DESIRED. PLEASE INDICATE, IN THE ORDER, THE QUANTITY AND THE ABOVE SPECIFICATIONS FOR EVERY KIND OF ABUTMENT. (Tools not included in the set)



Ref. BE42

1 EMPTY BOX FOR ASSORTED OT EQUATOR ABUTMENTS

### ALL ON FOUR OT BRIDGE SET WITH OT EQUATOR ABUTMENTS



THE ALL ON FOUR SET CONTAINS EVERYTHING NEEDED TO BUILD A FOUR-IMPLANT FRA-MEWORK WITH THE OT BRIDGE SYSTEM. THE FOUR OT EQUATOR ABUTMENTS MUST BE OF THE SAME CONNECTION AND HEIGHT. EXTRA-GRADE T-BASE WILL BE COMPATIBLE WITH THE CHO-SEN PLATFORM

### **ALL ON FOUR OT EQUATOR + EXTRAGRADE**

### SET CONTENT

- OT EQUATOR TITANIUM ABUTMENTS WITH TIN COATING TITANIUM EXTRAGRADE T-BASE

- PINK SEEGER WITH HANDLE (SOFT RETENTION)
  WHITE SEEGER WITH HANDLE (STANDARD RETENTION)
  TITANIUM SCREWS FOR EXTRAGRADE T-BASE

#### **ALL ON FOUR OT EQUATOR + EXTRAGRADE MINI SET CONTENT**

Ref. AO4M

Ref. AO4A

- OT EQUATOR TITANIUM ABUTMENTS WITH TIN COATING
- TITANIUM EXTRAGRADE T-BASE MINI PINK SEEGER WITH HANDLE (SOFT RETENTION) WHITE SEEGER WITH HANDLE (STANDARD RETENTION)
- TITANIUM SCREWS FOR MINI EXTRAGRADE TITANIUM T-BASE

### **ALL ON FOUR OT EQUATOR + EXTRAGRADE 15°**

### SET CONTENT

OT EQUATOR TITANIUM ABUTMENTS WITH TIN COATING TITANIUM EXTRAGRADE T-BASE INCLINED AT  $15^\circ$ 

- PINK SEEGER WITH HANDLE (SOFT RETENTION)
  WHITE SEEGER WITH HANDLE (STANDARD RETENTION)
- SPHERICAL TITANIUM SCREWS FOR 15° INCLINED T-BASE

### TITANIUM 4 T-BASE EXTRAGRADE SET



THE 4 T-BASE EXTRAGRADE SET, CONTAINS FOUR COMPLETE EXTRAGRADE T-BASE FOR MAKING A FRAMEWORK WITH THE OT BRIDGE SYSTEM.

**Ref. AO4SEG** 

Ø 3,5 mm

Ref. AO440SEG

Ø4mm

Ref. AO450SEG

Ø5mm

**4 T-BASE EXTRAGRADE** 

**SET CONTENT** 

- 4 4
- TITANIUM EXTRAGRADE T-BASE
  PINK SEEGER WITH HANDLE (SOFT RETENTION)
  WHITE SEEGER WITH HANDLE (STANDARD RETENTION)
  TITANIUM SCREWS FOR EXTRAGRADE T-BASE

**Ref. AO4MSEG** 

Ø 3,5 mm

Ref. AO4M40SEG

Ø4mm

**4 T-BASE EXTRAGRADE MINI** 

**SET CONTENT** 

- TITANIUM EXTRAGRADE T-BASE MINI • 4
- PINK SEEGER WITH HANDLE (SOFT RETENTION)
  WHITE SEEGER WITH HANDLE (STANDARD RETENTION)
- TITANIUM SCREWS FOR MINI EXTRAGRADE TITANIUM T-BASE

4 T-BASE EXTRAGRADE 15°

**Ref. AO4SFEG** 

Ø 3,5 mm

Ref. AO4EG15SF

Ø4 mm

**SET CONTENT** 

- TITANIUM EXTRAGRADE T-BASE INCLINED AT 15° PINK SEEGER WITH HANDLE (SOFT RETENTION) • 4
- WHITE SEEGER WITH HANDLE (STANDARD RETENTION)
  SPHERICAL TITANIUM SCREWS FOR 15° INCLINED T-BASE

### TITANIUM T-BASE FOR FIXED PROSTHESIS



**Ref. MTSEG** 

Ø 3.5 mm

Ref. MT40SEG

Ø4mm

**Ref. MT50SEG** 

Ø5mm

**EXTRAGRADE TITANIUM T-BASE** 

**SET CONTENT** 

- TITANIUM EXTRAGRADE T-BASE TO BE CEMENTED h9,5mm TITANIUM SCREW h2mm
- WHITE SEEGER WITH HANDLE PINK SEEGER WITH HANDLE



**Ref. MMSEG** 

Ø3.5 mm

Ref. MM40SEG

Ø4mm

MINI EXTRAGRADE TITANIUM T-BASE

SET CONTENT

- TITANIUM EXTRAGRADE T-BASE MINI TO BE CEMENTED h5,7mm TITANIUM SCREW FOR MINI EXTRAGRADE TITANIUM T-BASE
- WHITE SEEGER WITH HANDLE PINK SEEGER WITH HANDLE



**Ref. MTSFEG** 

Ø 3.5 mm

Ref. MTEG15SF

**EXTRAGRADE TITANIUM T-BASE INCLINED AT 15°** 

**SET CONTENT** 

- TITANIUM EXTRAGRADE T-BASE INCLINED AT 15° TO BE CEMENTED
- WHITE SEEGER WITH HANDLE
- PINK SEEGER WITH HANDLE SPHERICAL TITANIUM SCREW FOR 15° INCLINED T-BASE

### T-BASE FOR FIXED PROSTHESIS



Ref. CMEG40C

### **EXTRAGRADE CASTABLE T-BASE**

#### **SET CONTENT**

- EXTRAGRADE CASTABLE T-BASE WITH HOLE TITANIUM SCREW h2 mm
- WHITE SEEGER WITH HANDLE
- PINK SEEGER WITH HANDLE



Ref. MCEG15SF

### **EXTRAGRADE CASTABLE T-BASE WITHOUT HOLE**

#### SET CONTENT

- EXTRAGRADE CASTABLE T-BASE WITHOUT HOLE (ONLY AVAILABLE WITH Ø4mm)
- . 2 WHITE SEEGER WITH HANDLE
- PINK SEEGER WITH HANDLE



**Ref. CMTFP** 

h 1.8 mm

Ref. CMTB30FP

h3mm

TITANIUM T-BASE WITH LATERAL HOLES FOR TEMPORARY RESIN BRIDGE WITH WIRE-REINFORCEMENT

#### SET CONTENT

- TITANIUM T-BASE WITH LATERAL HOLES (ONLY AVAILABLE WITH Ø 3,5mm)
- TITANIUM SCREW h2 mm WHITE SEEGER WITH HANDLE
- PINK SEEGER WITH HANDLE



Ref. MT40GL

## TITANIUM T-BASE FOR ENDORAL WELDING

### **SET CONTENT**

- TITANIUM NORMAL T-BASE FOR INTRA ORAL WELDING h9,5 mm Ø4mm TITANIUM SCREW h2 mm
- WHITE SEEGER WITH HANDLE
- PINK SEEGER WITH HANDLE

### **HEALING ABUTMENTS**

Ø4mm



Ref. 159VAG

### **SCREW AND HEALING RING**

### SET CONTENT

- TITANIUM SCREW TITANIUM HEALING RING



Ref. 159EAG

### OT EQUATOR AND HEALING RING

### SET CONTENT

- TITANIUM + TIN OT EQUATOR SCREW FOR HEALING RING
- TITANIUM HEALING RING



Ref. 159VEA

Titanium + TiN OT EQUATOR screw (for self-extracting peek ring or titanium self-extracting healing ring)

### IMPRESSION COPINGS



Ref. 144TLE

• 1 TITANIUM OT EQUATOR LONG IMPRESSION COPING + TITANIUM LONG SCREW h11,5mm



Ref. 144TTE

• 1 TITANIUM OT EQUATOR IMPRESSION COPING + TITANIUM SCREW h6 mm



Ref. 144MTE

• 2 OT EQUATOR IMPRESSION COPINGS FOR INDIVIDUAL TRAY



Ref. TPM2

• 2 MINI OT EQUATOR IMPRESSION COPINGS h3,7mm

### LABORATORY ANALOGS



**Ref. 144AE** 

• 2 LABORATORY ANALOGS ø 4mm

**Ref. AI502** 

• 2 LABORATORY ANALOGS ø 5mm



Ref. 144AVC4 **Ref. 144AVC5** 

- 1 LABORATORY ANALOG WITH SCREW FOR CAD/CAM ø 4mm
- 1 LABORATORY ANALOG WITH SCREW FOR CAD/CAM ø 5mm

### OT EQUATOR SCAN ABUTMENT



Ref. 145SAE

• 1 OT EQUATOR TITANIUM SCAN ABUTMENT + TITANIUM SCREW

### **SCREWS**



Ref. 147VME

• 1 TITANIUM SCREW FOR EXTRAGRADE T-BASE (for screwdriver/stem for inclined Channel)



Ref. 146VMI

SPHERICAL TITANIUM SCREW FOR INCLINED T-BASE • 1 (for screwdriver/stem for inclined Channel)



Ref. 148VCM

TITANIUM SCREW FOR MINI EXTRAGRADE T-BASE (for screwdriver/stem for inclined Channel)



Ref. VC

• 1 TITANIUM SCREW FOR EXTRAGRADE T-BASE h2mm



Ref. VM

• 1 TITANIUM SCREW FOR EXTRAGRADE T-BASE h6mm

Ref. VXL

• 1 TITANIUM SCREW FOR EXTRAGRADE T-BASE h11,5mm

### **TOOLS**



**Ref. KFST** 

### SEEGER TOOL SET FOR DIRECT OT BRIDGE CONNECTION

### **SET CONTENT**

SEEGER TOOL BUR 1 (2.8xR0.05) SEEGER TOOL BUR 2 (2.6x30°)



Ref. 774CHE

**Ref. 775CLE\*** 

(\*Long version)

SQUARE SCREW DRIVER + HOLDER for tightening OT Equator (SQUARE TIP 1.25MM) (USABLE WITH THE MANUAL TORQUE WRENCH)



**Ref. 774HC** 

• 1 INTERCHANGEABLE HOLDER

### **FIXED PROSTHESIS OT BRIDGE**

### **TOOLS**



Ref. 760CSC

Ref. 760CS\* (\*Long version)

SCREWDRIVER/STEM FOR ANGULATED CHANNEL FOR TORQUE CONTROLLER CONTRA ANGLE (SQUARE 1.25MM) AND ADAPTER FOR MANUAL TORQUE WRENCH

Ref. 760ACC

SCREWDRIVER/STEM FOR ANGULATED CHANNEL ADAPTER

(USABLE WITH MANUAL TORQUE WRENCH)



**Ref. 760CE** 

OT EQUATOR SCREWDRIVER/STEM FOR TORQUE CONTROLLER CONTRA ANGLE AND FOR MANUAL TORQUE WRENCH ADAPTER READY FOR HOLDER (SQUARE 1,25mm)

Ref. 760CRD-US

RATCHET TORQUE CONTROL DEVICE (FOR OT EQUATOR, SPHERO BLOCK AND FLEX 15/35 Ncm STRENGHT - MAX 70 Ncm TORQUE.

Ref. 124ICP

• 1 UNIVERSAL HANDLE

**Ref.490IS** 

INSERTER/EXTRACTOR FOR SEEGER FOR FIXED PROSTHESIS WITH TESTER, HANDLE AND CASE

**Ref.190IS** 

• 1 INSERTER AND TESTER FOR FIXED PROSTHESIS SEEGER

### **MEDICAL STEEL T-BAR**



Ref. G80A

4 MEDICAL STEEL JOINTS

Ref. T20A

4 MEDICAL STEEL CYLINDERS

Ref. T20Ai

4 MEDICAL STEEL RETENTIVE CYLINDERS

### TITANIUM T-BAR



Ref. G80T

• 4 TITANIUM JOINTS

Ref. T20T

4 TITANIUM CYLINDERS

Ref. T20Ti

4 TITANIUM RETENTIVE CYLINDERS

### **TOOLS T-BAR**



**Ref. T-APLUS** 

2 T-BAR ADAPTERS

### **SPARE PARTS** -

SPAKE PA	KIS		1
64	Ref. RMTSFEG	• 1 EXTRAGRADE TITANIUM T-BASE INCLINED AT 15° Ø 3,5mm	
	Ref. RMTEG15SF	• 1 EXTRAGRADE TITANIUM T-BASE INCLINED AT15° Ø 4mm	
	Ref. RMTSEG	• 1 EXTRAGRADE TITANIUM T-BASE Ø 3,5mm	
	Ref. RMT40SEG	• 1 EXTRAGRADE TITANIUM T-BASE Ø 4mm	
	Ref. RMT50SEG	• 1 EXTRAGRADE TITANIUM T-BASE Ø 5mm	
<b>Q</b> .	Ref. RMMSEG	• 1 MINI EXTRAGRADE TITANIUM T-BASE Ø 3,5mm	
	Ref. RMM40SEG	• 1 MINI EXTRAGRADE TITANIUM T-BASE Ø 4mm	
<b>a</b>	Ref. GCMR	•1 CASTABLE SLEEVE FOR INCLINED TITANIUM T-BASE WITHOUT HOLE Ø 3,5mm	
1	Ref. GCMR4	• 1 CASTABLE SLEEVE FOR INCLINED TITANIUM T-BASE WITHOUT HOLE Ø4mm	
	Ref. GC	• 1 CASTABLE SLEEVE FOR TITANIUM T-BASE Ø 3,5mm	
0	Ref. GC4	•1 CASTABLE SLEEVE FOR TITANIUM T-BASE Ø 4mm	
	Ref. GC5	• 1 CASTABLE SLEEVE FOR TITANIUM T-BASE Ø 5mm	
0.	Ref. GMM	• 1 CASTABLE SLEEVE FOR MINI TITANIUM T-BASE Ø 3,5mm	
0	Ref. GMM4	• 1 CASTABLE SLEEVE FOR MINI TITANIUM T-BASE Ø 4mm	
Ele	Ref. MTFP	• 1 TITANIUM T-BASE WITH LATERAL HOLES FOR TEMPORARY RESIN BRIDGE WITH WIRE-REINFORCEMENT h 1,8 mm (ONLY AVAILABLE WITH Ø 3,5mm)	
	Ref. MTB30FP	• 1 TITANIUM T-BASE WITH LATERAL HOLES FOR TEMPORARY RESIN BRIDGE WITH WIRE-REINFORCEMENT h3 mm (ONLY AVAILABLE WITH Ø 3,5mm)	
	Ref. HR800	• 1 TITANIUM WIRE WITHOUT ELASTIC MEMORY Ø 1,0mm	
	Ref. MEG40C	• 1 EXTRAGRADE CASTABLE T-BASE Ø 4mm	
8	Ref. RMCEG15SF	• 1 EXTRAGRADE CASTABLE T-BASE WITHOUT HOLE (ONLY AVAILABLE WITH Ø 4,0mm)	
	Ref. SM6	• 6 WHITE SEEGER WITH HANDLE (STANDARD RETENTION)	
1,	Ref. SRM6	• 6 PINK SEEGER WITH HANDLE (SOFT RETENTION)	
	Ref. 022RGO	• 3 CASTABLE CONNECTION BAR	
-WEW	Ref. FST1	• 1 SEEGER TOOL BUR 1 (2.8xR0.05) FOR DIRECT OT BRIDGE CONNECTION	
M. James	Ref. FST2	• 1 SEEGER TOOL BUR 2 (2.6x30°) FOR DIRECT OT BRIDGE CONNECTION	
	Ref. FSS	• 1 BUR FOR SEEGER HOUSING POLISHING	
	Ref. OC	• 1 OT CEM (COMPOSITE CEMENT: 1 Two-component syringe , 10 AUTOMIX TIP5)	
111	Ref. TOC	• 10 AUTOMIX TIPS	

**Ref. 130** 

Ref. 131

### OT EQUATOR FOR IMPLANTS =



### **OT EQUATOR FOR IMPLANTS SET**

#### SET CONTENT

- OT EQUATOR TITANIUM + TIN ABUTMENT
- STAINLESS STEEL HOUSING PROTECTIVE DISK
- RETENTIVE CAPS
  - (1 EXTRA-SOFT, 1 SOFT, 1 STANDARD, 1 STRONG)



### OT EQUATOR AND SMARTBOX FOR IMPLANTS SET

### **SET CONTENT**

- OT EQUATOR TITANIUM + TIN ABUTMENT ANODIZED SMART BOX HOUSING WITH BLACK POSITIONING CAP
- PINK PROTECTIVE DISK
  - (1 EXTRA-SOFT, 1 SOFT, 1 STANDARD, 1 STRONG)



**Ref. 030** 

OT EQUATOR TITANIUM + TIN ABUTMENT

THE OT EQUATOR ATTACHMENTS ARE AVAILABLE FOR THE MAJOR IMPLANT BRANDS. UPON REQUEST WE CAN CUSTOMIZE THEM FOR ANY KIND OF IMPLANTS. WHEN ORDERING INDICATE: IMPLANT, DIAMETER AND CUFF HEIGHT. CUFF AVAILABLE TO 7mm.

### OT EQUATOR FOR BARS -



### **OT EQUATOR THREADED SLEEVE SYSTEM SET FOR BONDING**

#### SET CONTENT

- OT EQUATOR THREADED for titanium sleeve (1.6 mm THREAD)
- THREADED SLEEVES (1.6 mm THREAD) STAINLESS STEEL HOUSINGS • 2
- OT EOUATOR SPACERS
- RETENTIVE CAPS
  - (2 EXTRA-SOFT, 2 SOFT, 2 STANDARD, 2 PROCESSING)

### **ELASTIC SEEGER** PASSIVE BAR CONNECTIONS SYSTEM



Ref. 158ESA

**Ref. 160EQB** 

### **ELASTIC SEEGER SET**

- **SET CONTENT**
- CASTABLE CYLINDER FOR SEEGER
- SELF-EXTRACTING SEEGER TITANIUM SCREW FOR SELF-EXTRACTING SEEGER

### CAPS ASSORTMENT SET •



**Ref. 192ECE** 



### **CAPS ASSORTMENT SET** SET CONTENT

- STAINLESS STEEL HOUSING
- PROCESSING CAP
- RETENTIVE CAPS (1 EXTRA-SOFT, 1 SOFT, 1 STANDARD, 1 STRONG)

### **SMART BOX SET** SELF-ALIGNING OT EQUATOR HOUSING



### **SMART BOX SET**

### SET CONTENT

- SMART BOX HOUSING WITH BLACK POSITIONING CAP
- PINK PROTECTIVE DISK
- RETENTIVE CAPS (1 EXTRA-SOFT, 1 SOFT, 1 STANDARD, 1 STRONG)

### OT EQUATOR CASTABLE —





Ref. 092ECQ

### **OT EQUATOR CASTABLE SET**

### **SET CONTENT**

- SINGLE MALES STAINLESS STEEL HOUSINGS
- RETENTIVE CAPS (2 SOFT, 2 STANDARD)

### SPARE PARTS OT EQUATOR —

JI AILL I A	MIS OI EQUAIC	<b>)</b> (1		
	Ref. 140CEV	•4	RETENTIVE CAPS - VIOLET "STRONG" 2700g	
	Ref. 140CET	•4	RETENTIVE CAPS - CLEAR "STANDARD" 1800g	
	Ref. 140CER	•4	RETENTIVE CAPS - PINK "SOFT" 1200g	
	Ref. 140CEG	•4	RETENTIVE CAPS - YELLOW (EXTRA-SOFT 600g)	
	Ref. 140CEN	• 4	PROCESSING CAPS - BLACK (FOR LABORATORY USE ONLY)	
	Ref. 140CEDR8	• 4	RETENTIVE CAPS - AQUA CAPS UNDERSIZED INTERNAL DIAMETER 2,2mm (STRONG 1300g)	
	Ref. 140CEADR8	•4	RETENTIVE CAPS - ORANGE UNDERSIZED INTERNAL DIAMETER 2,2mm (VERY ELASTIC 350g)	
	Ref. 141CAE	• 2	STAINLESS STEEL HOUSINGS	
	Ref. 141CTE	• 2	TITANIUM HOUSINGS FUCHSIA ANODIZED	
	Ref. 330SBE	• 1	SMART BOX HOUSING (WITH BLACK POSITIONING CAP)	
	Ref. 335CSB	• 4	BLACK SMARTBOX (POSITIONING CAPS)	
<del></del>	Ref. 039SFE2	• 1	OT EQUATOR THREADED 2mm (DIRECT INTO 2mm THREAD)	
	Ref. 339SFE	• 1	OT EQUATOR THREADED 1.6mm (FOR TITANIUM SLEEVE)	
	Ref. 239GSF	• 1	TITANIUM SLEEVE (1.6mm THREAD)	
	Ref. 339SSE	• 1	OT EQUATOR SPACER	
	Ref. 159VEA	• 1	OT EQUATOR SEEGER (FOR SELF-EXTRACTING PEEK RING OR TITANIUM HEALING RING)	
*	Ref. 158VAE	• 1	TITANIUM SCREW FOR SELF-EXTRACTING SEEGER	
	Ref. 158SAE	• 1	SELF-EXTRACTING WHITE SEEGER (FOR BAR LOCKING)	
	Ref. 158CCE	• 1	CASTABLE CYLINDER height 2,5mm (FOR SEEGER)	
	Ref. 158CCAE	• 1	HIGH CASTABLE CYLINDER height 3,5mm (FOR SEEGER)	
3	Ref. 151SS	•2	CASTABLE OT EQUATOR MALES	

### **ACCESSORIES OT EQUATOR**



Ref. 144TLE • 1 OT EQUATOR TITANIUM TRANSFER + TITANIUM SCREW



Ref. 144TTE • 1 OT EQUATOR TITANIUM TRANSFER + TITANIUM SCREW



Ref. 144MTE • 2 IMPRESSION COPINGS (FOR INDIVIDUAL TRAY)



IMPRESSION COPINGS NORMAL SIZE - OT EQUATOR Ref. 044CAIN

(WITH INTERCHANGEABLE CAP FOR PICK UP IMPRESSION)



Ref. 144AE • 2 LABORATORY ANALOGS ø 4mm

• 2 LABORATORY ANALOGS ø 5mm **Ref. AI502** 

• 1 LABORATORY ANALOG WITH TITANIUM SCREW FOR CAD/CAM ø4mm **Ref. 144AVC4** 

**Ref. 144AVC5** LABORATORY ANALOG WITH TITANIUM SCREW FOR CAD/CAM ø5mm

Ref. 161SCME OT EQUATOR SPACERS FOR HOUSING

### OT EQUATOR TOOLS



Ref. 74AC01 • 1 PARALLELOMETER KEY "NORMAL SIZE"



OT FOUATOR SOLIARE SCREW DRIVER + HOLDER Ref. 774CHE (SQUARE 1.25MM) (USABLE WITH RATCHET TORQUE CONTROL DEVICE)



**Ref. 774HC** OT EQUATOR INTERCHANGEABLE HOLDER



SQUARE DRIVER CONNECTOR FOR TORQUE **Ref. 760CE** 

CONTROLLER CONTRA ANGLE (SQUARE 1.25mm)



MULTIUSE CAPS INSERTER/EXTRACTOR TOOL **Ref. 485IC** • 1



**Ref. 185SIS** METAL INSERTION TOOL FOR SEEGER (USE TOGETHER WITH 491EC)



MULTIUSE EXTRACTOR TOOL FOR CAPS

(WITH HOLE STEM HOLDER)



RATCHET TORQUE CONTROL DEVICE (FOR OT EQUATOR, SPHERO BLOCK AND FLEX 15/35 Ncm STRENGHT - MAX 70 Ncm TORQUE.

### SINGLE THREADED SPHERES WITH THREADED SLEEVE TO BE BONDED NORMAL - MICRO



Ref. 139KSFN

## SET SINGLE THREADED SPHERE NORMAL SIZE + TIN WITH SLEEVE TO BE BONDED

#### **SET CONTENT**

- 2 TITANIUM NORMAL SINGLE THREADED SPHERE (HEX 1,3mm THREAD 1,6mm)
- 2 TITANIUM BOND-IN THREADED SLEEVE
- 2 S/STEEL SPACER for the THREADED SLEEVE (NORMAL sphere)

## SET SINGLE THREADED SPHERE MICRO SIZE + TIN WITH SLEEVE TO BE BONDED

### Ref. 139KSFM SET CONTENT

- 2 TITANIUM MICRO SINGLE THREADED SPHERE (HEX 0.9mm THREAD 1.6mm)
- (HEX 0,9mm THREAD 1,6mm)
  2 TITANIUM BOND-IN THREADED SLEEVE
- 2 S/STEEL SPACER for the THREADED SLEEVE (MICRO sphere)

### **SPARE PARTS**



Ref. 239SFN

1 SINGLE THREADED SPHERE NORMAL SIZE + TIN FOR THREADED SLEEVE HEX 1,3mm THREAD 1,6mm



Ref. 239SFM

1 SINGLE THREADED SPHERE MICRO SIZE + TIN FOR THREADED SLEEVE HEX 0,9mm THREAD 1,6mm



Ref. 239GSF

1 TITANIUM BOND-IN THREADED SLEEVE FOR NORMAL AND MICRO - 1,6mm



Ref. 239SSN

• 1 STAINLESS STEEL SPACER FOR THREADED SLEEVE (NORMAL sphere)



Ref. 239SSM

• 1 STAINLESS STEEL SPACER FOR THREADED SLEEVE (Micro sphere)



**Ref. 039SFN2** 

1 SINGLE THREADED SPHERE NORMAL SIZE + TIN HEX 1,3mm THREAD 2mm



Ref. 039SFM2

 SINGLE THREADED SPHERE MICRO SIZE + TIN HEX 0,9mm THREAD 2mm



Ref. 760CBR

1 HEX DRIVER FOR THREADED SPHERE NORMAL HEX 1,27mm



Ref. 772CSF

1 HEX DRIVER FOR THREADED SPHERE MICRO HEX 0,9mm

### **OT CAP** IMPRESSION COPINGS



Ref. 044CAIN

2 IMPRESSION COPINGS NORMAL SIZE Ø2,5mm WITH INTERCHANGEABLE CAP



Ref. 044CAIM

• 2 IMPRESSION COPINGS MICRO SIZE Ø1,8mm WITH INTERCHANGEABLE CAP

### IMPLANTOLOGY TITANIUM ATTACHMENTS FOR OVERDENTURE —

SPHERO FLEX - SPHERO BLOCK

**Ref. 109** 



### **SPHERO FLEX TIN COATED**

### **SET CONTENT**

- TITANIUM ABUTMENT WITH SWIVEL SPHERE
- (SPHERE Ø 2,5mm) CAPS (Yellow Pink Clear) STAINLESS STEEL HOUSING
- PROTECTIVE DISK DIRECTIONAL RINGS

Ref. 002



### **SPHERO BLOCK NORMAL TIN COATED**

### **SET CONTENT**

- TITANIUM ABUTMENT
- (SPHERE Ø 2,5 mm) CAPS (Yellow Pink Clear)
- STAINLESS STEEL HOUSING PROTECTIVE DISK
- DIRECTIONAL RINGS

Ref. 003



### SPHERO BLOCK MICRO TIN COATED

#### **SET CONTENT**

- TITANIUM ABUTMENT (SPHERE Ø 1,8mm)
- CAPS (Yellow Pink Clear) STAINLESS STEEL HOUSING PROTECTIVE DISK • 3
- DIRECTIONAL RINGS

THE SPHERO FLEX, SPHERO BLOCK ATTACHMENTS ARE AVAILABLE FOR THE MAJOR IMPLANT BRANDS.UPON REQUEST WE CAN CUSTOMIZE THEM FOR ANY KIND OF IMPLANTS. WHEN (B.D.G) mm ORDERING, PLEASE ALWAYS INDICATE: IMPLANT BRAND, IMPLANT DIAMETER AND CUFF HEIGHT. FOR SPHERO FLEX ABUTMENT CUFF AVAILABLE FROM 0,5 TO 7mm.



### **CUFF HEIGHT MEASURER**

THIS TOOL IS USED TO DETERMINE THE **CUFF HEIGHT OF IMPLANTS WITH** INTERNAL OR EXTERNAL HEXAGON




### **SET CONTENT**

CUFF HEIGHT SLIDER GAUGE MEASURER



### **TOOLS SPHERO FLEX - SPHERO BLOCK**



**ELASTIC DOWEL** 

### ANTI-UNSCREWING SYSTEM WITH AN ELASTIC DOWEL APPLICABLE TO

SPHERO FLEX-BLOCK - OT REVERSE AND DIFFERENT SCREWS. THE ELASTIC DOWEL IS APPLICABLE ONLY UPON REQUESTTO ALL SCREWS WITH A NOT INFERIOR M. 1.8 to screws - Sphero Flex - Sphero Block - OT Reverse attachments, etc.

Extra charge for the dowel application



### ATTACHMENTS OT CAP & OT BOX - OT STRATEGY - OT BAR OT VERTICAL - OT TECNO - TITANIUM SINGLE SPHERE

MASTER STARTER SET FOR LABORATORY



#### **SET CONTENT**

#### OT CAP - OT BOX SYSTEM:

- 16 castable assorted pivots (NORMAL and micro size)
   4 castable spheres (NORMAL and micro size)
   2 titanium single spheres (NORMAL and micro size)
   + sphere holder + 2 sliding base
   4 castable OT CAP bars (NORMAL and micro size)
   2 complete castable OT CAP TECNO bars (NORMAL and micro size)

- castable OT BOX CLASSIC bars (upper + lower part) NORMAL and micro size

- castable OT BOX SPECIAL bars ( NORMAL and micro size) castable OT BOX connectors castable OT BOX MONO housings (NORMAL and micro size)
- plastic positioner rings (NORMAL and micro size)
  OT CAP retentive caps (NORMAL and micro size)
- Clear- standard retention PINK soft retention YELLOW extra soft retention - GREEN very elastic retention TITAN CAP titanium caps (NORMAL and micro size)
- BLACK CAPS- only for lab use (NORMAL and micro size)
- S/STEEL housings • 8
- ANALOGUES for plaster models (NORMAL and micro size)
- protective disks
- aluminium spacers (NORMAL and micro size)

### OT STRATEGY:

- castable spheres (2 standard and 2 high) castable STEADY
- retentive caps: 2 clear 2 pink 2 yellow for housing
- black caps for housing only for lab use
- s/steel housings
- plastic positioner rings
- retentive caps for duplication technique: CLEAR standard retention PINK soft retention YELLOW extra soft retention
- black caps for duplication technique only for lab use
- ANALOGUES for plaster models

### **OT BAR MULTIUSE:**

- castable bar
- gingival connector
- positioner clip A
- positioner clip B
- castable housings
- clips 2 PINK soft retention 2 YELLOW medium retention

### OT VERTICAL:

Ref. 005SKLUS

- castable male • 2
- castable STEADY
- clips 2 WHITE- standard retention 2 GREEN soft retention
- castable parallelometer keys
- ceramic bars

### TOOLS:

- Pincers 2 Reamer tools for caps (NORMAL and micro size)

- Reamer tools for Caps (NORMAL and micro size)
  Reamer tools for OT STRATEGY caps
  Parallelometer keys for OT CAP (NORMAL and micro size)
  Parallelometer key for OT STRATEGY size
  Parallelometer keys for OT CAP TECNO (NORMAL and micro size)
  Parallelometer keys for OT BAR MULTIUSE
  Blue Plastic "multiuse" insertion tool
- (for steel tool and for inserting the caps in the mouth)
- OT CAP insertion tool
- OT STRATEGY insertion tool
- OT BAR clips insertion tool
- OT VERTICAL clips insertion tool

### ATTACHMENTS OT CAP & OT BOX - OT STRATEGY- OT BAR OT VERTICAL - OT UNILATERAL - OT EOUATOR

**BASIC LABORATORY SET** 



### OT CAP - OT BOX SYSTEM:

- 16 castable assorted pivots (NORMAL and micro size)
- castable spheres NORMAL/micro
- castable OT CAP bars (NORMAL and micro size) castable OT BOX CLASSIC bars ( upper + lower part) NORMAL and micro size castable OT BOX SPECIAL bars ( NORMAL and micro size)
- castable OT BOX connectors
- castable OT BOX MONO housings (NORMAL and micro size)
- plastic positioner rings (NORMAL and micro size)
- OT CAP assorted retentive caps (NORMAL and micro size) PINK soft retention - YELLOW extra soft retention
- 13 BLACK CAPS- only for lab use (NORMAL and micro size)
  4 S/STEEL housings (NORMAL and micro size)

- castable spheres (2 standard and 2 high) castable STEADY
- pink retentive caps
- s/steel housings
- plastic positioner rings
- retentive caps for duplication technique: PINK soft retention YELLOW extra soft retention
- 4 black caps for duplication technique only for lab use

### **OT EQUATOR:**

- single males
- stainless steel housings
- retentive caps (2 soft, 2 standard)
- processing black caps

### **OT UNILATERAL:**

- unilateral castable attachment
- uni-box castable
- plastic positioner micro size
- ot cap micro size: 1 pink 1 black processing cap
- ot strategy (for duplication technique): 1 yellow 1 black processing cap
- connector

### **OT BAR MULTIUSE:**

- castable bar
- gengival connector
- positioner clip (A B)
- clips 2 PINK soft retention 2 YELLOW medium retention

### OT VERTICAL:

- castable male
- castable STEADY
- clips 2 WHITE- standard retention 2 GREEN soft retention
- castable parallelometer keys
- 2 ceramic bars

### TOOLS:

Ref. 005SKLBUS

- Tweezer
- Parallelometer keys for OT CAP NORMAL
- Parallelometer keys for OT CAP micro
- Parallelometer key for OT STRATEGY
- Parallelometer key for OT BAR MULTIUSE . 1 Blue Plastic "multiuse" insertion tool
- (for steel tool and for inserting the caps in the mouth)
- 1 OT CAP insertion tool
- OT STRATEGY insertion tool • 1
- OT BAR clips insertion tool
- OT VERTICAL clips insertion tool

### OT CAP AND OT BOX ATTACHMENTS



NORMAL

Ref. 092CAN

**Ref. 196BCN** 

Ref 197BCM

**OT CAP NORMAL/MICRO SET CONTENT** 

MICRO

• 4 CASTABLE SINGLE SPHERES PINK RETENTIVE CAPS

Ref. 092CAM

STAINLESS STEEL HOUSINGS

PLASTIC POSITIONERS



NORMAL

MICRO

**OT CAP NORMAL/MICRO "ECONOMIC"** 

**SET CONTENT** 

CASTABLE PRESECTIONED BAR
CASTABLE PRESECTIONED BEVELED BAR
CLEAR RETENTIVE CAPS (Standard retention)



NORMAL

Ref. 099BSN

OT CAP & OT BOX MONO NORMAL/MICRO

**SET CONTENT** 

CASTABLE PRESECTIONED BARS (1 straight + 1 beveled) • 2

MICRO CLEAR RETENTIVE CAPS
CASTABLE MONO OT BOX **Ref. 099BSM** 

PLASTIC POSITIONERS



NORMAL

OT BOX SPECIAL NORMAL/MICRO + CONNECTORS Ref. 058BSN

SET CONTENT

2 OT BOX SPECIAL BARS

4 PLASTIC POSITIONERS MICRO Ref. 058BSM

CONNECTORS



NORMAL

MICRO

OT BOX CLASSIC NORMAL/MICRO + CONNECTORS

**SET CONTENT** 

• 2 UPPER BARS LOWER BARS

PLASTIC POSITIONERS

CONNECTORS • 4

### OT STRATEGY ATTACHMENTS -



**Ref. 098SSS** 

**Ref. 153BCN** 

Ref. 153BCM

**OT STRATEGY** CAPS FOR DUPLICATION TECHNIQUE

SET CONTENT

CASTABLE STRATEGY MALES (2 STANDARD - 2 HIGH) CASTABLE STEADY

RETENTIVE CAPS FOR DUPLICATION TECHNIQUE



Ref. 098CAL

**OT STRATEGY** CAPS FOR STAINLESS STEEL HOUSING TECHNIQUE

SET CONTENT

CASTABLE STRATEGY MALES (2 STANDARD - 2 HIGH) RETENTIVE CAPS FOR HOUSING TECHNIQUE (2 SOFT, 2 STANDARD) STAINLESS STEEL HOUSINGS • 4 . 4

PLASTIC POSITIONERS

CASTABLE STEADY

### OT BAR MULTIUSE ATTACHMENTS =



**Ref. 0210BM** 

**OT BAR MULTIUSE + CONNECTORS** 

SET CONTENT

CASTABLE BARS POSITIONER CLIP A · 2

POSITIONER CLIP B

CASTABLE BOXES
RETENTIVE PINK CLIPS

RETENTIVE YELLOW CLIPS

RETENTIVE RED CLIPS

CONNECTORS

GINGIVAL CONNECTOR

### OT VERTICAL ATTACHMENTS



**Ref. 0710BV** 

### **OT VERTICAL**

#### **SET CONTENT**

- CASTABLE MALES • 4 • 2
- CASTABLE STEADY
  RETENTIVE WHITE CLIPS
  RETENTIVE GREEN CLIPS
- CASTABLE PARALLELOMETER KEYS + PIN

### **OT UNILATERAL** EXTRA CORONAL ATTACHMENTS WITH COMBINED FUNCTION

### Ref. 098UNL

### **OT UNILATERAL + CONNECTOR**

#### **SET CONTENT**

- UNILATERAL CASTABLE ATTACHMENT UNI-BOX CASTABLE
- PLASTIC POSITIONER MICRO SIZE PINK CAP MICRO SIZE

- STRATEGY YELLOW CAP (for duplication technique) BLACK CAP MICRO SIZE (for laboratory use only)
- STRATEGY BLACK CAP (for laboratory use only)
- CONNECTOR



Ref. 156BCU

• 1 UNILATERAL CASTABLE ATTACHMENT



Ref. 156CCU

• 1 UNI-BOX CASTABLE

### OT CAP TECNO =

### NORMAL

Ref. 093CTN

### **OT CAP TECNO NORMAL/MICRO SET CONTENT**

- Ref. 093CTM
- PRE-ANGULATED CASTABLE EXTENSIONS TITANIUM SINGLE THREADED SPHERES
- TITANIUM BOND-IN THREADED SLEEVES TRANSPARENT CAPS NORMAL/MICRO
- PINK CAPS NORMAL/MICRO YELLOW CAPS NORMAL/MICRO



Ref. 193BCA

• 2 PRE-ANGULATED CASTABLE EXTENSIONS

### **OT LOCK** LOCKING PIN TITANIUM AND CASTABLE



**Ref. 880CLT** 

### **OT LOCK SET TITANIUM**

### **SET CONTENT**

- LOCKING PIN SET
- HOUSING SHAPER CERAMIC PIN



Ref. 880CLR

### **OT LOCK SET TITANIUM ADJUSTABLE**

### **SET CONTENT**

- LOCKING PIN SET
- HOUSING SHAPER CERAMIC PIN
- CASTABLE SPACER RINGS

### OVERDENTURE ATTACHMENTS

### S.P.L. PIVOTS - COPING COVER

PROMOTIONAL SET FOR CLINIC



Ref. 011SPL

### **SET CONTENT**

- 24 CASTABLE ASSORTED PIVOTS MICRO SIZE
- 40 CASTABLE ASSORTED PIVOTS NORMAL SIZE
- 12 TITANIUM MICRO PIVOTS
- 12 TITANIUM NORMAL PIVOTS
- 16 NORMAL CAPS (Clear standard retention PINK soft retention)
- 16 MICRO CAPS (Clear standard retention PINK soft retention)
- 20 PROTECTIVE DISKS
- 2 REGULATING TOOLS (1 NORMAL and 1 micro size)
- BLUE PLASTIC "MULTIUSE" INSERTION TOOL • 1 (FOR THE INSERTION OF THE CAPS INSIDE THE METAL HOUSINGS)
- OT CAP insertion tool • 1
- NORMAL STAINLESS STEEL HOUSING MICRO STAINLESS STEEL HOUSING • 4
- 4
- 12 DIRECTIONAL RINGS (grey 0° green 7° red 14°)
- TWEEZER

#### Lengths available: 10 - 9 - 7mm

**Ref. 333PEK7** 

**Ref. 333PEK9** 

Ref. 333PEK10

### OT EQUATOR TITANIUM TIN COATED PIVOTS SET

### **LOW PROFILE HEAD PIVOTS**

SET CONTENT

- TITANIUM OT EQUATOR PIVOTS h7,9 or 10mm (of the same length chosen)
  OT EQUATOR ANODIZED TITANIUM HOUSINGS
  YELLOW EXTRA SOFT CAPS
- 2
- PINK SOFT CAPS CLEAR STANDARD CAPS
- TRANSPARENT PROTECTIVE DISKS

Lengths available: 10 - 9 - 7mm

**Ref. 336PTE7** 

**Ref. 336PTE9** 

Ref. 336PTE10

#### OT EQUATOR TITANIUM TIN COATED PIVOTS LOW PROFILE HEAD PIVOTS SET CONTENT

• 1 TITANIUM OT EQUATOR PIVOT h7,9 or 10mm (of the length chosen)






**Ref. 192ECE** 

### OT EQUATOR CAPS ASSORTMENT SET

### **SET CONTENT**

- STAINLESS STEEL HOUSING
- PROCESSING CAP
- RETENTIVE CAPS (1 EXTRA-SOFT, 1 SOFT, 1 STANDARD, 1 STRONG)
- PROTECTIVE DISK



Ref. 033PSF7R

Ref. 033PSF9R

Ref. 033PSF10R

**PIVOT FLEX TIN COATED** 

#### TITANIUM PIVOT WITH SWIVEL SPHERE - NORMAL size Adapted for COPING COVER SET CONTENT

- TITANIUM PIVOT WITH SWIVEL SPHERE
- STAINLESS STEEL HOUSING
- PINK CAPS NORMAL Size Soft retention
- PROTECTIVE DISK
- DIRECTIONAL RINGS

Available in lengths: 10 - 9 - 7mm NORMAL SPHERE ø2,5 MICRO SPHERE ø1,8

**Ref. 036PTN7** 

**Ref. 036PTN9** 

TITANIUM PIVOTS SET NORMAL SIZE **SET CONTENT** 

TITANIUM PIVOTS (Adapted for COPING COVER) • 5 Ref. 036PTN10

**Ref. 036PTM7** 

**Ref. 036PTM9** 

Ref. 036PTM10

TITANIUM PIVOTS SET MICRO SIZE

**SET CONTENT** 

TITANIUM PIVOTS (Adapted for COPING COVER)

NORMAL

Ref. 192ACN

MICRO

CAPS ASSORTMENT SET OT CAP NORMAL/MICRO SET CONTENT

Ref. 192ACM

STAINLESS STEEL HOUSING

PROCESSING CAP RETENTIVE CAPS (1 GREEN, 1 YELLOW, 1 PINK, 1 CLEAR)

PROTECTIVE DISK



Ref. 034PRK

**OT REVERSE 3 OVERDENTURE POSTS** FOR ROOTS, AVAILABLE IN ONE STANDARD SIZE

SET CONTENT

ROOT POSTS IN TITANIUM + TIN RETENTIVE MALE ATTACHMENT TITANIUM + NYLON

PLASTIC MANUAL TOOLS

STAINLESS STEEL HOUSINGS

PINK PROTECTIVE DISKS

Ø of sub-sphere pin: 2,2 2 1,8 1,7 1,4mm



Ref. 010PSP

**CASTABLE PIVOTS** FOR CANAL IMPRESSIONS SET CONTENT

20 CASTABLE PIVOTS NORMAL SIZE

Lenghts available: 14 - 12 - 10 - 9 - 7mm



Ref. 012PSM

**CASTABLE PIVOTS FOR CANAL IMPRESSIONS** 

SET CONTENT • 20 CASTABLE PIVOTS MICRO SIZE

Lenghts available: 10 - 9 - 7mm



Ref. A01MOG

**MOOSER BUR FOR PIVOTS OF 7 - 9 - 10mm** 



Ref. A03MOB

**MOOSER BUR FOR PIVOTS OF 12 - 14mm** 



**Ref. 121FD** 

**OT REVERSE 3** DIAMOND SHAPED BUR FOR ROOTS



**Ref. 100AD** 

**DIRECTIONAL RINGS FOR SPHERO FLEX AND BLOCK** COLOURED RINGS WITH INCLINATION OF 0° - 7'

**Ref. 485IC** 

• 1 EXTRACTOR/INSERTER TOOL FOR CAPS OT EQUATOR - NORMAL - MICRO

## SPARE PARTS OT CAP SYSTEM NORMAL Ø 2,5 mm

-	Ref. 055SCN	• 4	CASTABLE SINGLE SPHERES NORMAL SIZE	
Was .	Ref. 151BAN	• 2	CASTABLE PRESECTIONED BEVELED BARS NORMAL SIZE	
1	Ref. 150BPN	• 2	CASTABLE PRESECTIONED BARS NORMAL SIZE	
	Ref. 048CAN	• 6	SILVER CAPS NORMAL SIZE (EXTRA-RESILIENT) (ELASTIC AND GUMMY) - 350g	
	Ref. 048CON	• 6	GOLD CAPS NORMAL SIZE (EXTRA-RESILIENT) (SLIGHTLY ELASTIC) - 500g	
	Ref. 049PCN	•6	GREEN CAPS NORMAL SIZE (VERY ELASTIC RETENTION) - 350g	
	Ref. 060CRN AY	•6	YELLOW CAPS NORMAL SIZE (EXTRA SOFT RETENTION) - 500g	
	Ref. 040CRN SN	• 6	PINK CAPS NORMAL SIZE (SOFT RETENTION) - 900g	
	Ref. 040CRN	• 6	CLEAR CAPS NORMAL SIZE (STANDARD RETENTION) - 1300g	
	Ref. 040TCN	• 2	TITAN CAPS NORMAL SIZE	
	Ref. 043CLN	• 6	BLACK CAPS NORMAL SIZE (FOR LABORATORY USE ONLY)	
	Ref. 049PCNDR8	•6	ORANGE CAPS NORMAL SIZE (VERY ELASTIC RETENTION) (UNDERSIZED INTERNAL DIAMETER FOR WORN OUT SPHERE OR FOR Ø 2.25 SPHERE) - 350g	
	Ref. 060CRNAYDR8	•6	YELLOW CAPS NORMAL SIZE (EXTRA SOFT RETENTION) (UNDERSIZED INTERNAL DIAMETER FOR WORN OUT SPHERE OR FOR Ø 2.25 SPHERE) - 500g	
	Ref.040CRNSNDR8	•6	PINK CAPS NORMAL SIZE (SOFT RETENTION) (UNDERSIZED INTERNAL DIAMETER FOR WORN OUT SPHERE OR FOR Ø 2.25 SPHERE) - 900g	
	Ref. 040CRNDR8	•6	AQUA CAPS NORMAL SIZE (STRONG RETENTION) (UNDERSIZED INTERNAL DIAMETER FOR WORN OUT SPHERE OR FOR Ø 2.25 SPHERE)) - 1300g	
Ò	Ref. 042AMN	• 2	MONO OT BOX NORMAL SIZE (CASTABLE HOUSINGS FOR CAPS)	
	Ref. 041CAN	• 2	STAINLESS STEEL HOUSINGS NORMAL SIZE	
	Ref. 041CTN	• 2	TITANIUM HOUSINGS NORMAL SIZE	
	Ref. 101PPN	• 2	PLASTIC POSITIONERS NORMAL SIZE	
	Ref. 044PPN	• 2	NORMAL SPHERE ANALOGS	

## SPARE PARTS OT CAP SYSTEM MICRO Ø 1,8 mm

	Ref. 055SCM	• 4	CASTABLE SPHERES MICRO SIZE	
Sand I	Ref. 151BAM	• 2	CASTABLE PRESECTIONED BEVELED BARS MICRO SIZE	
No.	Ref. 150BPM	• 2	CASTABLE PRESECTIONED BARS MICRO SIZE	
	Ref. 048CAM	• 6	SILVER CAPS MICRO SIZE (EXTRA-RESILIENT) (ELASTIC AND GUMMY) - 200g	
	Ref. 048COM	•6	GOLD CAPS MICRO SIZE (EXTRA-RESILIENT) (SLIGHTLY ELASTIC) - 450g	
	Ref. 049PCM	•6	GREEN CAPS MICRO SIZE (VERY ELASTIC RETENTION) - 200g	
	Ref. 060CRM AY	•6	YELLOW CAPS MICRO SIZE (EXTRA SOFT RETENTION) - 450g	
	Ref. 040CRM SN	• 6	PINK CAPS MICRO SIZE (SOFT RETENTION) - 800g	
	Ref. 040CRM	•6	CLEAR CAPS MICRO SIZE (STANDARD RETENTION) - 1100g	
	Ref. 040TCM	• 2	TITAN CAPS MICRO SIZE	
	Ref. 043CLM	• 6	BLACK CAPS MICRO SIZE (FOR LABORATORY USE ONLY)	
	Ref. 049PCMDR8	•6	ORANGE CAPS MICRO SIZE (VERY ELASTIC RETENTION) (UNDERSIZED INTERNAL DIAMETER FOR WORN OUT SPHERE OR FOR Ø1,6 mm) - 200g	
	Ref. 040CRMDR8	• 6	AQUA CAPS MICRO SIZE (STRONG RETENTION) (UNDERSIZED INTERNAL DIAMETER FOR WORN OUT SPHERE OR FOR Ø1,6 mm) - 1100g	
6	Ref. 042AMM	• 2	MONO OT BOX MICRO SIZE (CASTABLE HOUSINGS FOR CAPS)	
	Ref. 041CAM	• 2	STAINLESS STEEL HOUSINGS MICRO SIZE	
	Ref. 041CTM	• 2	TITANIUM HOUSINGS MICRO SIZE	
	Ref. 101PPM	• 2	PLASTIC POSITIONERS MICRO SIZE	
	Ref. 044PPM	• 2	MICRO SPHERE ANALOGS	

### OT CAP ASSORTMENT SET



KIT NORMAL

Ref. 074KAN

KIT MICRO

Ref. 074KAM

**SET CONTENT** 

- 8
- 8
- 6
- CONTENT

  CLEAR CAPS STANDARD RETENTION
  PINK CAPS SOFT RETENTION
  YELLOW CAPS EXTRA SOFT RETENTION
  GREEN CAPS VERY ELASTIC RETENTION
  GOLD CAPS (EXTRA RESILIENT) SLIGHTLY ELASTIC
  SILVER CAPS (EXTRA RESILIENT) ELASTIC AND GUMMY
  LIGHT GREEN CAPS STANDARD (FOR WORN OUT SPHERE)
  ORANGE CAPS ELASTIC (FOR WORN OUT SPHERE)
  BLACK CAPS LABORATORY USE ONLY
  TITAN CAPS

- TITAN CAPS



### **ASSORTMENT SET NORMAL & MICRO**

### **SET CONTAINS**

- SET CONTAINS

   4+4 CLEAR CAPS STANDARD RETENTION

   4+4 PINK CAPS SOFT RETENTION

   2+2 YELLOW CAPS EXTRA SOFT RETENTION

   3+3 GREEN CAPS VERY ELASTIC RETENTION

   3+3 GOLD CAPS (EXTRA RESILIENT) SLIGHTLY ELASTIC

   3+3 SILVER CAPS (EXTRA RESILIENT) ELASTIC AND GUMMY

   2+2 LIGHT GREEN CAPS STANDARD (FOR WORN OUT SPHERE)

   2+2 ORANGE CAPS ELASTIC (FOR WORN OUT SPHERE)

   2+2 BLACK CAPS LABORATORY USE ONLY

   1+1 TITAN CAPS
- 1+1 TITAN CAPS

### OT CAP ASSORTED RETENTIVE CAPS -

Ref. 074KNM



### NORMAL

Ref. 064ACN

MICRO

SET CONTENT

- ORIENI
  GREEN CAPS VERY ELASTIC RETENTION
  YELLOW CAPS EXTRA SOFT RETENTION
  PINK CAPS SOFT RETENTION
  CLEAR CAPS STANDARD RETENTION
- 6
- Ref. 064ACM • 6

### OT CAP SYSTEM ATTACHMENTS •

Ref. 038STN

Ref. 038STM



NORMAL

TITANIUM SINGLE SPHERE NORMAL SIZE + TIN

FOR LASER SOLDERING - FOR BONDING

**SET CONTENT** 

- SINGLE TITANIUM SPHERES PINK CAPS SOFT RETENTION SPHERE HOLDER
- CASTABLE SLIDING BASES



MICRO

TITANIUM SINGLE SPHERE MICRO SIZE + TIN FOR LASER SOLDERING - FOR BONDING

SET CONTENT

- SINGLE TITANIUM SPHERES • 2
- PINK CAPS SOFT RETENTION SPHERE HOLDER
- CASTABLE SLIDING BASES

### OT CAP NORMAL-MICRO TOOLS -

	Ref. 485IC • 1 EXTRACTOR/INSERTER TOOL FOR CAPS OT EQUATOR - OT CAP NORMAL/MICRO	
Marie Contraction of the Contrac	Ref. 74AC01 • 1 PARALLELOMETER MANDREL NORMAL/OT EQUATOR	
	Ref. 74AC02 • 1 PARALLELOMETER MANDREL MICRO	
	Ref. 74AC03  • 1 PARALLELOMETER MANDREL TECNO NORMAL AND MICRO SIZE (ONLY FOR OT CAP TECNO)	
	Ref. 491EC • 1 MULTIUSE REMOVAL TOOL FOR RETENTIVE CAPS	
	Ref. 008MDS  • 1 SPHERE MEASURING GAUGE (TO CHECK THE SPHERICAL ATTACHMENT DIAMETERS: 1,6-1,8-2,25-2,5mm)	
	Ref. 080RCN • 1 REAMER TOOL FOR OT CAPS - NORMAL SIZE	
	Ref. 080RCM • 1 REAMER TOOL FOR OT CAPS - MICRO SIZE	
	Ref. 082ATN · 1 TOOL FOR TESTING CAP RETENTION - NORMAL SIZE	
	Ref. 083ATM • 1 TOOL FOR TESTING CAP RETENTION - MICRO SIZE	
NODMAN	MICDO ACCECCODIEC AND GO 2	

### NORMAL - MICRO ACCESSORIES AND Ø2.2mm =

	Ref. 155CC	• 6	CONNECTORS
	Ref. 187CRS	• 2	CONCAVE SPHERES AVAILABLE: Ø1.8 - 2.2 - 2.5mm
	Ref. 187CRE	• 2	CONCAVE OT EQUATOR
<b>?</b>	Ref. 189SRS	• 2	SOLID SPHERES AVAILABLE: Ø1.8mm
0	Ref. 100PD	• 10	PROTECTIVE DISKS (internal diameter 1,6mm)
9	Ref. 100PDR	• 10	PINK PROTECTIVE DISKS (internal diameter 2,9mm)
	Ref. 161SCMN  MICRO  Ref. 161SCMM	• 4	OT CAP NORMAL/MICRO SPACERS FOR HOUSING
	Ref. 038GCS	• 2	CASTABLE SLIDING BASES
	Ref. 044CAI22	• 2	IMPRESSION COPINGS Ø2,25mm SPHERES WITH INTERCHANGEABLE CAP
•	Ref. 044PP22	• 2	MICRO ANALOGS Ø2.2mm

### OT STRATEGY SPARE PARTS



Ref. 154PCS 4 OT STRATEGY STANDARD MALES



Ref. 054PLS 4 OT STRATEGY HIGH MALES



 2 CASTABLE STEADY **Ref. 154STY** 



STRATEGY YELLOW CAPS FOR DUPLICATION TECHNIQUE Ref. 047CSD AY (EXTRA SOFT RETENTION) - 500g



STRATEGY PINK CAPS FOR DUPLICATION TECHNIQUE Ref. 047CSD SN (SOFT RETENTION) - 900g



STRATEGY CLEAR CAPS FOR DUPLICATION TECHNIQUE Ref. 047CSD (STANDARD RETENTION) - 1300g



STRATEGY BLACK CAPS FOR DUPLICATION TECHNIQUE Ref. 047CLS (FOR LABORATORY USE ONLY)



Ref. 047CSDDR8

AQUA STRATEGY CAPS FOR DUPLICATION TECHNIQUE (UNDERSIZED INTERNAL DIAMETER FOR WORN SPHERE OR 1.7mm SPHERE)



Ref. 045CSM AY

 4 STRATEGY YELLOW CAPS FOR S/STEEL HOUSING (EXTRA SOFT RETENTION) - 500g (TO BE USED ONLY WITH HOUSINGS)



Ref. 045CSM SN

STRATEGY PINK CAPS FOR S/STEEL HOUSING (SOFT RETENTION) - 900g (TO BE USED ONLY WITH HOUSINGS)



Ref. 045CSM

STRATEGY CLEAR CAPS FOR S/STEEL HOUSING (STANDARD RETENTION) -1300g (TO BE USED ONLY WITH HOUSINGS)



Ref. 045CLS

STRATEGY BLACK CAPS FOR S/STEEL HOUSING (FOR LABORATORY USE ONLY)



Ref. 045CSMDR8

AQUA STRATEGY CAPS FOR S/STEEL HOUSING (UNDERSIZED INTERNAL DIAMETER FOR WORN SPHERE OR 1.7mm SPHERE)




Ref. 046CSM


2 STRATEGY STAINLESS STEEL HOUSINGS



**Ref. 1010TS** 

 2 STRATEGY PLASTIC POSITIONERS (FOR SOLDERING STAINLESS STEEL HOUSINGS)





### OT VERTICAL SPARE PARTS -





**Ref. 472ICV** • 1 INSERTION TOOL FOR OT VERTICAL CLIPS

OT REVERS	E 3 SPARE PARTS		
	<b>Ref. 120MR2</b> • 2	RETENTIVE MALES IN TITANIUM + NYLON (COLOR CLEAR/ STANDARD RETENTION)	
	<b>Ref. 041CAM</b> • 2	S/STEEL HOUSINGS MICRO SIZE	
Ĭ	<b>Ref. 121AAR</b> •2	INOX ROOTS ANALOGS	
	Ref. 121AM · 2	PLASTIC MANUAL TOOLS	
₽ [	Ref. 121PR · 2	PLASTIC POSITIONERS FOR DUPLICATION	
	Ref. 100PDR · 10	PINK PROTECTIVE DISKS	
	<b>Ref. 124IMR</b> •1	METAL INSERTING TOOL FOR OT REVERSE/3 MALE	
	<b>Ref. 72CBR</b> ·1	S/STEEL PARALLELOMETER KEY FOR CASTABLE BAR	
The state of the s	<b>Ref. 72CFR</b> ·1	BRASS PARALLELOMETER KEY FORTITANIUM FEMALE	

### **OT LOCK SPARE PARTS**



**Ref. 882CS** 

• 1 SLIDING KEY (ONLY FOR TITANIUM - CASTABLE LOCKING PIN)



**Ref. 882AF** 

• 1 FLEXIBLE RING



**Ref. 882CB** 

• 1 LOCKING KEY



Ref. 882CPT

• 1 TITANIUM PREFABRICATED HOUSING (ONLY FOR TITANIUM - TITANIUM ADJUSTABLE LOCKING PIN)



**Ref. 882PPC** 

• 1 CERAMIC PIN



Ref. 882CLSCT

• 1 ADJUSTABLE SLIDING KEY (ONLY FOR TITANIUM ADJUSTABLE LOCKING PIN)



**Ref. 882DCC** 

CASTABLE SPACER RINGS (ONLY FOR TITANIUM ADJUSTABLE LOCKING PIN)

### **OT LOCK ACCESSORIES**



**Ref. 882CG** 

• 1 CONICAL GUIDE



**Ref. 882CAS** 

• 1 UNLOCKING TOOL

### SPHERE REPAIR SETS •



Ø AVAILABLE: 1,8 - 2,25 - 2,5 Ref. 087CRS

**CONCAVE REPAIR SPHERE** 

### SET CONTENT

- CONCAVE SPHERES IN TITANIUM TIN COATED • 2
- PINK CAPS SOFT RETENTION INSERTING TOOL
- GAUGE AND STRIP HOLDER



Ref. 087CRE

### **CONCAVE REPAIR OT EQUATOR**

### **SET CONTENT**

- CONCAVE OT EQUATOR IN TITANIUM TIN COATED PINK CAPS SOFT RETENTION INSERTING TOOL
- 2 2
- GAUGE AND STRIP HOLDER



**Ref. 089SRS** 

### **SOLID REPAIR SPHERE**

### **SET CONTENT**

- SOLID SPHERES 1.8mm IN TITANIUM TIN COATED • 2
- PINK CAPS SOFT RETENTION PROTECTIVE DISKS
- KEY


### OT LEVER -

**UNIVERSAL LEVER** 



**Ref. 7610TL** 

• 2 OT LEVER FOR REMOVAL OF PROSTHESES AND ALIGNERS



### **OT LEVER BOX**



**Ref. BOL** 

### **SET CONTENT**

- 6 OT LEVER PACKAGES1 OT LEVER PACKAGING DISPLAY STAND

### OT CEM -

**COMPOSITE CEMENT** 



### **SET CONTENT**

- TWO-COMPONENT SYRINGE 5ml
- •10 AUTOMIX TIPS



**Ref. TOC** 

•10 AUTOMIX TIPS

### **BROKEN SCREW EXTRACTOR KIT**

**BROKEN SCREW EXTRACTOR SETS ARE** MANUFACTURED FOR MOST IMPLANT BRANDS



**Ref. 680** 



- SET CONTENT

  1 CLAW REAMER BUR
  REVERSE CUTTING BUR
  MANUAL EXTRACTOR
  CENTERING DEVICE
- **Ref. 680FS**
- REVERSE CUTTING BUR
- Ref. 680FA\*
- **CLAW REAMER BUR**

\*Specify type of implant and diameter

THE BROKEN SCREW EXTRACTOR SET MAY HAVE SPACE PROBLEMS IN THE REAR SECTION, WHICH COULD AFFECT ITS PERFORMANCE. IT'S IMPORTANT TO ANALYZE THE SITUATION BEFORE USING THE SET.



# **DIRECT SERVICES**

A telephone service is available during office hours, it furnishes direct answers for technical questions. You can address to Rhein83 for knowing the address of foreign distributors

Tel. +39 051 244510 - +39 051 244396 - Fax +39 051 245238

Rhein83 website www.rhein83.com

General informations info@rhein83.com

Sales department international.orders@rhein83.it

Technical support support@rhein83.com

Marketing department direzione@rhein83.it





Via ZAGO, 10/ABC 40128 - BOLOGNA (ITALY)

Tel. (+39) 051 244510 - (+39) 051 244396 Fax (+39) 051 245238







