



let's
keep
connected



REV 01 02/2019

 polydentia
swiss manufacture 

Polydentia SA
via Cantonale 47
6805 Mezzovico-Vira
Switzerland

phone: +41 91 946 29 48
fax: +41 91 946 32 03

info@polydentia.ch
www.polydentia.ch



my junior kit

the first
sectional
matrix system
for paediatric
dentistry

sectional
matrix
systems



 polydentia
swiss manufacture 



new!
the first
sectional matrix system
for paediatric dentistry

myjunior kit

myJunior kit is the first sectional matrix system for paediatric dentistry, especially developed to suit smaller primary teeth and young permanent teeth. It enables dentists to achieve predictable restorations and ideal contact surfaces. The kit offers child-friendly solutions that help saving chair-time and reduce the distress of children, allowing them to quickly return to their favourite activities.



myTines Junior
Autoclavable and replaceable plastic extremities that fit over the wedge and provide excellent matrix adaptation



Silicone Tubes
Soft extremities that provide extra grip and adaptation



myWedge
Hollow V-shaped, high-tech plastic wedge



Junior Quickmat Steel Matrices

use only with rubber dam

kit

myJunior Kit - REF 6330

- 1 pc. myClip Junior
- 1 pc. myRing Junior
- 5 pairs Cylindrical Silicone Tubes
- 5 pairs myTines Junior
- 40 pcs. LumiContrast Sectional Matrices, Premolar 0.04mm / 5mm
- 40 pcs. Junior LumiContrast Sectional Matrices, 0.04mm / 3.75mm
- 50 pcs. Quickmat Sectional Matrices, Premolar 0.04mm / 5mm
- 50 pcs. Junior Quickmat Sectional Matrices, 0.04mm / 3.75mm
- 50 pcs. Pinky Wood Wedges, 11mm (XS)
- 50 pcs. Wood Wedges, 12mm (S)
- 50 pcs. myWedges XS
- 50 pcs. myWedges S


refill

- | | | | |
|---|--|--|--|
|
REF 6335
myClip Junior,
2 pcs. |
REF 6336
myRing Junior,
2 pcs. |
REF 6333
myTines Junior,
15 pairs |
REF 5711
Silicone Tubes,
15 pairs |
|
REF 6708
Junior LumiContrast Sectional
Matrices, 0.04mm / 3.75mm
70 pcs. |
REF 6705
LumiContrast Sectional Matrices,
Premolar 0.04mm / 5mm
70 pcs. |
REF 5708
Junior Quickmat Sectional
Matrices, 0.04mm / 3.75mm
100 pcs. |
REF 5705
Quickmat Sectional Matrices,
Premolar 0.04mm / 5mm,
100 pcs. |
|
REF 5011
Pinky Wood Wedges,
11 mm (XS), 200 pcs. |
REF 5012
Wood Wedges,
12 mm (S), 100 pcs. |
REF 7101
myWedge,
Extrasmall (XS), 100 pcs. |
REF 7102
myWedge,
Small (S), 100 pcs. |

myClip Junior and myRing Junior are paediatric sectional matrix rings that feature minimally invasive, safe to use, bright-coloured solutions that allow easier restorative treatment on young teeth of low crown height. The junior sectional matrix rings combine stability with the right separation force enabling dentists to achieve ideal contact surfaces. In addition, the extra grip provided by the flexible extremities myTines Junior and rounded silicone tubes ensure a better adaptation of the sectional matrices, that helps preventing annoying composite flashes and overhangs and saving valuable chair-time.

myClipjunior

All in one paediatric sectional matrix ring with integrated forceps



ADVANTAGES

- All in one sectional matrix ring. No need for intimidating forceps
- Long-lasting ring, made out of Nickel Titanium
- Optimal separation force to achieve flawless contact surfaces
- Quick, easy, one-hand ring positioning
- Autoclavable and replaceable plastic myTines Junior extremities



"I regard myClip Junior as an optimal solution for paediatric treatments as it allowed me to achieve an ideal restoration with less efforts and with considerable chair-time savings."


Dr. Marina Papachroni
Greece



for more information


myRingjunior

Paediatric sectional matrix ring that fits a wide variety of clinical situations



ADVANTAGES

- Ideal size to fit on deciduous teeth in children's mouth
- Long-lasting ring made out of Nickel Titanium
- Optimal separation force to achieve flawless contact surfaces
- Excellent matrix adaptation and retention
- Autoclavable and replaceable plastic myTines Junior extremities



"Children restorative treatments won't be tricky business anymore thanks to myRing Junior and the paediatric solutions from Polydentia. Paediatric restorations are easy and efficient thanks to the small dimension of the sectional matrix ring and the junior matrices. The system ensures a perfect fit on milk teeth and prevent time-consuming composite flashes. Restorations are also more enjoyable as my little patients don't have to withstand prolonged restoration times."

Dr. Martina Sabová
Slovakia



for more information